



25404



25410



25408



25409

**Detector Supplies** for Waters HPLC & ACQUITY UPLC® Systems

25405

Description	Model #	Similar to Waters Part #	qty.	cat.#
Xenon Lamp (without holder or mirror)	470	—	ea.	25404
Xenon Lamp	474	—	ea.	25405
Deuterium Lamp, Long Life (2,000 hours)*	486	WAT080678	ea.	25410
Deuterium Lamp	996, 2996	WAT052586	ea.	25408
Deuterium Lamp	2487, 2488	WAS081142	ea.	25409

\*Standard lamps have nominal 1,000-hour life.

## Thomson SINGLE StEP® Filter Vials

**Sample filtration that's economical, eco-friendly, and fast!**



- Easy-to-use vials offer fast sample filtration and require only a squeeze of your fingers.
- Color-coded caps allow easy identification of 0.2 µm or 0.45 µm membranes in PVDF, PTFE, PES, or nylon.
- Pre-slit PTFE/silicone caps help eliminate broken autosampler needles and cored septa.
- Low dead volume units feature rugged polypropylene vial and insert with 450 µL loading capacity.
- Fit most standard 12x32 mm autosamplers, including UHPLC instruments.

See **page 410**.
[www.restek.com/singlestep](http://www.restek.com/singlestep)
