



Universal 10-32 PEEK Column Connectors and Plugs

Universal PEEK connectors allow easy installation of all $\frac{1}{16}$ " tubing, including stainless steel, up to 5,000 psi.

Description	qty.	cat.#
PEEK Column Connector (beige, round body)	10-pk.	25015
PEEK Column Plug (black)	10-pk.	25016
PEEK Finger-Tight Fittings (blue, flat-sided)	10-pk.	25324



Rheodyne® Style Nut and Ferrule

Replacement long nut for connecting stainless steel tubing to a Rheodyne® 6-port valve or other Rheodyne® part.

Description	qty.	cat.#
$\frac{1}{16}$ " Rheodyne Style Nut	10-pk.	25095
$\frac{1}{16}$ " Rheodyne Style Ferrule	10-pk.	25096



Zero Dead Volume Valco® Internal Union

Ends of tubing seat squarely at bottoms of fitting details. Made of 300-series stainless steel. For use with $\frac{1}{16}$ " OD tubing. Stainless steel ferrules included.

Description	Union Bore	Valco #	qty.	cat.#
Internal Union	0.15 mm	ZU1XC	ea.	20147
Internal Union	0.25 mm	ZU1C	ea.	20148
Internal Union	0.75 mm	ZU1	ea.	20149
Internal Union	$\frac{1}{16}$ "	ZU1T	ea.	20150



PEEK Union Connector

Quickly and reliably connect two pieces of $\frac{1}{16}$ " tubing. 0.3 mm union bore. End fittings included.

Description	qty.	cat.#
PEEK Union Connector $\frac{1}{16}$ "	2-pk.	25323



LC-MS Gas Generators From Restek

Affordable

Generators are cost-effective and quickly pay for themselves.

Reliable

Dependably produce a continuous supply of high-purity carrier and fuel gas.

Safe

Eliminate hazardous cylinders and free up valuable lab space.

See page 345.

www.restek.com/gas