



25336



25328

QuickSplit Post-Column Flow Splitters

Fluid resistor technology eliminates adjustments to capillary tubing for optimizing split ratio. Wide range of interchangeable resistors available.

Fixed Flow Splitters for HPLC & LC-MS

- Split ratio not affected by changes in viscosity or pressure.
- High operating pressure limit: 10,000 psi (68,948 kPa).
- Low dead volume—negligible effect on analyte bandwidth.
- Total flow: 0.1–5.0 mL/min.

Adjustable Flow Splitter

- Adjustable metering valve gives convenient control of split ratio.
- Split ratio not affected by changes in viscosity or pressure.
- High operating pressure limit: 5,000 psi (34,474 kPa).
- Low dead volume—negligible effect on analyte bandwidth.
- Total flow: 0.1–5.0 mL/min.

Description	Split Ratio	qty.	cat.#
Binary Fixed, Post Column	100:1	ea.	25326
	50:1	ea.	25327
	20:1	ea.	25328
	10:1	ea.	25329
	5:1	ea.	25330
Replacement Fixed Resistor Set, Post Column	100:1	ea.	25331
	50:1	ea.	25332
	20:1	ea.	25333
	10:1	ea.	25334
	5:1	ea.	25335
Adjustable, Post Column	5:1 to 100:1	ea.	25336
	1:1 to 20:1	ea.	25337
	50:1 to 1000:1	ea.	26416
	5:1 to 100:1	ea.	25338
Replacement Adjustable Resistor Set, Post Column*	1:1 to 20:1	ea.	25339

*QuickSplit post-column flow splitter serial number required for replacement resistor set to ensure compatibility. Every resistor is manufactured per serial number. The serial number is located on top of the splitter.



LCLocker LC Organizer




Deluxe BenchBooster Organizer



Mini pHPerch Storage Unit



LC 30-Column Storage Cabinet



Book Holders



Open Supply Bins, 13-bin unit



TopLoader Balance-Bank Storage Unit



Glove Box Dispensers

LC Organizers

Description	Dimensions	qty.	cat.#
LCLocker	24 x 12 x 6"	ea.	25149
BenchBooster	24 x 7 x 12"	ea.	25150
Mini pHPerch	13 x 12 x 6"	ea.	25147
TopLoader	12 x 12 x 7"	ea.	25148
30 Column Cabinet	17 x 15 x 3"	ea.	25159
Book Holder, 0.75" ID (Small)	4 1/2 x 6 x 1 1/4"	ea.	25151
Book Holder, 1.5" ID (Large)	4 1/2 x 6 x 1 3/4"	ea.	25152
Open Supply Bin	24 x 12 x 10" (4-Bin Unit)	ea.	25153
Open Supply Bin	12 x 16 x 10.5" (5-Bin Unit)	ea.	25154
Open Supply Bin	12 x 12 x 7.5" (13-Bin Unit)	ea.	25155
Glove Box Dispenser	5 1/8 x 11 x 4 1/4" (Single)	ea.	25156
Glove Box Dispenser	10 1/8 x 11 x 4 1/4" (Double)	ea.	25157
Glove Box Dispenser	15 1/8 x 11 x 4 1/4" (Triple)	ea.	25158

*Please note: Accessories and supplies in photographs are not included.