

**With Hub-Cap Bottle Tops:**  
secure fit keeps mobile phase  
lines where they belong

**Without Hub-Cap:**  
messy, loose-fitting  
tape

Bottle not included.

### Hub-Cap 4-Liter Bottle Tops

Hub-Cap bottle tops are a great way to neatly keep your mobile phase lines where they belong. Use them instead of plastic paraffin film, aluminum foil, or tape on your mobile phase reservoirs.



Description	qty.	cat.#
Hub-Cap (assembly of the bottle cap and plug)	kit	26541
Hub-Cap Multi-Pack	3-pk.	26542

### Hub-Cap Adaptors

Allow the use of the Opti-Cap® Top with 4-liter solvent bottles.

Description	qty.	cat.#
Hub-Cap Adaptor (allows use of the Opti-Cap w/4 L solvent bottles)	ea.	26538
Hub-Cap Adaptor Multi-Pack	3-pk.	26539
Hub-Cap Adaptor and Opti-Cap (for use w/GL-45 solvent bottles)	kit	26540

### Related Items and Replacement Parts for Hub-Cap 4-Liter Bottle Tops

Description	qty.	cat.#
Mobile Phase Sparge Filter, 2 µm	ea.	25311
Mobile Phase Sparge Filter, 10 µm	ea.	25312
PTFE Tubing, 1/8" OD x 0.094" (2.4 mm) ID x 3 m	3 m	25307
PTFE Tubing, 1/8" OD x 0.063" (1.6 mm) ID x 3 m	3 m	25306
PEEK Plug, 1/4"-28 threads	3-pk.	25319

### Opti-Cap® GL-45 Bottle Tops

The most economical way to helium-sparge and deliver LC mobile phases. Opti-Cap® tops fit all standard GL-45 bottles and have two 1/8" holes and one 1/16" hole for tubing. All three openings are designed for threaded PEEK plugs.



Description	qty.	cat.#
Opti-Cap (Cap and PEEK Plug)	ea.	25300
Opti-Cap Kit	kit	25301
Opti-Cap Kit with 1 L Bottle	kit	25302
Opti-Cap Kit with 2 L Bottle	kit	25303

### Opti-Cap® Adaptors

Allow the use of the Hub-Cap Top with GL-45 solvent bottles.

Description	qty.	cat.#
Opti-Cap Adaptor (allows use of Hub-Cap w/GL-45 solvent bottles)	ea.	27197
Opti-Cap Adaptor Multi-Pack	3-pk.	27198
Opti-Cap Adaptor w/Hub-Cap	kit	26551



### Related Items and Replacement Parts for Opti-Cap® GL-45 Bottle Tops

Description	qty.	cat.#
Mobile Phase Sparge Filter, 2 µm	ea.	25311
Mobile Phase Sparge Filter, 10 µm	ea.	25312
PTFE Tubing, 1/8" OD x 0.094" (2.4 mm) ID x 3 m	3 m	25307
PTFE Tubing, 1/8" OD x 0.063" (1.6 mm) ID x 3 m	3 m	25306
PEEK Plug, 1/4"-28 threads	3-pk.	25319
1 L Graduated Safety-Coated Bottle – GL-45 threads	ea.	25304
2 L Graduated Safety-Coated Bottle – GL-45 threads	ea.	25305

## Which Bottle Top Do I Need?

### I have a 4-Liter bottle...

**Option 1:** Hub-Cap Top (cat.# 26541)

or

**Option 2\*:** Opti-Cap® Top (cat.# 25300)  
plus Hub-Cap Adaptor (cat.# 26538)

\*Also available as a kit (cat.# 26540)

#### Option 1



#### Option 2\*



\*Also available as a kit  
(cat.# 26540)

### I have a GL-45 bottle...

**Option 1:** Opti-Cap® Top (cat.# 25300)

or

**Option 2\*:** Hub-Cap Top (cat.# 26541)  
plus Opti-Cap® Adaptor (cat.# 27197)

\*Also available as a kit (cat.# 26551)

#### Option 1



#### Option 2\*



\*Also available as a kit  
(cat.# 26551)