

Parker Balston® Nitrogen Gas Generators for LC-MS

- Turn compressed air into ultra-pure nitrogen (up to 99.5%).
- Flows from 1 to 44 L/min.
- Models N2-04, N2-14, N2-22, and N2-35 require no electricity.
- Safe, reliable, low maintenance.
- Maintenance kits include replacement filters.



CE
22130

CE
22129

Specifications	NitroFlow Lab	N2-04	N2-14	N2-22 or N2-22A	N2-35 or N2-35A
Maximum Nitrogen Flow Rate:	32 L/min	8 L/min	36 L/min	N2-22: 44 L/min N2-22A: 29 L/min	44 L/min
Nitrogen Purity:	99.50%	99%	95.0%–99.5%	99%	99%
Min/Max Inlet Pressure:	N/A	60 psig/145 psig	60 psig/145 psig	60 psig/145 psig	60 psig/145 psig
Electrical Requirements:	120 VAC/60 Hz	None	None	N2-22: None N2-22A: 120 VAC/60 Hz	N2-35: None N2-35A: 120 VAC/60 Hz
Dimensions:	27.6" h x 35.4" w x 12.2" d (70 cm x 90 cm x 31 cm)	11" h x 13" w x 16" d (27 cm x 34 cm x 41 cm)	50" h x 16" w x 16" d (127 cm x 41 cm x 41 cm)	50" h x 16" w x 16" d (127 cm x 41 cm x 41 cm)	50" h x 16" w x 16" d (127 cm x 41 cm x 41 cm)
Shipping Weight:	205 lb (93 kg)	43 lb (20 kg)	75 lb (34 kg)	N2-22: 101 lb (46 kg) N2-22A: 106 lb (48 kg)	N2-35: 115 lb (52 kg) N2-35A: 119 lb (54 kg)

Nitrogen Generators for LC-MS	Model #	qty.	cat.#
Nitrogen Generator for LC-MS	NitroFlow Lab Model, 32 L/min max. flow	ea.	22129
Nitrogen Generator for LC-MS	N2-04 Model for ELSD, 8 L/min max. flow	ea.	22130
Nitrogen Generator for LC-MS	N2-14 (general purpose) 36 L/min max. flow	ea.	20677
Nitrogen Generator for LC-MS	N2-22 Model, 44 L/min max. flow	ea.	22131
Nitrogen Generator for LC-MS	N2-22A Model, 29 L/min max. flow	ea.	22132
Nitrogen Generator for LC-MS	N2-35 Model, 44 L/min max. flow	ea.	22133
Nitrogen Generator for LC-MS			
Includes: Oxygen analyzer and audible alarm (monitors oxygen in nitrogen stream and signals high or low concentrations)	N2-35A Model, 44 L/min max. flow	ea.	22134
Maintenance Kits	Model #	qty.	cat.#
Maintenance Kit	for NitroFlow Lab	kit	22156
Includes: carbon filter			
6-Month Maintenance Kit*			
Includes: 1st stage prefilter, 2nd stage prefilter, and final filter (1 each)	for N2-14, N2-14A, 75-T2, 75-T20NA	kit	21648
6-Month Maintenance Kit With Carbon Filter*			
Includes: 1st stage prefilter, 2nd stage prefilter, hydrocarbon scrubber, and final filter (1 each)	for N2-14, N2-14A, 75-T2, 75-T20NA	kit	22135
Generator Maintenance Kit	for HPZA-3500, HPZA-7000, HPZA-18000, HPZA-30000, and 75-80 (Zero Air); N2-04 (Nitrogen); and TOC-1250 (TOC)	kit	21647
Includes: 1st and 2nd stage prefilters (1 each) and 1 final filter			

*The manufacturer recommends maintenance for this unit every 6 months.

International power cords are available. Contact Customer Service or your Restek® representative to order.



CE
20677

Need to clean up your nitrogen stream for LC-MS?

Check out our Restek® Super Clean® gas trapping system.

See **page 288**.

www.restek.com/gas

