



22292

Syringe for Hitachi LC Autosamplers

SGE

Volume	SGE Model	cat.#	qty.	Restek cat.#
500 µL, M10 X1 Thread	500C-HITACH1	007660	ea.	22292



22297

Syringes for PerkinElmer LC Autosamplers

SGE

Volume	SGE Model	cat.#	qty.	Restek cat.#
250 µL, 1/4-28 UNF Thread	250D-CX-GT	006995	ea.	22297
500 µL, 1/4-28 UNF Thread	500D-CX-GT	007995	ea.	22298

Syringe for Waters WISP LC Autosamplers

PTFE-tipped plungers

Hamilton

Volume	Hamilton Model	cat.#	qty.	Restek cat.#
250 µL, 1/4-28 UNF Thread	1725	80024	ea.	24529



22294

SGE

Volume	SGE Model	cat.#	qty.	Restek cat.#
250 µL, 1/4-28 UNF Thread	250D-WISP	006690	ea.	22294

Syringes for CTC LC Autosamplers

Hamilton



22746

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton Model	cat.#	qty.	Restek cat.#
10 µL	N	22s	2"/51 mm	3	701N	203073	ea.	22743
25 µL	N	22s	2"/51 mm	3	1702N	203075	ea.	22744
100 µL	N	22s	2"/51 mm	3	1710N Slim Line*	203077	ea.	22745
100 µL	N	22	2"/51 mm	3	1710N	203235	ea.	22746
250 µL	N	22	2"/51 mm	3	1725N	203079	ea.	22747
500 µL	N	22	2"/51 mm	3	1750N	203349	ea.	22748

*Barrel OD = 6.7 mm; all other 25 µL and 100 µL syringes have a 7.9 mm barrel OD.



22737

SGE

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE Model	cat.#	qty.	Restek cat.#
10 µL	F	22s	2"/51 mm	LC	10F-CTC-LC	002710	ea.	22737
100 µL*	R	22s	2"/51 mm	LC	100R-C/T-GT-LC	005330	ea.	22741
500 µL*	F	22s	2"/51 mm	LC	500F-CTC-GT-LC(0.41)	007720	ea.	22742

*Gas-tight syringe.

Guide to Needle Termination Codes

Hamilton:

(N) Cemented Needle

SGE:

(F) Fixed Needle

(R) Removable Needle

also available

More syringes!

See page 394.

