

CHROMTECH LAMP GUIDE

January 2016

BUCK (SEE CECIL/JENWAY)			
BLC-20	LH-94	D2	
CAMAG			
TLC II SCANNER	L-285	D2	
CARLO ERBA/THERMO SCIENTIFIC			
	L-280	D2	
CARY/VARIAN			
1,3,4,5E,100,300,400,500	L-230	D2	
1,3,4,5E,100,300	L-930	W	
50	L-1045	P/XE	
400	L-110	W	
219	L-205	D2	
219	L-925	W	
500	L-1015	W	
4000, 5000, 6000	* LH-240	D2	
CECIL (BUCK)			
CE 1020,2020,2040,2041,2501,2502	L-115	D2	
CE 1020,2020,2040,2041,2501,2502	L-845	W	
CE 1200, 1220	LH-20		
CE 1200, 1220	LH-20	D2	
CE 4000 thru 9000 SERIES	L-115	D2	
CE 4000 thru 9000 SERIES	L-845	W	
COLE-PALMER			
01953-00	LH-40	D2	
COMETRO TECHNOLOGY LTD			
6000 PUW			INQUIRE
DIONEX/LC PACKING/GYNKOTEK/THERMO (SEE THERMO SCIENTIFIC)			
D-STAR			
DDW-10,DLC-20,DVW-10	LH-94	D2	
DUPONT			
850	L-300	D2	
8800	L-280	D2	
DYAX (SEE BIOTAGE)			
ECOM			
LCD2073	LH-60	D2	
SAPPHIRE,OPAL,TOPAZ	LH-60	D2	
FLASH 06, FLASH 10	LH-60	D2	
LCD2071.2, LCD2071.3	L-735	HG	
* MUST RETURN HOLDER			
EG & G (SEE PERKIN-ELMER)			
EMD MILLIPORE/MERCK MILLIPORE (SEE HITACHI)			
L2455	LH-56	D2	
NOVA 30A,60A			INQUIRE
PHARO 100,300			INQUIRE
MILLIPORE WATER SYSTEM		HG	
ERMA/ERC			
7511, 7512, 7515, 7522, 7530, 7535	L-610	W	
ESA			
520, 522	LH-40	D2	
FINSTRUMENTS			
INQUIRE			
FISHER (SEE THERMO SCIENTIFIC)			
GBC			
AA LAMPS	L-200	D2	
CINTRA 906,908,916	L-200	D2	
CINTRA 916	L-840	W	
CINTRA 1155, 1255	LH-910	P/XE	
CINTRA 1200	L-280	D2	
CINTRA 1205,1206,1210	LH-40	D2	
CINTRA 5100	LH-82	D2	
GE HEALTHCARE/AMERSHAM BIOSCIENCES			
PHARMACIA/LKB/BIOCHROM			
213,821,402,141			INQUIRE
2151	LH-20	D2	
4050, 4051, 4053, 4050	L-210	D2	
4050, 4060	L-835	W	
AKTA SERIES/FPLC/PROTEIN PURIFICATION SYSTEM			
AKTA BASIC,ETAN,EXPLORER 100, 105	* LH-990	P/XE	
AKTA EXPRESS, CROSS FLOW	L-1010	HG	
AKTA FPLC	L-1010	HG	
AKTA MICROSYSTEM, PILOT	* LH-990	P/XE	
AKTA PURIFIER 10, 100	* LH-990	P/XE	
AKTA PURIFIER 10, 100	L-1010	HG	
AKTA PURE 25, 150	* LH-990		
GENEQUANT II	+++ LH-975	D2/W	
GENEQUANT 100, 1300			INQUIRE
GENEQUANT PRO	+++ LH-975	D2/W	
* Comes with Blocking Filter			
+++ MUST RE-USE HOLDER			