

CHROMTECH LAMP GUIDE

January 2016

| GE HEALTHCARE (cont.) | | | | |
|--|------------|---------|------|---------|
| AKTA SERIES/FPLC | | | | |
| NOVASPEC II | L-835 | W | | |
| NOVA SPEC III, PLUS | | W | | INQUIRE |
| RNA/DNA CALCULATOR | +++ LH-975 | D2/W | | |
| NANOVUE SERIES | | | | |
| | | P/XE | | INQUIRE |
| | | | | |
| ULTROSPEC SERIES | | | | |
| ULTROSPEC II,III,PLUS | L-210 | D2 | | |
| ULTROSPEC II,III,PLUS | L-835 | W | | |
| ULTROSPEC 500 PRO | LH-125 | D2/W | | |
| ULTROSPEC 500 PRO | L-935 | W | | |
| ULTROSPEC 1000 E, 1000 I | +++ LH-975 | D2/W | | |
| ULTROSPEC 1000 E | L-935 | W | | |
| ULTROSPEC 1100 PRO | +++ LH-975 | D2/W | | |
| ULTROSPEC 1100 PRO | L-935 | W | | |
| ULTROSPEC 2000, 3000, 3000 PRO | LH-125 | D2/W | | |
| ULTROSPEC 2000, 3000 | L-935 | W | | |
| ULTROSPEC 2100 PRO, 3100 PRO | L-640 | P/XE | | |
| ULTROSPEC 3300, 4000, 4300 | LH-125 | D2/W | | |
| ULTROSPEC 3300, 4000, 4300 | L-935 | W | | |
| ULTROSPEC 5300 PRO, 6300 PRO | LH-125 | D2/W | | |
| UV 900, UVIS 920 | * LH-990 | P/XE | | |
| UV-1 | * L-1030 | HG | | |
| UPC 900 | L-1010 | HG | | |
| | | | | |
| GENTECH (SEE BEIJING) | | | | |
| | | | | |
| GILFORD | | | | |
| 200, 220, 222, 230, 240, 250, 252, 260 | L-205 | D2 | | |
| 200, 220, 222, 230, 240, 250, 252, 260 | L-875 | W | | |
| 2000, 2400, 2600 | L-205 | D2 | | |
| 2000, 2400, 2600 | L-875 | W | | |
| FLUORO IV | L-801 | XE | | |
| RESPONSE I | L-210 | D2 | | |
| RESPONSE I | L-900 | W | | |
| STASAR I, II, III, SYVA S-3 | L-705 | W | | |
| | | | | |
| +++ MUST RE-USE HOLDER | | | | |
| * Comes with Blocking Filter | | | | |
| * MUST RETURN SPENT LAMP & HOLDER | | | | |
| | | | | |
| GILSON | | | | |
| 115, 116, 117 | | LH-150 | D2 | |
| 118UV, 119 | | LH-151 | D2 | |
| 122 | | L-800 | XE | |
| 151, 152, 155, 156 | | LH-151 | D2 | |
| 160, 170, 171, 172 | | | | INQUIRE |
| GRACE/ALLTECH | | | | |
| 200, 201, 205 | | LH-40 | D2 | |
| 200 | | LH-760 | W | |
| 305 | | LH-910 | P/XE | |
| REVELERIS UV1, UV2 | | LH-40 | D2 | |
| GROTON | | | | |
| FD100 | | L-890 | P/XE | |
| GYNKOTEK (SEE DIONEX) | | | | |
| | | | | |
| HACH | | | | |
| DR/2700,2800,3800,3900 | | INQUIRE | W | |
| DR/2000, 2010 | | | | |
| DR/4000U,DR/4000V | | L-730 | W | |
| DR/4000U | | LH-110 | D2 | |
| 48000, 48100 | | L-730 | W | |
| DR/5000, 6000 | | LH-78 | D2 | |
| DR/5000, 6000 | | L-780 | W | |
| 57806-18 | | | W | INQUIRE |
| ASTRO | | | | INQUIRE |
| HARVARD BIOSCIENCE (SEE BIOCHROM) | | | | |
| | | | | |
| HITACHI | | | | |
| 100-10,20,30,40,60,80 | | L-280 | D2 | |
| 101, 102, 111 | | L-430 | D2 | |
| 124 | | L-210 | D2 | |
| 150-20 | | LH-51 | D2 | |
| 200, 300, 330, 340 | | LH-51 | D2 | |
| 5440 | | L-800 | XE | |
| 655A-22, -23 | | LH-51 | D2 | |
| F-1000,1050,1080,1200,2000,2500 | | L-800 | XE | |
| F-3000,3010,4000,4010,4500,7000 | | L-800 | XE | |
| GENE SPEC I | | LH-58 | D2 | |
| GENE SPEC III, V | | | | INQUIRE |
| L-2400,2420,2450,2455 | | LH-56 | D2 | |
| L-2420,2450,2455 | | L-835 | W | |
| | | | | |