

CHROMTECH LAMP GUIDE

January 2016

JY (JOBIN YVON)/SPEX/HORIBA			
1934D		P/XE	INQUIRE
FLUOROLOG-2/ F111AI,F112AI		XE	INQUIRE
FLUOROLOG-2/F111XI,F112XI	L-855	XE	
FLUOROLOG-2/F1T11I,F212I,F222I	L-855	XE	
FLASHLAB			INQUIRE
FLUOROG-3,FLUOROLOG TAU-3	L-855	XE	
FLUOROMAX-2,-3,-4	L-955	XE	
SYSTEM 5000	L-60	P/XE	
KNAUER/SONNTEK			
AZURA			
2.1S, 2.1S FIBER OPTIC	A5193	D2	INQUIRE
2.1 L UVD, MWD, DAD	A5193	D2	INQUIRE
6.1 L	AZL01	D2	INQUIRE
6.1 L	AZC1L	W	INQUIRE
BLUE SHADOW			
10D, 20D, 40D, 50D	A5193		INQUIRE
SMARTLINE			
200	A4142	HG	INQUIRE
2500, 2600	A4071	D2	INQUIRE
2500, 2600	A4073	W	INQUIRE
2520	A5193	D2	INQUIRE
2550	A5194	D2	INQUIRE
2550	A5195	W	INQUIRE
2800, 2850	A4447	D2	INQUIRE
2800, 2850	A4448	W	INQUIRE
WELLCHROM			
200	A4142	HG	INQUIRE
2000, 2001	A4071	D2	INQUIRE
2001	A4073	W	INQUIRE
2500, 2501, 2600, 2650, 2800	A4071	D2	INQUIRE
2500, 2501, 2600	A4073	W	INQUIRE
2800	A4448	W	INQUIRE
NON-SERIES			
A59200, A59201	A59210	XE	INQUIRE
KONTRON/BIOTEK/NORTH STAR SCIENTIFIC			
440	LHA-34	Lamp Ass'y	
440	LH-34	D2 Lamp only	
720	L-205	D2	
810, 820, 860, 900, 930, 940	L-205	D2	
820, 860, 900, 930, 940	L-110	W	
UVIKON	L-205	D2	
UVIKON XS, XL			INQUIRE
ALL SPECTROPHOTOMETERS	L-205	D2	
KRATOS/SCHOEFFEL (ALSO SEE ABI /PERKIN-ELMER)			
FS950		HG	INQUIRE
FS970, 980	LH-98	D2	
SF757, SF759A, SF769, SF770, SF773	LH-98	D2	
SF783, SF783A, 785, 785A	LH-98	D2	
VISION	LH-40	D2	
SF769, SF770, FS970	LHA-1	Lamp Ass'y	
ALL OTHERS	LHA-3	Lamp Ass'y	
LABALLIANCE/SCIENTIFIC SYSTEMS (SSI)			
200, 201, 500, 504, 506, V520, 525	LH-40	D2	
200, 201, 500, 504, 506, V520, 525	LH-760	W	
1200/1500 SERIES			
AGG1314B	LH-14	D2	
AGG1315D	LH-13	D2	
AGG1321A	L-870	P/XE	
LC PACKINGS (SEE DIONEX)			
LDC (SEE THERMO SCIENTIFIC)			
LIFE (SEE ABI)			
LINEAR/OEM/PRIVATE LABEL			
(SEE ALL TECH, BAS, BIO-RAD, COLE-PALMER, DIONEX, FISHER, GBC, JSCO, LABALLIANCE			
LDC/THERMO, MICHROM BIORESOURCES, PERSEPTIVE BIOSYSTEMS, RAININ			
SPECTRA-PHYSICS, THERMO, THERMO SCIENTIFIC, VARIAN)			
LKB (SEE GE HEALTHCARE)			
MCPHERSON			
FL-750/A, FL-750/B	LH-101	XE	
FL-750/A, FL-750/B	LH-100	D2	
FL-750/A, FL-750/B	LH-102	XE-HG	