

# CHROMTECH LAMP GUIDE

January 2016

MERCK (SEE HITACHI)			
METERTECH (SEE BARNSTEAD/TURNER)			
METROHM		INQUIRE	
MICHROM BIORESOURCES/BRUKER			
UV-200	LH-40	D2	
MD1, 2, 3, 4	LH-25	D2	
MICROTECH SCIENTIFIC, INC.			
LINEAR MT200/200CE	LH-40	D2	
LINEAR MT201, 205	LH-40	D2	
MILTON-ROY (SEE THERMO SCIENTIFIC/BAUSCH&LOMB)			
MOLECULAR DEVICES			
ANALYST AD, HT		XE	INQUIRE
EMAX, EMAX PLUS		W	
FLEX STATION, I , II, II 384	L-870	P/XE	
FLEX STATION III	L-775	P/XE	
OPTIMAX	L-60	P/XE	
SPECTRAMAX 190, 340 PC, 384 PLUS	L-60	P/XE	
SPECTRAMAX 250	L-870	P/XE	
SPECTRAMAX 340	L-65	P/XE	
SPECTRAMAX M2/M2E,M3,M4,M5/M5E	L-775	P/XE	
SPECTRAMAX ANALYST GT	L-870	P/XE	
SPECTRAMAX GEMINI,EM,XPS,XS	L-870	P/XE	
SPECTRAMAX i3, PARADIGM			INQUIRE
SPECTRAMAX PLUS	L-60	P/XE	
THERMO MAX	L-120	W	\$80
UV MAX, VMAX		W	INQUIRE
NICOLET (SEE THERMO SCIENTIFIC/ATI)			
OCEAN OPTICS			
D1000			INQUIRE
DH2000			INQUIRE
PX-2		P/XE	INQUIRE
OLIS			
RSM-1000	L-1000	X/E	
PARAYTEC			
CLARITY 17, 620, 1000			INQUIRE
TDA200 (100D)	L-870		

PERKIN-ELMER/EG&G			
HPLC UV DETECTORS			
785A	LH-98	D2	\$395
INTEGRAL 4000	L-210	D2	\$295
LC-55, 55B, 65, 75, 85B, 95	L-210	D2	\$295
LC-135, 235	L-210	D2	\$295
LC-200	LH-98	D2	
LC-200		W	INQUIRE
LC-90, 250, 290, 295			INQUIRE
LC-480	LH-92	D2	
FLEXAR			
UV/VIS LC	LH-98	D2	
UHPLC,PDA, FLUORESCENCE			INQUIRE
200/DAD, PHASE I, PHASE II	LH-97	D2	
200/DAD	L-725	W	
200/FLUORO	L-800	XE	
200/UV-VIS/785A	LH-98	D2	
200/UV-VIS/785A		W	INQUIRE
200EP			INQUIRE
FLUOROMETERS			
200, 203, 204 SERIES	L-800	XE	
512, 650 SERIES	L-800	XE	
HTS 7000	INQUIRE	P/XE	
LS-3,4,5,30,40,50,50B,55		P/XE	INQUIRE
LC-240		P/XE	INQUIRE
MPF-2A, 3, 4	L-800	XE	
MPF-43, 43A, 44, 44A, 44B	L-800	XE	
SPECTROPHOTOMETERS			
46, 55, 57	L-200	D2	
55		W	INQUIRE
220,320,330,340	LH-51	D2	
3840	L-200	D2	
550, 551, 555	L-240	D2	
552A,553,559A,570,571,572,575,576	L-200	D2	
556, 557	LH-51	D2	
BIO	LH-92	D2	
BIO	LH-845	W	
EZ150, EZ201, EZ210	LH-51	D2	
EZ150, EZ201, EZ210	L-840	W	