

CHROMTECH LAMP GUIDE

January 2016

SHIMADZU (cont.)				SHIMADZU (cont.)			
HPLC UV DETECTORS (cont.)				SPECTROPHOTOMETERS (cont.)			
SPD-M1A			INQUIRE	UV-2501PC, 2550, 2600, 2700	L-265	D2	
SPD-M6A, SPD-M10A		D2	INQUIRE	UV-2501PC, 2550, 2600, 2700	L-110	W	
SPD-M6A, SPD-M10A	L-935	W		UV-3000			INQUIRE
SPD-M10AV		D2	INQUIRE	UV-3100, 3100PC, 3100PC	L-265	D2	
SPD-M10AV	L-245	W		UV-3101PC	L-110	W	
SPD-M10AVP	LH-88	D2	\$495	UV-3600, 3700	L-265	D2	
SPD-M10AVP		W	INQUIRE	UV-3600	L-110	W	
SPD-M20A, M20AV	LH-87	D2					
SPD-M20A	L-245	W		ATOMIC ABSORPTION SPECTROMETERS			
SPD-M20AV	L-830	W		AA-625,630,640,660,670,680U	L-265	D2	
SPD-M30A			INQUIRE	AA-6200	L-265	D2	
ELSD-LTII(2), ELSD-LT			INQUIRE	AA-6501, 6601, 6701, 6800			INQUIRE
FLUOROMETERS				AA-6650			INQUIRE
RF-10A, 10AXL	L-800	XE					
RF-10AVP			INQUIRE	SHODEX			
RF-20A, 20AXS	* L-795	XE		RI71,72,73,74			INQUIRE
RF-530, 535, 540, 551	L-800	XE		RI101			INQUIRE
RF-1501, 5000, 5301	L-800	XE		SICKMAIHAK			
SPECTROPHOTOMETERS				GM 31, 32	L-1005	D2	
BIOSPEC-MINI, BIOSPEC-NANO		P/XE	INQUIRE				
BIOSPEC 1601, MULTISPEC 1501	L-265	D2		SIEVERS/IONICS INSTR. BUSINESS GROUP/GE ANALYTICAL INST.			
BIOSPEC 1601, MULTISPEC 1501	L-110	W		400, 800	L-700	HG	
MULTISPEC 1501	L-935	W		900, 5310C, PPT	L-605	HG	
CS-920, 930, 9301	L-265	D2					
CS-930		W					
SERIES 160	L-110	W		SKALAR			
UV-120, 150, 160A, 160U, 163	L-265	D2		SA6711	LH-20	D2	
UV-120, 160A, 160U	L-110	W		SA6712	LH-22	D2	
UV-120-01, UV-120-02	L-265	D2		6750	LH-22	D2	
UV-120-01, UV-120-02		W		SONNTEK (SEE KNAUER)			
UV-140-02, UV-150-02	L-265	D2					
UV-240, 250, 260, 265	L-265	D2		SPECTRA-PHYSICS (SEE THERMO SCIENTIFIC)			
UV-240, 260, 265	L-110	W					
UV-1200, 1201, 1201PC, 1203, 1240	L-265	D2		SPECTRONICS (SEE THERMO SCIENTIFIC)			
UV-1200, 1201, 1203, 1240, 1240V	L-245	W					
UV-1501, 1600, 1601, 1601PC, 1603, 1650	L-265	D2		SPECTRONICS CORP			
UV-1501, 1600, 1601, 1601PC, 1603, 1650	L-110	W		TRANSILLUMINATORS			
UV-1700, 1800	L-265	D2		UV-B	INQUIRE	HG	
UV-1700, 1800	L-245	W					
UV-2100, 2101, 2401PC, 2450	L-265	D2		SPEX (SEE JY /JOBIN YVON)			
UV-2101PC, 2401, 2450	L-110	W					
* MUST RETURN SPENT LAMP & HOLDER				SPECTROVISION (SEE GROTON/DIONEX)			