

# CHROMTECH LAMP GUIDE

January 2016

STRATAGENE/AGILENT				THERMO SCIENTIFIC (cont.)			
MX 3000, 3005 P	INQUIRE	W		<b>BAUSCH &amp; LOMB (cont.)</b>			
				BIOMATE 5	LH-54	D2	
				BIOMATE 5	L-910	W	
SPECTRAFLOW 500	LH-40	D2		DE 50	L-225	D2	
SPECTRAFLOW 501	LH-25	D2		EPSILON	L-915	W	
SPECTRAFLOW 600 DAD			INQUIRE	EVOLUTION 60			INQUIRE
				EVOLUTION 6000	LH-54	D2	
				EVOLUTION 6000	L-910	W	
<b>SKYCAM/LASERCHROM</b>				GENESYS 2, 5		D2	INQUIRE
S3210 UV, S3210 DAD	INQUIRE	D2		GENESYS 2, 5	L-915	W	
S3210 UV, S3210 DAD	INQUIRE	W		GENESYS 6,10,10BIO,10UV,10UV-S		P/XE	INQUIRE
				GENESYS 8	LH-54	D2	
				GENESYS 8	L-910	W	
<b>TECAN</b>				GENESYS 10, 20	L-915	W	
GENIOS PLUS/GENIOUS PRO	L-1025	P/XE		SPECTRONIC 20, 20 PLUS, 20D, 20D PLUS	L-905	W	
INFINITE F200, M200 NANOQUANT	L-1025	P/XE		SPECTRONIC 21,21DV,21MV,21UVD,21UVM	L-945	W	
INFINITE F200, M200 PRO	L-1025	P/XE		SPECTRONIC 21D, 21DUV	L-225	D2	
INFINITE F500, M1000			INQUIRE	SPECTRONIC 70, 100	L-920	W	
SAFIRE II, ULTRA EVOLUTION	L-870			SPECTRONIC 200	L-225	D2	
SPECTRO FLUOR PLUS	L-1025	P/XE		SPECTRONIC 301, 401	L-915	W	
SUNRISE			INQUIRE	SPECTRONIC 501, 601		W	INQUIRE
ULTRA	L-870	P/XE	\$695	SPECTRONIC 601,700,701,710	L-225	D2	
<b>TEKMAR-DOHRMAN (SEE TELEDYNE)</b>				SPECTRONIC 700,701,710	L-920	W	
				SPECTRONIC 1001,1201,2000	L-225	D2	
				SPECTRONIC 1001, 1201	L-810	W	
<b>TELEDYNE/TEKMAR INSTRUMENTS</b>				SPECTRONIC 2000	L-920	W	
PHOENIX 8000			INQUIRE	SPECTRONIC 3000	L-270	D2	
TOX ROUND GRID LAMP			INQUIRE				
OZONE GENERATOR LAMP 7			INQUIRE	<b>DIONEX</b>			
				<b>DIONEX ICS SERIES</b>			
<b>THERMO SCIENTIFIC/DIONEX/LC PACKING/GYNKOTEK</b>				DAD, DAD-3000, MWD, MWD-3000	LH-45	D2	INQUIRE
<b>AMINCO</b>				VWD-064377,070220,070221	LH-45	D2	INQUIRE
AB-2 (SERIES 2)		XE	INQUIRE				
SLM 2000	L-280	D2		<b>DIONEX ULTIMATE 3000 SERIES/SYSTEMS</b>			
SLM 3000	L-270	D2		MWD-064654,VWD-064377	LH-45	D2	INQUIRE
SLM 4800, 8000, 8100	L-855	XE		DAD-3000(RS),MWD-3000(RS),VWD 3000	LH-45	D2	INQUIRE
SPF-500			INQUIRE	VWD3100 (RS),VWD 3400(RS)	LH-45	D2	INQUIRE
<b>ATI</b>				VDM-2	L-235	D2	
240, 250	LHA-34	D2 Lamp Ass'y		ULT-UVD, UVD-3000	LH-25	D2	
240, 250	LH-34	Lamp Only		ULT-UVD	LHA-23	W	
<b>BARNSTEAD/TURNER/METERTECH</b>				RF 2000	L-800	XE	
110, 111, 112	L-185	HG		<b>ENVIRONMENTAL INSTRUMENTS/SO2 ANALYZER</b>			
SP-850	L-835	W		43A, 43CT, 43CH, 43i(TLE)	L-40	P/XE	
SP-830+, SP-870+	L-190	W		580B	L-40	P/XE	
SP-890, SP 8001, M965			INQUIRE				
<b>BAUSCH &amp; LOMB/MILTON ROY</b>							
AQUAMATE	LH-54	D2					
AQUAMATE, AQUAMATE UV	L-910	W					
BIOMATE 3, 3-S		P/XE	INQUIRE				