

CHROMTECH LAMP GUIDE

January 2016

VARIAN/CARY/RAININ/AGILENT (cont.)				WATERS (cont.)			
CARY (cont.)							
100, 210, 219, 220, 300, 400, 500	L-230	D2		996	LH-33	D2	
100, 210, 219, 220, 300	L-930	W		2420, 2424 ELSD	L-120	W	
400	L-110	W		2475	LH-925	XE	
500	L-1015	W		2487, 2488	LH-36X2	D2	
2050	L-290	D2		2996	LH-33		
2200, 2300, 2400	L-205	D2		ACQUITY UPLC			
2550	L-295	D2		FLR DETECTOR		HG/XE	INQUIRE
ECLIPSE	L-1045	P/XE		PDA 2998, TUV 2489		D2	INQUIRE
9070, 9075	LH-910	P/XE		LC MODULE 1/34551 SYSTEM	LH-31	D2	
AA200,240,275,280,300,400,600,640	L-230	D2		* MUST RETURN SPENT LAMP & HOLDER			
AA775,800,875,975,1275,1475,3000	L-230	D2					
DMS70,80,90,100,200,300	L-230	D2		YOUNG LIN			
FL-1, 2	L-305	D2		YL 9120, 9160			INQUIRE
LC-5000, 5500	L-305	D2		ZEISS			
PRO STAR SERIES				MCS 501			INQUIRE
PDA-2	LH-82	D2					
PRO STAR, PRO STAR 320	LH-950	P/XE					
PRO STAR 310, 9050 **	INQUIRE	D2					
PRO STAR 325, 335	LH-325	D2					
PRO STAR 330	LH-56	D2					
PRO STAR 340, 345	LH-40	D2					
PRO STAR 340, 345	LH-760	W					
PRO STAR 360	LH-910	P/XE					
PRO STAR 363	L-800	XE					
SPECTRA A10,20,30,40,640,880	L-230	D2					
STAR 3050	L-305	D2					
UV-1, 1A	LH-950	P/XE					
UV50,100,200	L-305	D2					
UV-C,DII,M	LH-40	D2					
VISTA 5000,5500	L-305	D2					
** RE-USE HOLDER							
WATERS							
401, 403, 404	INQUIRE						
420	L-770	HG					
440, 441	L-820	HG					
450	LHA-1	D2 Lamp Ass'y					
	LH-1	Lamp Only					
470, 474	L-800	XE					
480, 481	L-205	D2					
484	LH-30	D2					
486	LH-31	D2					
490E	L-890	P/XE					
990, 991, 994	LHA-34	D2 Lamp Ass'y					
	LH-34	Lamp Only					

LEGEND

D2	DEUTERIUM
HG	MERCURY
L	LAMP ONLY
LH	LAMP & HOLDER
LHA	LAMP HOLDER ASSEMBLY
NA	SODIUM
P/XE	PULSED XENON
W	TUNGSTEN
XE	XENON
X2	EXTENDED WARRANTY
ZN	ZINC

Lamp types and prices subject to change without prior notice. Sonntek assumes no liability for representations made herein. All information is accurate to the best of our knowledge. Majority of light sources listed are genuine Sonntek parts made to meet or exceed original equipment manufacturer's specifications. Sonntek is not affiliated with any companies mentioned; company or product names may be trademarks of their respective owners. 01/01/2016.