

Restek Replacement Parts for Beckman HPLC Systems

Description	Model #	Similar to Beckman part #	qty.	cat.#	price
Inlet Check Valve Assembly	110A&B, 112, 114M, 116, 118, 125, 126, 127, 128 Pumps	240720	ea.	25439	
Outlet Check Valve Assembly	110A&B, 112, 114M, 116, 118, 125, 126, 127, 128 Pumps	240721	ea.	25440	
Inlet Check Valve Cartridge for Check Valve Assembly	110A&B, 112, 114M, 116, 118, 125, 126, 127, 128 Pumps	240620	ea.	25441	
Outlet Check Valve Cartridge for Check Valve Assembly	110A&B, 112, 114M, 116, 118, 125, 126, 127, 128 Pumps	240621	ea.	25442	
Guide Sleeve for Graphite Plunger Guide	110 Series Pumps	243713	ea.	25444	
Outlet Check Filter Frit Assembly	110 Series Pumps	240619	ea.	25445	
Piston Guide Assembly	110 Series and 112 Pumps	243045	ea.	25446	
Plunger	110 Series and 112 Pumps	243053	ea.	25447	
Plunger Seal	110 Series Pumps	887138	ea.	25448	
Sapphire Plunger	114M, 116, 118, 125, 126, 127, 128 Pumps	240714	ea.	25450	
Plunger Seal	114M, 116, 118, 126, 127, 128 Pumps	241037, 237162	ea.	25451	
Pump Seal, Gold**	114, 116, 125, 126, 127, 128 Pumps	241037	ea.	25452	
Plunger Wash Seal	125, 126, 127, 128 Pumps	728772, 728772	ea.	25453	
Deuterium Lamp	DU60, 62, 64, 65	596791	ea.	25454	

**Ultra-high molecular weight polyethylene (UHMWPE); increased resistance to abrasion, best for buffers.

Restek Replacement Parts for Hitachi HPLC Systems

Description	Model #	Similar to Hitachi part #	qty.	cat.#	price
SS Check Valve Cartridges (1/16" Ruby Ball)	655/6000/6200	AN0-085	2-pk.	25455	
Inlet Check Valve Assembly	655/6000/6200	AN0-0833, 885-1330	ea.	25456	
Outlet Check Valve Assembly	655/6000/6200	AN0-0834, 885-1331	ea.	25457	
Plunger Assembly	655, 6000, 6200, 7100	810-1033	ea.	25458	
Plunger Seal, Black*	L655A, 6000, 6200, 7100	655-1080	ea.	25459	
Plunger Seal, Gold**	655, 6000, 6200, 7100	655-1080	ea.	25460	
Inlet Check Valve Assembly	L-7100	AN0-0836	ea.	25461	
Outlet Check Valve Assembly	L-7100	AN0-0837	ea.	25462	
Rotor Assembly for Dilutor Valve	AS-7200 Autosampler AS-7200,	810-3095	ea.	25463	
Rotor Seal Kit	AS-7250 Injection Valve	AN0-0818	ea.	25464	
Deuterium Lamp, Prealigned	L4000, L4200, L4250, L7400	AN0-0424	ea.	25465	

*Graphite-filled Teflon; best for organic solvents.

**Ultra-high molecular weight polyethylene (UHMWPE); increased resistance to abrasion, best for buffers.

Restek Replacement Parts for PerkinElmer HPLC Systems

Description	Model #	Similar to PE part #	qty.	cat.#	price
Inlet Check Valve	PE 200, 4, 250, 400, 410, 620	2540177	ea.	25438	
Outlet Check Valve	PE 200, 4, 250, 400, 410, 620	2540197	ea.	25437	
Standard Pump Seal Kit	PE 200	N2910383	kit	26425	
High Pressure Seal Kit	PE Series 200 Pumps-1, 2, 3, 3B, 4, 10, 250, 400, 410, 620, Int. 4000 Pump	2540275	ea.	25433	
Standard High Pressure Seal, Gold**	PE Series 200 Pumps-1, 2, 3, 3B, 4, 10, 250, 400, 410, 620, Int. 4000 Pump	9907324	ea.	25434	
Low Pressure Seal, Gold**	PE Series 200 Pumps	9907339	ea.	25435	
Needle	PE Series 200 Pumps	N2930023	ea.	25432	
Deuterium Lamp	PE Lambda-2, 5, 7, 8, 10, 11, 12, 14, 15, 16, 17, 18, 19, 20, 40, 800, 900	B0160917	ea.	25436	
Deuterium Lamp	PE 200/785A	N2920149	ea.	25431	

**Ultra-high molecular weight polyethylene (UHMWPE); increased resistance to abrasion, best for buffers.

