

## Restek Replacement Parts for Shimadzu HPLC Systems

Description	Model #	Similar to Shimadzu part #	qty.	cat. #	price
Inlet Check Valve	LC-6A, LC-10AS	228-12353-91	ea.	25287	
Inlet Check Valve	LC-600, LC-9A, LC-10AD	228-18522-91	ea.	25295	
Inlet Check Valve	LC-10ADvp	228-39093-92	ea.	24984	
Inlet Check Valve	LC-10AT, LC-10ATvp	228-32166-91	ea.	26521	
Inlet Check Valve	LC-20AD/B	228-48249-91	ea.	26426	
Outlet Check Valve	LC-6A, LC-10AS	228-09054-93	ea.	25288	
Outlet Check Valve Rebuild Kit	LC-6A, LC-10AS	228-11200-91	2-pk.	25289	
Outlet Check Valve	LC-600, LC-9A, LC-10AD, LC-10AT	228-18522-92	ea.	25282	
Outlet Check Valve	LC-10ADvp, LC-10ATvp	228-34976-91	ea.	24983	
Outlet Check Valve	LC-20AD/B	228-45705-91	ea.	26427	
Plunger Seal, Gold**	LC-6A	228-11999-00	ea.	25285	
Plunger Seal, Polyethylene, Gold**	LC-10AS, LC-10AT, LC-10ATvp	228-21975-00	ea.	25290	
Plunger Seal	LC-600, LC-9A, LC-10AD	228-18745-00	ea.	25293	
Plunger Seal, Black*	LC-10ADvp, LC 20AD	228-35146-00	ea.	24980	
Plunger Seal, Gold**	LC-10ADvp	228-32628-00	ea.	24981	
Plunger Seal, Black*	SIL-10ADvp, LC-10ATvp	228-35145-00	ea.	24985	
Plunger Rinse Seal	LC-10AS	228-28499-00	ea.	25292	
Sapphire Plunger	LC-10AS	228-17019-93	ea.	25291	
Sapphire Plunger	LC-600, LC-9A, LC-10AD	228-18523-91	ea.	25294	
Needle Seal	SIL-10A, 10XL, 10ADvp	228-33355-04	ea.	25468	
Rotor Seal	SIL-10ADvp	228-21217-97	ea.	24986	
Rotor Seal Assembly	SIL-10A, 10AXL	228-21217-91	ea.	25469	
Stator Assembly	SIL-10A, 10AXL	228-21220-91	ea.	25470	
Syringe, 500µL	SIL-10A, 10AXL	228-25237-04	ea.	25471	
Plunger Assembly, Ceramic	LC-10ADvp	228-35601-91	ea.	25472	
Plunger Assembly, Ceramic	LC-10ATvp	228-35009-92	ea.	25473	
Plunger Assembly, Sapphire	LC-10ADvp	228-35601-92	ea.	24982	
Deuterium Lamp	SPD-10A, 10AV, SPD-20A	228-34016-02	ea.	25284	

\*Graphite-filled Teflon; best for organic solvents.

\*\*Ultra-high molecular weight polyethylene (UHMWPE); increased resistance to abrasion, best for buffers.



26521



## Catch the Buzz

Sign up for Restek's e-newsletter, *The Buzz*  
[www.restek.com/buzz](http://www.restek.com/buzz)

