



25474



25482



25484



25476

25280
(both sides shown)

25407

25408

25404

Restek Replacement Parts for Thermo Scientific HPLC Systems

Description	Model #	Similar to SP/TSP			
		part #	qty.	cat.#	price
Inlet Check Valve Assembly	SP8800 & P-Series Pumps	A3495-010	ea.	25474	
Outlet Check Valve Assembly	SP8800 Series Pumps	A3490-010	ea.	25475	
Piston	SP8800 & P-Series Pumps	A3102-010	ea.	25476	
Back-up Seal	SP8800 & P-Series Pumps	A2963-010	ea.	25477	
Plunger Seal, Gold**	SP8800 & P-Series Pumps	A2962-010	ea.	25478	
Superseal					
Check Valve and Transducer Assembly	P-Series Pumps	A3990-010	ea.	25479	
Kel-F Washer	P-Series Pumps	A2973-010	ea.	25480	
Rotor Seal Assembly, Rheodyne 7010	TSP AS100, 300, 1000, 3000, 3500, 8875, and 8880 Autosamplers	7010-039	ea.	25481	
Syringe Assembly, 250µL	TSP AS100, 300, 1000, 3000, 3500, 8875, and 8880 Autosamplers	A3588-020	ea.	25482	
Syringe, 500µL	8800 Series	3301-0100	ea.	25483	
Lamp, UV	Linear UV-200, 203, 204, 205, 206, and UV 100, 150, 1000, and 2000 Detectors	9551-0023	ea.	25484	

Description	Model #	Similar to TSP			
		part #	qty.	cat.#	price
Check Valve Cartridge	LDC Constametric Pumps	900946	ea.	25485	
Sapphire Plunger	LDC Constametric Pumps	801306	ea.	25486	
Plunger Seal Kit, Gold**	LDC Constametric Pumps	31-36-00754	kit	25487	
Plunger Seal, Black*	LDC Constametric Pumps	206129001	ea.	25488	
Plunger Seal, Gold**	LDC Constametric Pumps	206156001	ea.	25489	
Lamp, Deuterium	LDC SM-I, II, III, 3000, 3100, 3100X, and 4000 Detectors	108035	ea.	25490	

*Graphite-filled Teflon; best for organic solvents.

**Ultra-high molecular weight polyethylene (UHMWPE); increased resistance to abrasion, best for buffers.

Restek Rheodyne® Style Replacement Parts for Waters HPLC Systems

Description	Model #	Similar to Rheodyne			
		part #	qty.	cat.#	price
Vespel Rotor Seal	1090, 7000, 7010, 7040, 7067	7010-039	ea.	25279	
Vespel Rotor Seal	7125, 7126, 7725, 7725i, 9725	7125-047	ea.	25280	
Isolation Seal	7010	7010-015	ea.	25281	

For easy access to Rheodyne fittings, use a ValvTool Wrench. See page 337.

Restek Replacement Parts for Waters Detectors

Description	Model #	Similar to Waters			
		part #	qty.	cat.#	price
Window Gasket	484, 486, 490	80335	ea.	25423	
Lamp Side Gasket	484, 486, 490	80336	ea.	25424	
Quartz Lens	486	80687	ea.	25427	
Xenon Lamp (without holder or mirror)	470	—	ea.	25404	
Xenon Lamp	474	—	ea.	25405	
Deuterium Lamp (UV/Vis)	480, 481	99499	ea.	25403	
Deuterium Lamp (UV/Vis)	484	80357	ea.	25406	
Deuterium Lamp (UV/Vis)	486	80678	ea.	25407	
Deuterium Lamp	996, 2996	WAT052586	ea.	25408	
Deuterium Lamp	2487	WAS081142	ea.	25409	
Deuterium Lamp, long life (2,000 hours)*	486	80678	ea.	25410	

*Standard lamps have nominal 1,000-hour life.

also available

LC autosampler syringes

See page 349.

