

## Restek Replacement Parts for Waters HPLC Systems

Description	Model #	Similar to Waters part #	qty.	cat.#	price
<b>Preventive Maintenance Kit</b>					
Includes: sparge diffuser, filter insert, compression screws, SS ferrule, battery for 2690/717, 250µL WISP syringe, seal wash plunger seal kits (2), wash tube seal kits (4), 2690 seal pack rebuild kit (steel bodies not included), 2690 head plunger seal kits (2), solvent reservoir 20µm filters (4), Alliance check valve cartridges (2), Alliance plunger assemblies (2), 2690 face seals (4)					
	Alliance 2690, 2695	WAT270944	kit	25143	
<b>Preventive Maintenance Kit</b> Includes: PerformancePLUS cartridges (2), sparge diffusers (4), Super Seals (2), solvent reservoir 20µm filters (2), sapphire plungers (2), reference valve rebuild kit, inlet manifold kit					
	600 Pump	WAT052675	kit	25144	
<b>Preventive Maintenance Kit</b> Includes: 717 seal pack with needle, filter insert, 250µL WISP syringe					
	717 Autosampler	WAT052669	kit	25145	
<b>Preventive Maintenance Kit</b>					
Includes: 616/326/625 plunger assemblies (2), pump seal kit, sparge diffusers (4), solvent reservoir 20µm filters (4), 616 cartridge assemblies (4)					
	616 Pump	WAT052672	kit	25146	
<b>Preventive Maintenance Kit</b> Includes: PerformancePLUS check valves (2), sparge diffuser, solvent 10µm filter, sapphire plungers (2), plunger seals (2)					
	515 Pump	WAT052587	kit	26519	
<b>Preventive Maintenance Kit</b>					
Includes: plungers, check valves, seals, filters					
	1525 Pump	201000114	kit	26430	
Inlet Check Valve Assembly	M6KA, 501, 510, 515, 590, 600E	33679, 25214	ea.	25360	
Inlet Check Valve Housing	M6KA, 501, 510, 515, 590, 600E	25203	ea.	25361	
Inlet Check Valve Rebuild Kit	M6KA, 501, 510, 515, 590, 600E	60495	2-pk.	25362	
Outlet Check Valve Assembly (Actuator Style)	M6KA, 501, 510, 515, 590, 600E	WAT025028	ea.	25363	
Outlet Check Valve Housing (Actuator Style)	M6KA, 501, 510, 515, 590, 600E	25212	ea.	25364	
Outlet Check Valve Rebuild Kit (Actuator Style)	M6KA, 501, 510, 515, 590, 600E	26016	2-pk.	25365	
Outlet Check Valve Assembly (Ball & Seat Style)	M6KA, 501, 510, 515, 590, 600E	25216	ea.	25366	
Outlet Check Valve Housing (Ball & Seat Style)	M6KA, 501, 510, 515, 590, 600E	25207	ea.	25367	
Outlet Check Valve Rebuild Kit (Ball & Seat Style)	M6KA, 501, 510, 515, 590, 600E	26014	2-pk.	25368	
Inlet Check Valve Assembly, 225µL (Extended Flow)	M6KA, 501, 510, 515, 590, 600E	60307	ea.	25369	
PerformancePLUS Check Valve Cartridge	M6KA, 501, 510, 515, 590, 600E	700000254	2-pk.	25370	
Check Valve Rebuild Kit (Extended Flow)	M6KA, 501, 510, 515, 590, 600E	88223	2-pk.	25371	
PerformancePLUS Check Valve Housing	M6KA, 501, 510, 515, 590, 600E		ea.	25372	
Check Valve Cartridges	Alliance	WAT270941	2-pk.	25373	
Super Seal for Analytical Heads	M6KA, 501, 510, 515, 590, 600E	22946, 22934	ea.	25374	
Plunger Seal, Gold for Analytical Heads**	M6KA, 501, 510, 515, 590, 600E	22934	ea.	25375	
Plunger Seal, Black for graphite-filled Teflon*	M6KA, 501, 510, 515, 590, 600E	26613	ea.	25378	
Plunger Seal, Black for EF Heads*	510, 590, 600E	26644	ea.	25379	
Plunger Seal, Gold for EF Heads**	510, 590, 600E	26644	ea.	25380	
Seal Wash Plunger Seal	Alliance	WAT271018	2-pk.	25386	
Head Plunger Seal Kit	Alliance	WAT270938	2-pk.	25387	
Head Plunger Seal Kit (Black)	Alliance	WAT271066	2-pk.	25388	
Head Plunger Seal	Acquity, nanoAcquity	700002599	2-pk.	26428	
Insert Seal Parts Kit	M6KA, 501, 510, 515, 590, 600E	60012	kit	25389	
Sapphire Plunger	M6KA, 510, 590, 600	25656	ea.	25381	
Sapphire Plunger (Extended Flow)	510, 590, 600E	60304	ea.	25382	
Sapphire Plunger Assembly Kit	M515	WAS207069	ea.	25384	
Sapphire Plunger	616, 625, 626	31788	ea.	25420	
Sapphire Plunger	Alliance	WAT270959	ea.	25385	
Draw-Off Tube Assembly	510, 515, 600, 610	WAT060476	ea.	25392	
Inlet Manifold Kit	M45, 501, 510, 590, 600E	60448	kit	25412	
Ferrule, Stainless Steel	515	22330	ea.	25417	
Face Seals Replacement Kit	Alliance 2690	WAT270939	4-pk.	26547	
Wash Face Seal	Alliance 2690	WAT271017	ea.	25428	
Wash Tube Seal Kit	Alliance 2690	WAT270940	4-pk.	25429	
Wash Seal	Acquity, nanoAcquity	700002598	2-pk.	26429	
Rebuild Kit, Reference Valve	600 Pump	25746	ea.	25492	
Compression Screws, Stainless Steel	Alliance 2690		2-pk.	25493	
Battery	Alliance 2690, 717, 486, 484		ea.	25494	
Seal Pack Rebuild Kit, without Seal Wash Tube Assembly	Alliance 2690	WAT271019	ea.	25495	
717 Seal Pack (Rebuilt) with Needle	717 Autosampler	4559	ea.	25496	
Seal Kit, Pump	616 Pump, 626, 625	31790	ea.	25497	
Cartridge Assembly, Stainless Steel	616 Pump	24960	ea.	25498	

\*Graphite-filled Teflon; best for organic solvents.;

\*\*Ultra-high molecular weight polyethylene (UHMWPE); increased resistance to abrasion, best for buffers.



Kits are convenient and economical, compared to ordering individual items!

