

QUICK  
CONNECTS

# SBLSeries

## Quick Connects

The SBL Series connector system is so intuitive that users just automatically know how to operate it the minute they pick it up. SBL Series quick connect fittings are cost effective and designed to meet your needs for a reusable component for small-bore applications in all fluid delivery areas. The SBL Series range is available in nylon, in two contrasting colors and fits 1/8" (3.2 mm) and 3/16" (4.8 mm) ID tubing for in-line applications.

**SBL SERIES:**  
Female Open Flow

**SBLF013 (-1 -2 -8B -CM030)**  
Side Latch Quick Connect Female 500 Series Barb, 1/8" (3.2 mm) Tubing ID



**SBLF035 (-1 -2 -8B -CM030)**  
Side Latch Quick Connect Female 500 Series Barb, 3/16" (4.8 mm) Tubing ID



**SBL SERIES:**  
Male Open Flow

**SBLM013 (-001 -002 -004 -005) Nylon, Buna-N O-Ring**  
Side Latch Quick Connect Male 500 Series Barb, 1/8" (3.2 mm) Tubing ID



-001 Gray Nylon  
-002 Blue Nylon

**SBLM035 (-001 -002 -004 -005) Nylon, Buna-N O-Ring**  
Side Latch Quick Connect Male 500 Series Barb, 3/16" (4.8 mm) Tubing ID



-001 Gray Nylon  
-002 Blue Nylon

**SBLMFGB013 (-001)**  
Bag of 50 Side Latch Quick Connect Male and Female Sets, 500 Series Barb 1/8" (3.2 mm) Tubing ID, Grey Male and Blue Female Nylon



**SBLMFGB013 (-002)**  
Bag of 50 Side Latch Quick Connect Male and Female Sets, 500 Series Barb 1/8" (3.2 mm) Tubing ID, Blue Male and Grey Female Nylon



**SBLMFGB035 (-003)**  
Bag of 50 Side Latch Quick Connect Male and Female Sets, 500 Series Barb 3/16" (4.8 mm) Tubing ID, Grey Male and Blue Female Nylon



**SBLMFGB035 (-004)**  
Bag of 50 Side Latch Quick Connect Male and Female Sets, 500 Series Barb 3/16" (4.8 mm) Tubing ID, Blue Male and Grey Female Nylon



### Need it Fast?

You can order any of our 4,000+ parts online now with no order \$ minimums, 24 hours a day, 7 days a week!  
[www.nordsonmedical.com](http://www.nordsonmedical.com)