


Ordering Information

Screw Top, Standard Opening Vials		
Part No.	Description	Unit
LBSV012C	2ml Clear vial, 8-425 screw top	100/pk
LBSV032C	2ml Clear vial, 8-425 screw top, graduated with writing area	100/pk
LBSV012A	2ml Amber vial, 8-425 screw top	100/pk
LBSV032A	2ml Amber vial, 8-425 screw top, graduated with writing area	100/pk
Screw caps - 8-425, polypropylene plastic for standard screw top vial		
LBSV012	Black screw cap with hole, for 2ml 8-425 screw top vial	100/pk
LBSV032	Blue screw cap with hole, for 2mL 8-425 screw top vial	100/pk
Septa for 8-425 Screw caps		
LBSV02RS	White PTFE/red silicone septa, for 2ml 8-425 screw top vial	10000/pk
LBSV02SS	Red PTFE/white silicone septa, for 2ml 8-425 screw top vial	10000/pk
Screw caps + Septa for standard screw top vial		
LBSV012CRS	White PTFE/red silicone septa + Black screw cap with hole, for 2ml 8-425 screw top vial	100/pk
LBSV012CSS	Red PTFE/white silicone septa + Black screw cap with hole, for 2ml 8-425 screw top vial	100/pk



Screw Top, Wide Opening Vials		
Part No.	Description	Unit
LBSV002C	2ml Clear vial, 9-425 screw top	100/pk
LBSV022C	2ml Clear vial, 9-425 screw top, graduated with writing area	100/pk
LBSV002A	2ml Amber vial, 9-425 screw top	100/pk
LBSV022A	2ml Amber vial, 9-425 screw top, graduated with writing area	100/pk
Screw caps - 9-425, polypropylene plastic for Wide Opening screw top vial		
LBSV002	Blue screw cap with hole, for 2ml 9-425 screw top vial	100/pk
Septa for 9-425 Screw caps		
LBSV03RS	White PTFE/red silicone septa, for 2ml 9-425 screw top vial	10000/pk
LBSV03SS	Red PTFE/white silicone septa, for 2ml 9-425 screw top vial	10000/pk
LBSV23SS	Blue PTFE/white silicone septa, Pre-slit, for 2mL 9-425 screw top vial	10000/pk
Screw caps + Septa for Wide Opening screw top vial		
LBSV002CRS	White PTFE/red silicone septa + Blue screw cap with hole, for 2ml 9-425 screw top vial	100/pk
LBSV102CSS	Red PTFE/white silicone septa + Blue screw cap with hole, for 2ml 9-425 screw top vial	100/pk
LBSV222CSS	Blue PTFE/white silicone septa, Pre-slit + Blue screw cap with hole, for 2ml 9-425 screw top vial	100/pk

