

Technical Parameter:

Parameters	CA (Sterile)			PES (Sterile)		
	13mm	25mm	33mm	13mm	25mm	33mm
Membrane material	CA	CA	CA	PES	PES	PES
Housing material	PP	PES	PP	PP	PES	PP
Filtration area (cm ²)	0.65	3.90	4.60	0.65	3.90	4.60
Pore Size(μm)	0.22 / 0.45	0.22 / 0.45	0.22 / 0.45	0.22 / 0.45	0.22 / 0.45	0.22 / 0.45
Holdup volume (μl)	<10	<30	<55	<10	<30	<55
Sample volume (ml)	<12	<100	<140	<12	<100	<140
Maximum Operating Temperature	110°C	110°C	110°C	90°C	90°C	90°C
Maximum Operating Pressure (psi)	50	95	120	50	95	120
Applicable pH value	4-8	4-8	4-8	1-14	1-14	1-14

Order Information:

PES Sterile Syringe Filter, Gama Sterile, Individually Package		
Item NO.	Discription	package
SFPES013022S	Sterile PES Syring Filters, 0.22(μm), 13(mm)	100
SFPES013045S	Sterile PES Syring Filters, 0.45(μm), 13(mm)	100
SFPES025022S	Sterile PES Syring Filters, 0.22(μm), 25(mm)	100
SFPES025045S	Sterile PES Syring Filters, 0.45(μm), 25(mm)	100
SFPES033022S	Sterile PES Syring Filters, 0.22(μm), 33(mm)	100
SFPES033045S	Sterile PES Syring Filters, 0.45(μm), 33(mm)	100

CA Sterile Syringe Filter, Gama Sterile, Individually Package		
Item NO.	Discription	package
SFCA013022S	Sterile CA Syringe Filters, 0.22(μm), 13(mm)	100
SFCA013045S	Sterile CA Syringe Filters, 0.45(μm), 13(mm)	100
SFCA025022S	Sterile CA Syringe Filters, 0.22(μm), 25(mm)	100
SFCA025045S	Sterile CA Syringe Filters, 0.45(μm), 25(mm)	100
SFCA033022S	Sterile CA Syringe Filters, 0.22(μm), 33(mm)	100
SFCA033045S	Sterile CA Syringe Filters, 0.45(μm), 33(mm)	100