

MS SimplePure Syringe Filter

Item#	Description	/Pkt
Nylon66 Syringe Filter (Polyamide-66, 13mm/25mm)		
SFNY013022N	Nonsterile Nylon Syringe Filters, Pore:0.22(μm), Housing Diameter:13(mm)	100
SFNY013045N	Nonsterile Nylon Syringe Filters, Pore:0.45(μm), Housing Diameter:13(mm)	100
SFNY025022N	Nonsterile Nylon Syringe Filters, Pore:0.22(μm), Housing Diameter:25(mm)	100
SFNY025045N	Nonsterile Nylon Syringe Filters, Pore:0.45(μm), Housing Diameter:25(mm)	100
PTFE Syringe Filter (Polytetrafluoroethylene, Hydrophobic, 13mm/25mm)		
SFPTFE013022NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.22(μm), Housing Diameter:13(mm)	100
SFPTFE013045NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.45(μm), Housing Diameter:13(mm)	100
SFPTFE025022NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.22(μm), Housing Diameter:25(mm)	100
SFPTFE025045NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.45(μm), Housing Diameter:25(mm)	100
PTFE Syringe Filter (Polytetrafluoroethylene, Hydrophilic, 13mm/25mm)		
SFPTFE013022NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.22(μm), Housing Diameter:13(mm)	100
SFPTFE013045NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.45(μm), Housing Diameter:13(mm)	100
SFPTFE025022NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.22(μm), Housing Diameter:25(mm)	100
SFPTFE025045NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.45(μm), Housing Diameter:25(mm)	100
PES Syringe Filter (13mm/25mm)		
SFPES013022N	Nonsterile PES Syringe Filters, Pore:0.22(μm), Housing Diameter:13(mm)	100
SFPES013045N	Nonsterile PES Syringe Filters, Pore:0.45(μm), Housing Diameter:13(mm)	100
SFPES025022N	Nonsterile PES Syringe Filters, Pore:0.22(μm), Housing Diameter:25(mm)	100
SFPES025045N	Nonsterile PES Syringe Filters, Pore:0.45(μm), Housing Diameter:25(mm)	100
PVDF Syringe Filter (Polyvinylidene Fluoride, 13mm/25mm)		
SFPVDF013022N	Nonsterile Hydrophobic PVDF Syringe Filters, Pore:0.22(μm), Housing Diameter:13(mm)	100
SFPVDF013045N	Nonsterile Hydrophobic PVDF Syringe Filters, Pore:0.45(μm), Housing Diameter:13(mm)	100
SFPVDF025022N	Nonsterile Hydrophobic PVDF Syringe Filters, Pore:0.22(μm), Housing Diameter:25(mm)	100
SFPVDF025045N	Nonsterile Hydrophobic PVDF Syringe Filters, Pore:0.45(μm), Housing Diameter:25(mm)	100
GF Syringe Filter (Glass Fiber, 13mm/25mm)		
SFGF013070N	Nonsterile GF Syringe Filters, Pore:0.7(μm), Housing Diameter:13(mm)	100
SFGF013100N	Nonsterile GF Syringe Filters, Pore:1.0(μm), Housing Diameter:13(mm)	100
SFGF025070N	Nonsterile GF Syringe Filters, Pore:0.7(μm), Housing Diameter:25(mm)	100
SFGF025100N	Nonsterile GF Syringe Filters, Pore:1.0(μm), Housing Diameter:25(mm)	100
CA Syringe Filter (Cellulose Acetate, 13mm/25mm)		
SFCA013022N	Nonsterile CA Syringe Filters, Pore:0.22(μm), Housing Diameter:13(mm)	100
SFCA013045N	Nonsterile CA Syringe Filters, Pore:0.45(μm), Housing Diameter:13(mm)	100
SFCA025022N	Nonsterile CA Syringe Filters, Pore:0.22(μm), Housing Diameter:25(mm)	100
SFCA025045N	Nonsterile CA Syringe Filters, Pore:0.45(μm), Housing Diameter:25(mm)	100
MCE Syringe Filter (Mixed Cellulose Ester, 13mm/25mm)		
SFMCE013022N	Nonsterile MCE Syringe Filters, Pore:0.22(μm), Housing Diameter:13(mm)	100
SFMCE013045N	Nonsterile MCE Syringe Filters, Pore:0.45(μm), Housing Diameter:13(mm)	100
SFMCE025022N	Nonsterile MCE Syringe Filters, Pore:0.22(μm), Housing Diameter:25(mm)	100
SFMCE025045N	Nonsterile MCE Syringe Filters, Pore:0.45(μm), Housing Diameter:25(mm)	100
PP Syringe Filter (Polypropylene, 13mm/25mm)		
SFPP013022N	Nonsterile PP Syringe Filters, Pore:0.22(μm), Housing Diameter:13(mm)	100
SFPP013045N	Nonsterile PP Syringe Filters, Pore:0.45(μm), Housing Diameter:13(mm)	100
SFPP025022N	Nonsterile PP Syringe Filters, Pore:0.22(μm), Housing Diameter:25(mm)	100
SFPP025045N	Nonsterile PP Syringe Filters, Pore:0.45(μm), Housing Diameter:25(mm)	100

MS Syringe Filters Product List 2014

PTFE Syringe Filters : ALL Sizes

SFPTFE013022NG	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFPTFE013045NG	Nonsterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFPTFE025022NG	Nonsterile PTFE Syringe Filters,Pore: 0.22(µm), Diameter:25(mm), GF Prefilter	100
SFPTFE025045NG	Nonsterile PTFE Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), GF Prefilter	100
SFPTFE030045NG	Nonsterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:30(mm), GF Prefilter	100
SFPTFE050022NG	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:50(mm), GF Prefilter	12
SFPTFE050045NG	Nonsterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:50(mm), GF prefilter	12
SFPTFE013022SG	Sterile PTFE Syringe Filters, Pore:0.22(µm),Diameter: 13(mm), GF Prefilter	100
SFPTFE013045SG	Sterile PTFE Syringe Filters,Pore: 0.45(µm), Diameter:13(mm), GF Prefilter	100
SFPTFE025022SG	Sterile PTFE Syringe Filters, Pore:0.22(µm),Diameter: 25(mm), GF Prefilter	100
SFPTFE025045SG	Sterile PTFE Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), GF Prefilter	100
SFPTFE013022NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.22(µm), Diameter:13(mm)	500
SFPTFE013045NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.45(µm), Diameter:13(mm)	500
SFPTFE025010NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.22(µm),Diameter:25(mm)	200
SFPTFE025022NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.22(µm),Diameter:25(mm)	200
SFPTFE025045NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.45(µm), Diameter:25(mm)	200
SFPTFE025500NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:5.0(µm),Diameter:25(mm)	200
SFPTFE033022NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.22(µm), Diameter:33(mm)	200
SFPTFE033045NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.45(µm), Diameter:33(mm)	200
SFPTFE050022NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.22(µm), Diameter:50(mm)	12
SFPTFE050045NL	Nonsterile Hydrophilic PTFE Syringe Filters, Pore:0.45(µm), Diameter:50(mm)	12
SFPTFE013010NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.1(µm),Diameter:13(mm)	500
SFPTFE013022NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.22(µm), Diameter:13(mm)	500
SFPTFE013045NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.45(µm),Diameter:13(mm)	500
SFPTFE013100NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:1.0(µm), Diameter:13(mm)	500
SFPTFE013500NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:5.0(µm),Diameter:13(mm)	500
SFPTFE017022NB	Nonsterile Hydrophobic PTFE Syringe Filter, Pore:0.22(µm), Diameter:17(mm)	500
SFPTFE017045NB	Non-sterile Hydrophobic PTFE Syringe Filter, Pore:0.45(µm), Diameter:17(mm)	500
SFPTFE025010NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.1(µm), Diameter:25(mm)	200
SFPTFE025022NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.22(µm), Diameter:25(mm)	200
SFPTFE025045NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.45(µm),Diameter:25(mm)	200
SFPTFE025100NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:1.0(µm),Diameter:25(mm)	200
SFPTFE025200NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:2.0(µm),Diameter:25(mm)	200
SFPTFE025300NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:3.0(µm),Diameter:25(mm)	200
SFPTFE025500NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:5.0(µm),Diameter:25(mm)	200
SFPTFE030022NB	Nonsterile Hydrophobic PTFE Syringe Filter, Pore:0.22(µm), Diameter:30(mm)	200
SFPTFE030045NB	Nonsterile Hydrophobic PTFE Syringe Filter, Pore:0.45(µm), Diameter:30(mm)	200
SFPTFE030100NB	Nonsterile Hydrophobic PTFE Syringe Filter, Pore:0.45(µm), Diameter:30(mm)	200
SFPTFE033022NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.22(µm),Diameter:33(mm)	200
SFPTFE033045NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.45(µm),Diameter:33(mm)	200
SFPTFE033100NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:1.0(µm),Diameter:33(mm)	200
SFPTFE050022NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.22(µm),Diameter:50(mm)	12
SFPTFE050045NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:0.45(µm), Diameter:50(mm)	12
SFPTFE050100NB	Nonsterile Hydrophobic PTFE Syringe Filters, Pore:1.0(µm), Diameter:50(mm)	12
SFPTFE013022SL	Sterile Hydrophilic PTFE Syringe Filters, Pore:0.22(µm), Diameter:13(mm)	200
SFPTFE013045SL	Sterile Hydrophilic PTFE Syringe Filters, Pore:0.45(µm), Diameter:13(mm)	200
SFPTFE025022SL	Sterile Hydrophilic PTFE Syringe Filters, Pore:0.22(µm), Diameter:25(mm)	200
SFPTFE025045SL	Sterile Hydrophilic PTFE Syringe Filters, Pore:0.45(µm), Diameter:25(mm)	200
SFPTFE033022SL	Sterile Hydrophilic PTFE Syringe Filters, Pore:0.22(µm),Diameter:33(mm)	200
SFPTFE033045SL	Sterile Hydrophilic PTFE Syringe Filters, Pore:0.45(µm),Diameter:33(mm)	200

MS Syringe Filters Product List 2014

SFPTFE013022SB	Sterile Hydrophobic PTFE Syringe Filters, Pore:0.22(µm), Diameter:13(mm)	
SFPTFE013045SB	Sterile Hydrophobic PTFE Syringe Filters, Pore:0.45(µm), Diameter:13(mm)	200
SFPTFE025022SB	Sterile Hydrophobic PTFE Syringe Filters, Pore:0.22(µm),Diameter:25(mm)	200
SFPTFE025045SB	Sterile Hydrophobic PTFE Syringe Filters, Pore:0.45(µm), Diameter:25(mm)	200
SFPTFE033022SB	Sterile Hydrophobic PTFE Syringe Filters, Pore:0.22(µm), Diameter:33(mm)	200
SFPTFE033045SB	Sterile Hydrophobic PTFE Syringe Filters, Pore:0.45(µm),Diameter:33(mm)	200
SFPTFE050022SB	Sterile Hydrophobic PTFE Syringe Filters, Pore:0.22(µm), Diameter:50(mm)	10
SFPTFE013022NP	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPTFE013045NP	Nonsterile PTFE Syringe Filters,Pore: 0.45(µm), Diameter:13(mm), PP Prefilter	100
SFPTFE025022NP	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPTFE025045NP	Nonsterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP Prefilter	100
SFPTFE050022NP	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFPTFE050045NP	Nonsterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:50(mm), PP Prefilter	12
SFPTFE013022SP	Sterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPTFE013045SP	Sterile PTFE Syringe Filters, Pore:0.45(µm),Diameter: 13(mm), PP prefilter	100
SFPTFE025022SP	Sterile PTFE Syringe Filters,Pore: 0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPTFE025045SP	Sterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP Prefilter	100
PVDF Syringe Filters :: ALL Sizes		
SFPVDF013022NG	Nonsterile PVDF Syringe Filter, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFPVDF013045NG	Nonsterile PVDF Syringe Filter,Pore: 0.45(µm),Diameter: 13(mm), GF Prefilter	100
SFPVDF025022NG	Nonsterile PVDF Syringe Filter, Pore:0.22(µm),Diameter: 25(mm), GF Prefilter	100
SFPVDF025045NG	Nonsterile PVDF Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFPVDF050022NG	Nonsterile PVDF Syringe Filter, Pore:0.22(µm), Diameter:50(mm), GF Prefilter	12
SFPVDF050045NG	Nonsterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:50(mm), GF Prefilter	12
SFPVDF013022SG	Sterile PVDF Syringe Filter,Pore: 0.22(µm), Diameter:13(mm), GF Prefilter	100
SFPVDF013045SG	Sterile PVDF Syringe Filters, Pore:0.45(µm),Diameter: 13(mm), GF Prefilter	100
SFPVDF025022SG	Sterile PVDF Syringe Filter,Pore: 0.22(µm),Diameter: 25(mm), GF Prefilter	100
SFPVDF025045SG	Sterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFPVDF013010N	Nonsterile PVDF Syringe Filter, Pore:0.1(µm),Diameter:13(mm)	500
SFPVDF013022N	Nonsterile PVDF Syringe Filter, Pore:0.22(µm), Diameter:13(mm)	500
SFPVDF013045N	Nonsterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:13(mm)	500
SFPVDF017022N	Nonsterile PVDF Syringe Filter, Pore:0.22(µm), Diameter:17(mm)	500
SFPVDF017045N	Nonsterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:17(mm)	500
SFPVDF025010N	Nonsterile PVDF Syringe Filter, Pore:0.1(µm), Diameter:25(mm)	200
SFPVDF025022N	Nonsterile PVDF Syringe Filter, Pore:0.22(µm),Diameter:25(mm)	200
SFPVDF025045N	Nonsterile PVDF Syringe Filter, Pore:0.45(µm),Diameter:25(mm)	200
SFPVDF025200N	Nonsterile PVDF Syringe Filter, Pore:2.0(µm),Diameter:25(mm)	200
SFPVDF025300N	Nonsterile PVDF Syringe Filter, Pore:3.0(µm),Diameter:25(mm)	200
SFPVDF030010N	Nonsterile PVDF Syringe Filter, Pore:0.1(µm), Diameter:30(mm)	100
SFPVDF030022N	Nonsterile PVDF Syringe Filter, Pore:0.22(µm), Diameter:30(mm)	200
SFPVDF030045N	Nonsterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:30(mm)	200
SFPVDF033022N	Nonsterile PVDF Syringe Filter, Pore:0.22(µm),Diameter:33(mm)	200
SFPVDF033045N	Nonsterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:33(mm)	200
SFPVDF050022N	Nonsterile PVDF Syringe Filter, Pore:0.22(µm),Diameter:50(mm)	12
SFPVDF050045N	Nonsterile PVDF Syringe Filter, Pore:0.45(µm),Diameter:50(mm)	12
SFPVDF013022S	Sterile PVDF Syringe Filter, Pore:0.22(µm),Diameter:13(mm)	200
SFPVDF013045S	Sterile PVDF Syringe Filter, Pore:0.45(µm),Diameter:13(mm)	200
SFPVDF025022S	Sterile PVDF Syringe Filter, Pore:0.22(µm),Diameter:25(mm)	200
SFPVDF025045S	Sterile PVDF Syringe Filter, Pore:0.45(µm),Diameter:25(mm)	200
SFPVDF033022S	Sterile PVDF Syringe Filter, Pore:0.22(µm),Diameter:33(mm)	200

MS Syringe Filters Product List 2014

SFPVDF033045S	Sterile PVDF Syringe Filter, Pore:0.45(µm),Diameter:33(mm)	
SFPVDF013022NP	Nonsterile PVDF Syringe Filter,Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPVDF013045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFPVDF025022NP	Nonsterile PVDF Syringe Filter,Pore: 0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPVDF025045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(µm),Diameter: 25(mm), PP Prefilter	100
SFPVDF050022NP	Nonsterile PVDF Syringe Filter, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFPVDF050045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:50(mm), PP Prefilter	12
SFPVDF013022SP	Sterile PVDF Syringe Filter,Pore: 0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPVDF013045SP	Sterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFPVDF025022SP	Sterile PVDF Syringe Filter,Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPVDF025045SP	Sterile PVDF Syringe Filter,Pore: 0.45(µm), Diameter:25(mm), PP Prefilter	100
PES Syringe Filters : : ALL Sizes		
SFPES013002NG	Nonsterile PES Syring Filters,Pore: 0.02(um), Diamter:13(mm), GF Prefilter	100
SFPES013022NG	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF prefilter	100
SFPES013045NG	Nonsterile PES Syringe Filters, Pore:0.45(µm),Diameter: 13(mm), GF Prefilter	100
SFPES025002NG	Nonsterile PES Syringe Filters, Pore:0.02(µm), Diameter:25(mm), GF Prefilter	100
SFPES025022NG	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm), GF Prefilter	100
SFPES025045NG	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFPES050022NG	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:50(mm), GF Prefilter	12
SFPES050045NG	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:50(mm), GF prefilter	12
SFPES013022SG	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFPES013045SG	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFPES025022SG	Sterile PES Syringe Filters, Pore:0.02(µm), Diameter:25(mm), GF Prefilter	100
SFPES025045SG	Sterile PES Syringe Filters, Pore:0.02(µm), Diameter:25(mm), GF Prefilter	100
SFPES013010N	Nonsterile PES Syringe Filters, Pore:0.1(µm), Diameter:13(mm)	500
SFPES013022N	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm)	500
SFPES013045N	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:13(mm)	500
SFPES017022N	Nonsterile PES Syringe Filter, Pore:0.22(µm), Diameter:17(mm)	500
SFPES017045N	Nonsterile PES Syringe Filter, Pore:0.45(µm), Diameter:17(mm)	500
SFPES025010N	Nonsterile PES Syringe Filters, Pore:0.1(µm), Diameter:25(mm)	200
SFPES025022N	Nonsterile PES Syringe Filters, Pore:0.22(µm),Diameter:25(mm)	200
SFPES025045N	Nonsterile PES Syringe Filters, Pore:0.45(µm),Diameter:25(mm)	200
SFPES025200N	Nonsterile PES Syringe Filters, Pore:2.0(µm), Diameter:25(mm)	200
SFPES030022N	Nonsterile PES Syringe Filter, Pore:0.22(µm), Diameter:30(mm)	200
SFPES030045N	Nonsterile PES Syringe Filter, Pore:0.45(µm), Diameter:30(mm)	200
SFPES030200N	Nonsterile PES Syringe Filter, Pore:2.0(µm), Diameter:30(mm)	200
SFPES033022N	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:33(mm)	200
SFPES033045N	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:33(mm)	200
SFPES050022N	Nonsterile PES Syringe Filters, Pore:0.22(µm),Diameter:50(mm)	12
SFPES050045N	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:50(mm)	12
SFPES013022S	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm)	200
SFPES013045S	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:13(mm)	200
SFPES025010S	Sterile PES Syringe Filters, Pore:0.1(µm), Diameter:25(mm)	100
SFPES025022S	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm)	200
SFPES025045S	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm)	100
SFPES030022S	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:30(mm)	100
SFPES033022S	Sterile PES Syringe Filters, Pore:0.22(µm),Diameter:33(mm)	200
SFPES033045S	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:33(mm)	200
SFPES013022NP	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100

MS Syringe Filters Product List 2014

SFPES013045NP	Nonsterile PES Syringe Filters,Pore: 0.45(µm), Diameter:13(mm), PP Prefilter	
SFPES025022NP	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPES025045NP	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP prefilter	100
SFPES050022NP	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFPES050045NP	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:50(mm), PP Prefilter	12
SFPES013022SP	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPES013045SP	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP prefilter	100
SFPES025022SP	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP prefilter	100
SFPES025045SP	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP prefilter	100
Nylon Syringe Filters : : ALL Sizes		
SFNY013010NG	Nonsterile Nylon Syringe Filters, Pore:0.10(µm), Diameter:13(mm), GF Prefilter	100
SFNY013022NG	Nonsterile Nylon Syringe Filters,Pore: 0.22(µm), Diameter: 13(mm), GF Prefilter	100
SFNY013045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFNY025022NG	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm), GF Prefilter	100
SFNY025045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFNY025500NG	Nonsterile Nylon Syringe Filters, Pore:5.0(µm), Diameter:25(mm), GF Prefilter	100
SFNY030045NG	Nonsterile Nylon Syringe Filter, Pore:0.45(µm), Diameter:30(mm),GF Prefilter	100
SFNY050022NG	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:50(mm), GF Prefilter	12
SFNY050045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:50(mm), GF Prefilter	12
SFNY013022SG	Sterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFNY013045SG	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFNY025022SG	Sterile Nylon Syringe Filters, Pore:0.22(µm),Diameter: 25(mm), GF Prefilter	100
SFNY025045SG	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SF1322N1000	Non-sterile Nylon Syringe Filter, Pore:0.22(µm),Diameter:13(mm)	1000
SF1345N1000	Non-sterile Nylon Syringe Filter, Pore:0.45(µm),Diameter:13(mm)	1000
SFNY004020N	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:4(mm)	100
SFNY013010N	Nonsterile Nylon Syringe Filters, Pore:0.1(µm), Diameter:13(mm)	500
SFNY013022N	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:13(mm)	500
SFNY013045N	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm)	500
SFNY013080N	Nonsterile Nylon Syringe Filters, Pore:0.8(µm), Diameter:13(mm)	500
SFNY013100N	Nonsterile Nylon Syringe Filters, Pore:1.0(µm), Diameter:13(mm)	500
SFNY013300N	Nonsterile Nylon Syringe Filters, Pore:3.0(µm), Diameter:13(mm)	500
SFNY013500N	Nonsterile Nylon Syringe Filters, Pore:5.0(µm), Diameter:13(mm)	500
SFNY017022N	Nonsterile Nylon Syringe Filter, Pore:0.22(µm), Diameter:17(mm)	500
SFNY017045N	Nonsterile Nylon Syringe Filter, Pore:0.45(µm), Diameter:17(mm)	500
SFNY025022N	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm)	200
SFNY025045N	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:25(mm)	200
SFNY025080N	Nonsterile Nylon Syringe Filters, Pore:0.8(µm), Diameter:25(mm)	200
SFNY025100N	Nonsterile Nylon Syringe Filters, Pore:1.0(µm), Diameter:25(mm)	200
SFNY025200N	Non-sterile Nylon Syringe Filters, Pore:2.0(µm), Diameter:25(mm)	200
SFNY025300N	Nonsterile Nylon Syringe Filters, Pore:3.0(µm), Diameter:25(mm)	200
SFNY025500N	Nonsterile Nylon Syringe Filters, Pore:5.0(µm), Diameter:25(mm)	200
SFNY030022N	Nonsterile Nylon Syringe Filter, Pore:0.22(µm), Diameter:30(mm)	200
SFNY030045N	Nonsterile Nylon Syringe Filter, Pore:0.45(µm), Diameter:30(mm)	200

MS Syringe Filters Product List 2014

SFNY033022N	Nonsterile Nylon Syringe Filters, Pore:0.22(µm),Diameter:33(mm)	
SFNY033045N	Nonsterile Nylon Syringe Filters, Pore:0.45(µm),Diameter:33(mm)	200
SFNY050022N	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:50(mm)	12
SFNY050045N	Nonsterile Nylon Syringe Filters, Pore:0.45µm, Diameter:50(mm)	12
SFNY013010S	Sterile Nylon Syringe Filters, Pore:0.1(µm), Diameter:13(mm)	200
SFNY013022S	Sterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:13(mm)	200
SFNY013045S	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm)	200
SFNY013080S	Sterile Nylon Syringe Filters, Pore:0.8(µm), Diameter:13(mm)	200
SFNY013100S	Sterile Nylon Syringe Filters, Pore:1.0(µm), Diameter:13(mm)	200
SFNY013300S	Sterile Nylon Syringe Filters, Pore:3.0(µm), Diameter:13(mm)	200
SFNY013500S	Sterile Nylon Syringe Filters, Pore:5.0(µm), Diameter:13(mm)	200
SFNY025010S	Sterile Nylon Syringe Filters, Pore:0.1(µm), Diameter:25(mm)	200
SFNY025022S	Sterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm)	200
SFNY025045S	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:25(mm)	200
SFNY025080S	Sterile Nylon Syringe Filters, Pore:0.8(µm), Diameter:25(mm)	200
SFNY025100S	Sterile Nylon Syringe Filters, Pore:1.0(µm), Diameter:25(mm)	200
SFNY025200S	Sterile Nylon Syringe Filters, Pore:2.0(µm), Diameter:25(mm)	200
SFNY025300S	Sterile Nylon Syringe Filters, Pore:3.0(µm), Diameter:25(mm)	200
SFNY025500S	Sterile Nylon Syringe Filters, Pore:5.0(µm), Diameter:25(mm)	200
SFNY033022S	Sterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:33(mm)	200
SFNY033045S	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:33(mm)	200
SFNY013022NP	Nonsterile Nylon Syringe Filters,Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFNY013045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFNY025022NP	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFNY025045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), PP Prefilter	100
SFNY050022NP	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFNY050045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:50(mm), PP Prefilter	12
SFNY013022SP	Sterile Nylon Syringe Filters, Pore:0.22(µm),Diameter: 13(mm), PP Prefilter	100
SFNY013045SP	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFNY025022SP	Sterile Nylon Syringe Filters,Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFNY025045SP	Sterile Nylon Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), PP Prefilter	100
Pre Filter Syringe Filters : ALL Sizes		
SFCA025045NG	Nonsterile CA Syringe Filter, Pore:0.45(µm), Diameter:25(mm), GF prefilter	200
SFNY013010NG	Nonsterile Nylon Syringe Filters, Pore:0.10(µm), Diameter:13(mm), GF Prefilter	100
SFNY013022NG	Nonsterile Nylon Syringe Filters,Pore: 0.22(µm), Diameter: 13(mm), GF Prefilter	100
SFNY013022NP	Nonsterile Nylon Syringe Filters,Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFNY013022SG	Sterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFNY013022SP	Sterile Nylon Syringe Filters, Pore:0.22(µm),Diameter: 13(mm), PP Prefilter	100
SFNY013045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFNY013045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFNY013045SG	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFNY013045SP	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFNY025022NG	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm), GF Prefilter	100
SFNY025022NP	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100

MS Syringe Filters Product List 2014

Membrane Solutions

SFNY025022SG	Sterile Nylon Syringe Filters, Pore:0.22(µm),Diameter: 25(mm), GF Prefilter	
SFNY025022SP	Sterile Nylon Syringe Filters,Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFNY025045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFNY025045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), PP Prefilter	100
SFNY025045SG	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFNY025045SP	Sterile Nylon Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), PP Prefilter	100
SFNY025500NG	Nonsterile Nylon Syringe Filters, Pore:5.0(µm), Diameter:25(mm), GF Prefilter	100
SFNY030045NG	Nonsterile Nylon Syringe Filter, Pore:0.45(µm), Diameter:30(mm),GF Prefilter	100
SFNY050022NG	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:50(mm), GF Prefilter	12
SFNY050022NP	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFNY050045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:50(mm), GF Prefilter	12
SFNY050045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:50(mm), PP Prefilter	12
SFPES013022NG	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF prefilter	100
SFPES013022NP	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPES013022SG	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFPES013022SP	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPES013045NG	Nonsterile PES Syringe Filters, Pore:0.45(µm),Diameter: 13(mm), GF Prefilter	100
SFPES013045NP	Nonsterile PES Syringe Filters,Pore: 0.45(µm), Diameter:13(mm), PP Prefilter	100
SFPES013045SG	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFPES013045SP	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP prefilter	100
SFPES025022NG	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm), GF Prefilter	100
SFPES025022NP	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPES025022SG	Sterile PES Syringe Filters, Pore:0.02(µm), Diameter:25(mm), GF Prefilter	100
SFPES025022SP	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP prefilter	100
SFPES025045NG	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFPES025045NP	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP prefilter	100
SFPES025045SG	Sterile PES Syringe Filters, Pore:0.02(µm), Diameter:25(mm), GF Prefilter	100
SFPES025045SP	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP prefilter	100
SFPES050022NG	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:50(mm), GF Prefilter	12
SFPES050022NP	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFPES050045NG	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:50(mm), GF prefilter	12
SFPES050045NP	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:50(mm), PP Prefilter	12
SFPTFE013022NG	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFPTFE013022NP	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPTFE013022SG	Sterile PTFE Syringe Filters, Pore:0.22(µm),Diameter: 13(mm), GF Prefilter	100
SFPTFE013022SP	Sterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPTFE013045NG	Nonsterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFPTFE013045NP	Nonsterile PTFE Syringe Filters,Pore: 0.45(µm), Diameter:13(mm), PP Prefilter	100
SFPTFE013045SG	Sterile PTFE Syringe Filters,Pore: 0.45(µm), Diameter:13(mm), GF Prefilter	100
SFPTFE013045SP	Sterile PTFE Syringe Filters, Pore:0.45(µm),Diameter: 13(mm), PP prefilter	100
SFPTFE025022NG	Nonsterile PTFE Syringe Filters,Pore: 0.22(µm), Diameter:25(mm), GF Prefilter	100
SFPTFE025022NP	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPTFE025022SG	Sterile PTFE Syringe Filters, Pore:0.22(µm),Diameter: 25(mm), GF Prefilter	100
SFPTFE025022SP	Sterile PTFE Syringe Filters,Pore: 0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPTFE025045NG	Nonsterile PTFE Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), GF Prefilter	100

MS Syringe Filters Product List 2014

Membrane Solutions

SFPTFE025045NP	Nonsterile PTFE Syringe Filters, Pore:0.45(μm), Diameter:25(mm), PP Prefilter	
SFPTFE025045SG	Sterile PTFE Syringe Filters, Pore:0.45(μm),Diameter: 25(mm), GF Prefilter	100
SFPTFE025045SP	Sterile PTFE Syringe Filters, Pore:0.45(μm), Diameter:25(mm), PP Prefilter	100
SFPTFE050022NG	Nonsterile PTFE Syringe Filters, Pore:0.22(μm), Diameter:50(mm), GF Prefilter	12
SFPTFE050022NP	Nonsterile PTFE Syringe Filters, Pore:0.22(μm), Diameter:50(mm), PP Prefilter	12
SFPTFE050045NG	Nonsterile PTFE Syringe Filters, Pore:0.45(μm), Diameter:50(mm), GF prefilter	12
SFPTFE050045NP	Nonsterile PTFE Syringe Filters, Pore:0.45(μm), Diameter:50(mm), PP Prefilter	12
SFPVDF013022NG	Nonsterile PVDF Syringe Filter, Pore:0.22(μm), Diameter:13(mm), GF Prefilter	100
SFPVDF013022NP	Nonsterile PVDF Syringe Filter,Pore:0.22(μm), Diameter:13(mm), PP Prefilter	100
SFPVDF013022SG	Sterile PVDF Syringe Filter,Pore: 0.22(μm), Diameter:13(mm), GF Prefilter	100
SFPVDF013022SP	Sterile PVDF Syringe Filter,Pore: 0.22(μm), Diameter:13(mm), PP Prefilter	100
SFPVDF013045NG	Nonsterile PVDF Syringe Filter,Pore: 0.45(μm),Diameter: 13(mm), GF Prefilter	100
SFPVDF013045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(μm), Diameter:13(mm), PP Prefilter	100
SFPVDF013045SG	Sterile PVDF Syringe Filters, Pore:0.45(μm),Diameter: 13(mm), GF Prefilter	100
SFPVDF013045SP	Sterile PVDF Syringe Filter, Pore:0.45(μm), Diameter:13(mm), PP Prefilter	100
SFPVDF025022NG	Nonsterile PVDF Syringe Filter, Pore:0.22(μm),Diameter: 25(mm), GF Prefilter	100
SFPVDF025022NP	Nonsterile PVDF Syringe Filter,Pore: 0.22(μm), Diameter:25(mm), PP Prefilter	100
SFPVDF025022SG	Sterile PVDF Syringe Filter,Pore: 0.22(μm),Diameter: 25(mm), GF Prefilter	100
SFPVDF025022SP	Sterile PVDF Syringe Filter,Pore:0.22(μm), Diameter:25(mm), PP Prefilter	100
SFPVDF025045NG	Nonsterile PVDF Syringe Filters, Pore:0.45(μm), Diameter:25(mm), GF Prefilter	100
SFPVDF025045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(μm),Diameter: 25(mm), PP Prefilter	100
SFPVDF025045SG	Sterile PVDF Syringe Filter, Pore:0.45(μm), Diameter:25(mm), GF Prefilter	100
SFPVDF025045SP	Sterile PVDF Syringe Filter,Pore: 0.45(μm), Diameter:25(mm), PP Prefilter	100
SFPVDF050022NG	Nonsterile PVDF Syringe Filter, Pore:0.22(μm), Diameter:50(mm), GF Prefilter	12
SFPVDF050022NP	Nonsterile PVDF Syringe Filter, Pore:0.22(μm), Diameter:50(mm), PP Prefilter	12
SFPVDF050045NG	Nonsterile PVDF Syringe Filter, Pore:0.45(μm), Diameter:50(mm), GF Prefilter	12
SFPVDF050045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(μm), Diameter:50(mm), PP Prefilter	12
SFCA025045NG	Nonsterile CA Syringe Filter, Pore:0.45(μm), Diameter:25(mm), GF prefilter	200
SFNY013010NG	Nonsterile Nylon Syringe Filters, Pore:0.10(μm), Diameter:13(mm), GF Prefilter	100
SFNY013022NG	Nonsterile Nylon Syringe Filters,Pore: 0.22(μm), Diameter: 13(mm), GF Prefilter	100
SFNY013022NP	Nonsterile Nylon Syringe Filters,Pore:0.22(μm), Diameter:13(mm), PP Prefilter	100
SFNY013022SG	Sterile Nylon Syringe Filters, Pore:0.22(μm), Diameter:13(mm), GF Prefilter	100
SFNY013022SP	Sterile Nylon Syringe Filters, Pore:0.22(μm),Diameter: 13(mm), PP Prefilter	100
SFNY013045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(μm), Diameter:13(mm), GF Prefilter	100
SFNY013045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(μm), Diameter:13(mm), PP Prefilter	100
SFNY013045SG	Sterile Nylon Syringe Filters, Pore:0.45(μm), Diameter:13(mm), GF Prefilter	100
SFNY013045SP	Sterile Nylon Syringe Filters, Pore:0.45(μm), Diameter:13(mm), PP Prefilter	100
SFNY025022NG	Nonsterile Nylon Syringe Filters, Pore:0.22(μm), Diameter:25(mm), GF Prefilter	100
SFNY025022NP	Nonsterile Nylon Syringe Filters, Pore:0.22(μm), Diameter:25(mm), PP Prefilter	100
SFNY025022SG	Sterile Nylon Syringe Filters, Pore:0.22(μm),Diameter: 25(mm), GF Prefilter	100
SFNY025022SP	Sterile Nylon Syringe Filters,Pore:0.22(μm), Diameter:25(mm), PP Prefilter	100
SFNY025045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(μm), Diameter:25(mm), GF Prefilter	100
SFNY025045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(μm),Diameter: 25(mm), PP Prefilter	100
SFNY025045SG	Sterile Nylon Syringe Filters, Pore:0.45(μm), Diameter:25(mm), GF Prefilter	100
SFNY025045SP	Sterile Nylon Syringe Filters, Pore:0.45(μm),Diameter: 25(mm), PP Prefilter	100

MS Membrane Filters Product List 2014

SFNY025500NG	Nonsterile Nylon Syringe Filters, Pore:5.0(µm), Diameter:25(mm), GF Prefilter	
SFNY030045NG	Nonsterile Nylon Syringe Filter, Pore:0.45(µm), Diameter:30(mm),GF Prefilter	100
SFNY050022NG	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:50(mm), GF Prefilter	12
SFNY050022NP	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFNY050045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:50(mm), GF Prefilter	12
SFNY050045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:50(mm), PP Prefilter	12
SFPES013022NG	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF prefilter	100
SFPES013022NP	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPES013022SG	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFPES013022SP	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPES013045NG	Nonsterile PES Syringe Filters, Pore:0.45(µm),Diameter: 13(mm), GF Prefilter	100
SFPES013045NP	Nonsterile PES Syringe Filters,Pore: 0.45(µm), Diameter:13(mm), PP Prefilter	100
SFPES013045SG	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFPES013045SP	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP prefilter	100
SFPES025022NG	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm), GF Prefilter	100
SFPES025022NP	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPES025022SG	Sterile PES Syringe Filters, Pore:0.02(µm), Diameter:25(mm), GF Prefilter	100
SFPES025022SP	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP prefilter	100
SFPES025045NG	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFPES025045NP	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP prefilter	100
SFPES025045SG	Sterile PES Syringe Filters, Pore:0.02(µm), Diameter:25(mm), GF Prefilter	100
SFPES025045SP	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP prefilter	100
SFPES050022NG	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:50(mm), GF Prefilter	12
SFPES050022NP	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFPES050045NG	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:50(mm), GF prefilter	12
SFPES050045NP	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:50(mm), PP Prefilter	12
SFPTFE013022NG	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFPTFE013022NP	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPTFE013022SG	Sterile PTFE Syringe Filters, Pore:0.22(µm),Diameter: 13(mm), GF Prefilter	100
SFPTFE013022SP	Sterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPTFE013045NG	Nonsterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFPTFE013045NP	Nonsterile PTFE Syringe Filters,Pore: 0.45(µm), Diameter:13(mm), PP Prefilter	100
SFPTFE013045SG	Sterile PTFE Syringe Filters,Pore: 0.45(µm), Diameter:13(mm), GF Prefilter	100
SFPTFE013045SP	Sterile PTFE Syringe Filters, Pore:0.45(µm),Diameter: 13(mm), PP prefilter	100
SFPTFE025022NG	Nonsterile PTFE Syringe Filters,Pore: 0.22(µm), Diameter:25(mm), GF Prefilter	100
SFPTFE025022NP	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPTFE025022SG	Sterile PTFE Syringe Filters, Pore:0.22(µm),Diameter: 25(mm), GF Prefilter	100
SFPTFE025022SP	Sterile PTFE Syringe Filters,Pore: 0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPTFE025045NG	Nonsterile PTFE Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), GF Prefilter	100
SFPTFE025045NP	Nonsterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP Prefilter	100
SFPTFE025045SG	Sterile PTFE Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), GF Prefilter	100
SFPTFE025045SP	Sterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP Prefilter	100
SFPTFE050022NG	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:50(mm), GF Prefilter	12
SFPTFE050022NP	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFPTFE050045NG	Nonsterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:50(mm), GF prefilter	12

MS Syringe Filters Product List 2014

SFPTFE050045NP	Nonsterile PTFE Syringe Filters, Pore:0.45(μm), Diameter:50(mm), PP Prefilter	
SFPVDF013022NG	Nonsterile PVDF Syringe Filter, Pore:0.22(μm), Diameter:13(mm), GF Prefilter	100
SFPVDF013022NP	Nonsterile PVDF Syringe Filter,Pore:0.22(μm), Diameter:13(mm), PP Prefilter	100
SFPVDF013022SG	Sterile PVDF Syringe Filter,Pore: 0.22(μm), Diameter:13(mm), GF Prefilter	100
SFPVDF013022SP	Sterile PVDF Syringe Filter,Pore: 0.22(μm), Diameter:13(mm), PP Prefilter	100
SFPVDF013045NG	Nonsterile PVDF Syringe Filter,Pore: 0.45(μm),Diameter: 13(mm), GF Prefilter	100
SFPVDF013045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(μm), Diameter:13(mm), PP Prefilter	100
SFPVDF013045SG	Sterile PVDF Syringe Filters, Pore:0.45(μm),Diameter: 13(mm), GF Prefilter	100
SFPVDF013045SP	Sterile PVDF Syringe Filter, Pore:0.45(μm), Diameter:13(mm), PP Prefilter	100
SFPVDF025022NG	Nonsterile PVDF Syringe Filter, Pore:0.22(μm),Diameter: 25(mm), GF Prefilter	100
SFPVDF025022NP	Nonsterile PVDF Syringe Filter,Pore: 0.22(μm), Diameter:25(mm), PP Prefilter	100
SFPVDF025022SG	Sterile PVDF Syringe Filter,Pore: 0.22(μm),Diameter: 25(mm), GF Prefilter	100
SFPVDF025022SP	Sterile PVDF Syringe Filter,Pore:0.22(μm), Diameter:25(mm), PP Prefilter	100
SFPVDF025045NG	Nonsterile PVDF Syringe Filters, Pore:0.45(μm), Diameter:25(mm), GF Prefilter	100
SFPVDF025045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(μm),Diameter: 25(mm), PP Prefilter	100
SFPVDF025045SG	Sterile PVDF Syringe Filter, Pore:0.45(μm), Diameter:25(mm), GF Prefilter	100
SFPVDF025045SP	Sterile PVDF Syringe Filter,Pore: 0.45(μm), Diameter:25(mm), PP Prefilter	100
SFPVDF050022NG	Nonsterile PVDF Syringe Filter, Pore:0.22(μm), Diameter:50(mm), GF Prefilter	12
SFPVDF050022NP	Nonsterile PVDF Syringe Filter, Pore:0.22(μm), Diameter:50(mm), PP Prefilter	12
SFPVDF050045NG	Nonsterile PVDF Syringe Filter, Pore:0.45(μm), Diameter:50(mm), GF Prefilter	12
SFPVDF050045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(μm), Diameter:50(mm), PP Prefilter	12

SuperPure HPLC Syringe Filters : All Sizes - All Materials

SFNY017022N	SuperPure HPLC, Nylon66, 0.22um,17mm	100/pk
SFNY017045N	SuperPure HPLC, Nylon66, 0.45um,17mm	100/pk
SFNY030022N	SuperPure HPLC, Nylon66, 0.22um,30mm	100/pk
SFNY030045N	SuperPure HPLC, Nylon66, 0.45um,30mm	100/pk
SFPES017022N	SuperPure HPLC, PES, 0.22um,17mm	100/pk
SFPES017045N	SuperPure HPLC, PES, 0.45um,17mm	100/pk
SFPES030022N	SuperPure HPLC, PES, 0.22um,30mm	100/pk
SFPES030045N	SuperPure HPLC, PES, 0.45um,30mm	100/pk
SFMCE017022N	SuperPure HPLC, MCE, 0.22um,17mm	100/pk
SFMCE017045N	SuperPure HPLC, MCE , 0.45um,17mm	100/pk
SFMCE030022N	SuperPure HPLC, MCE, 0.22um,30mm	100/pk
SFMCE030045N	SuperPure HPLC, MCE , 0.45um,30mm	100/pk
SFPTFE017022NB	SuperPure HPLC, PTFE, 0.22um,17mm	100/pk
SFPTFE017045NB	SuperPure HPLC, PTFE , 0.45um,17mm	100/pk
SFPTFE030022NB	SuperPure HPLC, PTFE, 0.22um,30mm	100/pk
SFPTFE030045NB	SuperPure HPLC, PTFE, 0.45um,30mm	100/pk
SFPVDF017022N	SuperPure HPLC, PVDF, 0.22um,17mm	100/pk
SFPVDF017045N	SuperPure HPLC, PVDF, 0.45um,17mm	100/pk
SFPVDF030022N	SuperPure HPLC, PVDF, 0.22um,30mm	100/pk
SFPVDF030045N	SuperPure HPLC, PVDF, 0.45um,30mm	100/pk

SteriBio Syringe Filters : All Sizes - All Material		
SFCA013022S	Sterile CA Syringe Filter, Pore:0.22(µm), Diameter:13(mm)	200
SFCA013045S	Sterile CA Syringe Filter, Pore:0.45(µm), Diameter:13(mm)	200
SFCA025022S	Sterile CA Syringe Filter, Pore:0.22(µm),Diameter:25(mm)	200
SFCA025045S	Sterile CA Syringe Filter, Pore:0.45(µm), Diameter:25(mm)	200
SFCA025500S	Sterile CA Syringe Filter, Pore:5.0(µm), Diameter:25(mm)	100
SFCA033022S	Sterile CA Syringe Filter, Pore:0.22(µm), Diameter:33(mm)	200
SFCA033045S	Sterile CA Syringe Filter, Pore:0.45(µm), Diameter:33(mm)	200
SFCA033500S	Sterile CA Syringe Filter, Pore:5.0(µm),Diameter:33(mm)	200
SFPES013022S	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm)	200
SFPES013045S	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:13(mm)	200
SFPES025010S	Sterile PES Syringe Filters, Pore:0.1(µm), Diameter:25(mm)	100
SFPES025022S	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm)	200
SFPES025045S	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm)	100
SFPES033022S	Sterile PES Syringe Filters, Pore:0.22(µm),Diameter:33(mm)	200
SFPES033045S	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:33(mm)	200
MS G-MP MultiLayer Syringe Filter : 25mm Diam; 75psi max, PP Body		
10um Glass Fibre; 1.0um PP; 0.7um Glass Fibre / media layer		
SFNY025022NM	Nylon, 0.22um	100
SFNY025045NM	Nylon, 0.45um	100
SFPVDF025022NM	PVDF, 0.22um	100
SFPVDF025045NM	PVDF, 0.45um	100
SFPTFE025022NM	PTFE, 0.22um	100
SFPTFE025045NM	PTFE, 0.45um	100
SFPES025022NM	PES, 0.22um	100
SFPES025045NM	PES, 0.45 um	100
SFPP025022NM	PP, 0.22um	100
SFPP025045NM	PP, 0.45um	100
SFCA025022NM	CA, 0.22um	100
SFCA025045NM	CA, 0.45um	100
Reusable Syringe Filter Holder		
SFRPP13	Reusable Syringe Filter, Polypropylene Housing,Diameter:13(mm)	10
SFRSS13	Reusable Syringe Filter, Stainless Steel Housing,Diameter:13(mm)	10
SFRPP25	Reusable Syringe Filter, Polypropylene Housing,Diameter:25(mm)	10
SFRSS25	Reusable Syringe Filter, Stainless Steel Housing,Diameter:25(mm)	10
SFRPP35	Reusable Syringe Filter, Polypropylene Housing,Diameter:35(mm)	10
SFRPP50	Reusable Syringe Filter, Polypropylene Housing,Diameter:50(mm)	10

MS Membrane Filters Product List 2014

Item#	Type	Description	Packaging
Nylon Membrane Filter			
MFNY013022	Pore size: 0.22um	Nylon Membrane Filter, Pore:0.22(μm), Diameter:13(mm)	400
MFNY025022		Nylon Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFNY047022		Nylon Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFNY090022		Nylon Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFNY142022		Nylon Membrane Filter, Pore:0.22(μm), Diameter:142(mm)	50
MFNY293022		Nylon Membrane Filter, Pore:0.22(μm), Diameter:293(mm)	25
MFNY013045	Pore size: 0.45um	Nylon Membrane Filter, Pore:0.45(μm), Diameter:13(mm)	400
MFNY025045		Nylon Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200
MFNY047045		Nylon Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFNY090045		Nylon Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
MFNY142045		Nylon Membrane Filter, Pore:0.45(μm), Diameter:142(mm)	50
MFNY293045		Nylon Membrane Filter, Pore:0.45(μm), Diameter:293(mm)	25
MFNY013080	Pore size: 0.8um	Nylon Membrane Filter, Pore:0.8(μm), Diameter:13(mm)	400
MFNY025080		Nylon Membrane Filter, Pore:0.8(μm), Diameter:25(mm)	200
MFNY047080		Nylon Membrane Filter, Pore:0.8(μm), Diameter:47(mm)	200
MFNY090080		Nylon Membrane Filter, Pore:0.8(μm), Diameter:90(mm)	200
MFNY142080		Nylon Membrane Filter, Pore:0.8(μm), Diameter:142(mm)	200
MFNY293080		Nylon Membrane Filter, Pore:0.8(μm), Diameter:293(mm)	200
MFNY013100	Pore size: 1.2um	Nylon Membrane Filter, Pore:1.2(μm), Diameter:13(mm)	400
MFNY025100		Nylon Membrane Filter, Pore:1.2(μm), Diameter:25(mm)	200
MFNY047100		Nylon Membrane Filter, Pore:1.2(μm), Diameter:47(mm)	200
MFNY090100		Nylon Membrane Filter, Pore:1.2(μm), Diameter:90(mm)	100
MFNY142100		Nylon Membrane Filter, Pore:1.2(μm), Diameter:142(mm)	50
MFNY293100		Nylon Membrane Filter, Pore:1.2(μm), Diameter:293(mm)	25
MFNY013300	Pore size: 3.0um	Nylon Membrane Filter, Pore:3.0(μm), Diameter:13(mm)	400
MFNY025300		Nylon Membrane Filter, Pore:3.0(μm), Diameter:25(mm)	200
MFNY047300		Nylon Membrane Filter, Pore:3.0(μm), Diameter:47(mm)	200
MFNY090300		Nylon Membrane Filter, Pore:3.0(μm), Diameter:90(mm)	100
MFNY142300		Nylon Membrane Filter, Pore:3.0(μm), Diameter:142(mm)	50
MFNY293300		Nylon Membrane Filter, Pore:3.0(μm), Diameter:293(mm)	25
MFNY013500	Pore size: 5.0um	Nylon Membrane Filter, Pore:5.0(μm), Diameter:13(mm)	400
MFNY025500		Nylon Membrane Filter, Pore:5.0(μm), Diameter:25(mm)	200
MFNY047500		Nylon Membrane Filter, Pore:5.0(μm), Diameter:47(mm)	200
MFNY090500		Nylon Membrane Filter, Pore:5.0(μm), Diameter:90(mm)	100
MFNY142500		Nylon Membrane Filter, Pore:5.0(μm), Diameter:142(mm)	50
MFNY293500		Nylon Membrane Filter, Pore:5.0(μm), Diameter:293(mm)	25
Hydrophobic PTFE Membrane Filter			
MFPTFE013022B	Pore size: 0.22um	Hydrophobic PTFE Membrane Filter, Pore:0.22(μm), Diameter:13(mm)	400
MFPTFE025022B		Hydrophobic PTFE Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFPTFE047022B		Hydrophobic PTFE Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFPTFE090022B		Hydrophobic PTFE Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFPTFE142022B		Hydrophobic PTFE Membrane Filter, Pore:0.22(μm), Diameter:142(mm)	50
MFPTFE013045B	Pore size: 0.45um	Hydrophobic PTFE Membrane Filter, Pore:0.45(μm), Diameter:13(mm)	400
MFPTFE025045B		Hydrophobic PTFE Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200

MS Membrane Filters Product List 2014

MFPTFE047045B		Hydrophobic PTFE Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFPTFE090045B		Hydrophobic PTFE Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
MFPTFE142045B		Hydrophobic PTFE Membrane Filter, Pore:0.45(μm), Diameter:142(mm)	50
MFPTFE025010B	Pore size: 0.1μm	Hydrophobic PTFE Membrane Filter, Pore:0.1(μm), Diameter:25(mm)	200
MFPTFE047010B		Hydrophobic PTFE Membrane Filter, Pore:0.1(μm), Diameter:47(mm)	200
MFPTFE090010B		Hydrophobic PTFE Membrane Filter, Pore:0.1(μm), Diameter:90(mm)	100
MFPTFE142010B		Hydrophobic PTFE Membrane Filter, Pore:0.1(μm), Diameter:142(mm)	50
MFPTFE025100B	Pore size: 1.2μm	Hydrophobic PTFE Membrane Filter, Pore:1.2(μm), Diameter:25(mm)	200
MFPTFE047100B		Hydrophobic PTFE Membrane Filter, Pore:1.2(μm), Diameter:47(mm)	200
MFPTFE090100B		Hydrophobic PTFE Membrane Filter, Pore:1.2(μm), Diameter:90(mm)	100
MFPTFE142100B		Hydrophobic PTFE Membrane Filter, Pore:1.2(μm), Diameter:142(mm)	50
MFPTFE025300B	Pore size: 3.0μm	Hydrophobic PTFE Membrane Filter, Pore:3.0(μm), Diameter:25(mm)	200
MFPTFE047300B		Hydrophobic PTFE Membrane Filter, Pore:3.0(μm), Diameter:47(mm)	200
MFPTFE090300B		Hydrophobic PTFE Membrane Filter, Pore:3.0(μm), Diameter:90(mm)	100
MFPTFE142300B		Hydrophobic PTFE Membrane Filter, Pore:3.0(μm), Diameter:142(mm)	50
MFPTFE025300B	Pore size: 5.0μm	Hydrophobic PTFE Membrane Filter, Pore:3.0(μm), Diameter:25(mm)	200
MFPTFE047500B		Hydrophobic PTFE Membrane Filter, Pore:5.0(μm), Diameter:47(mm)	200
MFPTFE090500B		Hydrophobic PTFE Membrane Filter, Pore:5.0(μm), Diameter:90(mm)	100
MFPTFE142500B		Hydrophobic PTFE Membrane Filter, Pore:5.0(μm), Diameter:142(mm)	50
Hydrophilic PTFE Membrane Filter			
MFPTFE013022L	Pore size: 0.22μm	Hydrophilic PTFE Membrane Filter, Pore:0.22(μm), Diameter:13(mm)	400
MFPTFE025022L		Hydrophilic PTFE Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFPTFE047022L		Hydrophilic PTFE Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFPTFE090022L		Hydrophilic PTFE Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFPTFE013045L	Pore size: 0.45μm	Hydrophilic PTFE Membrane Filter, Pore:0.45(μm), Diameter:13(mm)	400
MFPTFE025045L		Hydrophilic PTFE Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200
MFPTFE047045L		Hydrophilic PTFE Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFPTFE090045L		Hydrophilic PTFE Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
CA Membrane Filter			
MFCA013022	Pore size: 0.22μm	CA Membrane Filter, Pore:0.22(μm), Diameter:13(mm)	400
MFCA025022		CA Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFCA047022		CA Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFCA090022		CA Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFCA142022		CA Membrane Filter, Pore:0.22(μm), Diameter:142(mm)	50
MFCA013045	Pore size: 0.45μm	CA Membrane Filter, Pore:0.45(μm), Diameter:13(mm)	400
MFCA025045		CA Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200
MFCA047045		CA Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFCA090045		CA Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
MFCA142045		CA Membrane Filter, Pore:0.45(μm), Diameter:142(mm)	50
MFCA013080	Pore size: 0.8μm	CA Membrane Filter, Pore:0.8(μm), Diameter:13(mm)	400
MFCA025080		CA Membrane Filter, Pore:0.8(μm), Diameter:25(mm)	200
MFCA047080		CA Membrane Filter, Pore:0.8(μm), Diameter:47(mm)	200
MFCA090080		CA Membrane Filter, Pore:0.8(μm), Diameter:90(mm)	100
MFCA142080		CA Membrane Filter, Pore:0.8(μm), Diameter:142(mm)	50
CN Membrane Filter			

MS Membrane Filters Product List 2014

MFCN047022	Pore size: 0.22um	CN Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFCN047045		CN Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFCN047080		CN Membrane Filter, Pore:0.8(μm), Diameter:47(mm)	200
MCE Membrane Filter			
MFMC013022	Pore size: 0.22um	MCE Membrane Filter, Pore:0.22(μm), Diameter:13(mm)	400
MFMC025022		MCE Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFMC037022		MCE Membrane Filter, Pore:0.22(μm), Diameter:37(mm)	200
MFMC047022		MCE Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFMC090022		MCE Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFMC142022		MCE Membrane Filter, Pore:0.22(μm), Diameter:142(mm)	50
MFMC013045	Pore size: 0.45um	MCE Membrane Filter, Pore:0.45(μm), Diameter:13(mm)	400
MFMC025045		MCE Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200
MFMC037045		MCE Membrane Filter, Pore:0.45(μm), Diameter:37(mm)	200
MFMC047045		MCE Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFMC090045		MCE Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
MFMC142045		MCE Membrane Filter, Pore:0.45(μm), Diameter:142(mm)	50
MFMC013080	Pore size: 0.8um	MCE Membrane Filter, Pore:0.8(μm), Diameter:13(mm)	400
MFMC025080		MCE Membrane Filter, Pore:0.8(μm), Diameter:25(mm)	200
MFMC037080		MCE Membrane Filter, Pore:0.8(μm), Diameter:37(mm)	200
MFMC047080		MCE Membrane Filter, Pore:0.8(μm), Diameter:47(mm)	200
MFMC090080		MCE Membrane Filter, Pore:0.8(μm), Diameter:90(mm)	100
MFMC142080		MCE Membrane Filter, Pore:0.8(μm), Diameter:142(mm)	50
MFMC013100	Pore size: 1.2um	MCE Membrane Filter, Pore:1.2(μm), Diameter:13(mm)	400
MFMC025100		MCE Membrane Filter, Pore:1.2(μm), Diameter:25(mm)	200
MFMC037100		MCE Membrane Filter, Pore:1.2(μm), Diameter:37(mm)	200
MFMC047100		MCE Membrane Filter, Pore:1.2(μm), Diameter:47(mm)	200
MFMC090100		MCE Membrane Filter, Pore:1.2(μm), Diameter:90(mm)	100
MFMC142100		MCE Membrane Filter, Pore:1.2(μm), Diameter:142(mm)	50
MFMC013300	Pore size: 3.0um	MCE Membrane Filter, Pore:3.0(μm), Diameter:13(mm)	400
MFMC025300		MCE Membrane Filter, Pore:3.0(μm), Diameter:25(mm)	200
MFMC037300		MCE Membrane Filter, Pore:3.0(μm), Diameter:37(mm)	200
MFMC047300		MCE Membrane Filter, Pore:3.0(μm), Diameter:47(mm)	200
MFMC090300		MCE Membrane Filter, Pore:3.0(μm), Diameter:90(mm)	100
MFMC142300		MCE Membrane Filter, Pore:3.0(μm), Diameter:142(mm)	50
MCE Membrane Filter, White, Gridded, Non-sterile			
MFMC047022GWN	Pore size: 0.22um	MCE Gridded Membrane Filter, White, 0.22(μm), 47(mm)	200
MFMC047045GWN	Pore size: 0.45um	MCE Gridded Membrane Filter, White, 0.45(μm), 47(mm)	200
MFMC047080GWN	Pore size: 0.80um	MCE Gridded Membrane Filter, White, 0.80(μm), 47(mm)	200
MCE Membrane Filter, Black, Gridded, Non-sterile			
MFMC047022GBN	Pore size: 0.22um	MCE Gridded Membrane Filter, Black, 0.22(μm), 47(mm)	200
MFMC047045GBN	Pore size: 0.45um	MCE Gridded Membrane Filter, Black, 0.45(μm), 47(mm)	200
MFMC047080GBN	Pore size: 0.80um	MCE Gridded Membrane Filter, Black, 0.80(μm), 47(mm)	200
PES Membrane Filter			
MFPE013022	Pore size: 0.22um	PES Membrane Filter, Pore:0.22(μm), Diameter:13(mm)	400
MFPE025022	Pore size: 0.22um	PES Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFPE047022	Pore size: 0.22um	PES Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200

MS Membrane Filters Product List 2014

MFPE090022	Pore size: 0.22m	PES Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFPE142022	Pore size: 0.22m	PES Membrane Filter, Pore:0.22(μm), Diameter:142(mm)	50
MFPE013045	Pore size: 0.45m	PES Membrane Filter, Pore:0.45(μm), Diameter:13(mm)	400
MFPE025045		PES Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200
MFPE047045		PES Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFPE090045		PES Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
MFPE142045		PES Membrane Filter, Pore:0.45(μm), Diameter:142(mm)	50
MFPE013010	Pore size: 0.1m	PES Membrane Filter, Pore:0.1(μm), Diameter:13(mm)	400
MFPE025010		PES Membrane Filter, Pore:0.1(μm), Diameter:25(mm)	200
MFPE047010		PES Membrane Filter, Pore: 0.1(μm), Diameter: 47(mm)	200
MFPE090010		PES Membrane Filter, Pore:0.1(μm), Diameter:90(mm)	100
PVDF Membrane Filter			
MFPVDF013022	Pore size: 0.22m	PVDF Membrane Filter, Pore: 0.22(μm), Diameter: 13(mm)	400
MFPVDF025022		PVDF Membrane Filter, Pore: 0.22(μm), Diameter: 25(mm)	200
MFPVDF047022		PVDF Membrane Filter, Pore: 0.22(μm), Diameter: 47(mm)	200
MFPVDF090022		PVDF Membrane Filter, Pore: 0.22(μm), Diameter: 90(mm)	100
MFPVDF142022		PVDF Membrane Filter, Pore: 0.22(μm), Diameter: 142(mm)	50
MFPVDF013045	Pore size: 0.45m	PVDF Membrane Filter, Pore: 0.45(μm), Diameter: 13(mm)	400
MFPVDF025045		PVDF Membrane Filter, Pore: 0.45(μm), Diameter: 25(mm)	200
MFPVDF047045		PVDF Membrane Filter, Pore: 0.45(μm), Diameter: 47(mm)	200
MFPVDF090045		PVDF Membrane Filter, Pore: 0.45(μm), Diameter: 90(mm)	100
MFPVDF142045		PVDF Membrane Filter, Pore: 0.45(μm), Diameter: 142(mm)	50
PP Membrane Filter			
MFPP025022	Pore size: 0.22m	PP Membrane Filter, Pore:0.22(μm), Diameter:25(mm)	200
MFPP047022		PP Membrane Filter, Pore:0.22(μm), Diameter:47(mm)	200
MFPP090022		PP Membrane Filter, Pore:0.22(μm), Diameter:90(mm)	100
MFPP142022		PP Membrane Filter, Pore:0.22(μm), Diameter:142(mm)	50
MFPP250022		PP Membrane Filter, Pore:0.22(μm), Diameter:250(mm)	50
MFPP293022		PP Membrane Filter, Pore:0.22(μm), Diameter:293(mm)	25
MFPP025045	Pore size: 0.45m	PP Membrane Filter, Pore:0.45(μm), Diameter:25(mm)	200
MFPP047045		PP Membrane Filter, Pore:0.45(μm), Diameter:47(mm)	200
MFPP090045		PP Membrane Filter, Pore:0.45(μm), Diameter:90(mm)	100
MFPP142045		PP Membrane Filter, Pore:0.45(μm), Diameter:142(mm)	50
MFPP250045		PP Membrane Filter, Pore:0.45(μm), Diameter:250(mm)	50
MFPP293045		PP Membrane Filter, Pore:0.45(μm), Diameter:293(mm)	25

MS Nylon Mesh Filters Product List 2014

Nylon Mesh Filters			
MENY090041	41um mesh	Nylon mesh filter,90mm, 41µm,100Pcs/pack	100
MENY303041	41um mesh	Nylon mesh filter, 30cm×3m,41µm, 1Pc/Pack	1
MENY025060	60um mesh	Nylon mesh filter,25mm,60µm,100Pcs/pack	100
MENY047060	60um mesh	Nylon mesh filter,47mm,60µm,100Pcs/pack	100
MENY090060	60um mesh	Nylon mesh filter,90mm, 60µm,100Pcs/pack	100
MENY303060	60um mesh	Nylon mesh filter, 30cm×3m,60µm, 1Pc/Pack	1
MENY025080	80um mesh	Nylon mesh filter,25mm,80µm,100Pcs/pack	100
MENY047080	80um mesh	Nylon mesh filter,47mm,80µm,100Pcs/pack	100
MENY090080	80um mesh	Nylon mesh filter,90mm, 80µm,100Pcs/pack	100
MENY303080	80um mesh	Nylon mesh filter, 30cm×3m,80µm, 1Pc/Pack	1
MENY025100	100um mesh	Nylon mesh filter,25mm,100µm,100Pcs/pack	100
MENY047100	100um mesh	Nylon mesh filter,47mm,100µm,100Pcs/pack	100
MENY090100	100um mesh	Nylon mesh filter,90mm, 100µm,100Pcs/pack	100
MENY303100	100um mesh	Nylon mesh filter, 30cm×3m,100µm, 1Pc/Pack	1
MENY025120	120um mesh	Nylon mesh filter,25mm,120µm,100Pcs/pack	100
MENY047120	120um mesh	Nylon mesh filter,47mm,120µm,100Pcs/pack	100
MENY090120	120um mesh	Nylon mesh filter,90mm, 120µm,100Pcs/pack	100
MENY303120	120um mesh	Nylon mesh filter, 30cm×3m,120µm, 1Pc/Pack	1
MENY025140	140um mesh	Nylon mesh filter,25mm,140µm,100Pcs/pack	100
MENY047140	140um mesh	Nylon mesh filter,47mm,140µm,100Pcs/pack	100
MENY090140	140um mesh	Nylon mesh filter,90mm, 140µm,100Pcs/pack	100
MENY303140	140um mesh	Nylon mesh filter, 30cm×3m,140µm, 1Pc/Pack	1
MENY025160	160um mesh	Nylon mesh filter,25mm,160µm,100Pcs/pack	100
MENY047160	160um mesh	Nylon mesh filter,47mm,160µm,100Pcs/pack	100
MENY090160	160um mesh	Nylon mesh filter,90mm,160µm,100Pcs/pack	100
MENY303160	160um mesh	Nylon mesh filter, 30cm×3m,160µm,100Pcs/pack	100
MENY025180	180um mesh	Nylon mesh filter,25mm,180µm,100Pcs/pack	100
MENY047180	180um mesh	Nylon mesh filter,47mm,180µm,100Pcs/pack	100
MENY090180	180um mesh	Nylon mesh filter,90mm,180µm,100Pcs/pack	100
MENY303180	180um mesh	Nylon mesh filter, 30cm×3m,180µm,100Pcs/pack	100

Note: 1) All of the price listed below are **Ex.work price**.

4) Membrane Solutions reserves itself the right to modify prices and content of this document without previous notice.

Cat.No.	Membrane Material	Filter Pore Size(μm)	Qty(pcs/cs)	Case Size(cm)
Filter System (Membrane Diameter: 50mm, Funnel Capacity : 250ml,				
Receiver Bottle Capacity : 250ml)				
VFPPE122250	PES	0.22	12	33*33*30
VFPPE145250		0.45	12	
VFPMCE122250	MCE	0.22	12	
VFPMCE145250		0.45	12	
VFPCA122250	CA	0.22	12	
VFPCA145250		0.45	12	
VFPNY122250	Nylon	0.22	12	
VFPNY145250		0.45	12	
VFPVDF122250	PVDF	0.22	12	
VFPVDF145250		0.45	12	
Filter System (Membrane Diameter : 50mm, Funnel Capacity : 250ml,				
Receiver Bottle Capacity : 500ml)				
VFPPE122500	PES	0.22	12	33*33*30
VFPPE145500		0.45	12	
VFPMCE122500	MCE	0.22	12	
VFPMCE145500		0.45	12	
VFPCA122500	CA	0.22	12	
VFPCA145500		0.45	12	
VFPNY122500	Nylon	0.22	12	
VFPNY145500		0.45	12	
VFPVDF122500	PVDF	0.22	12	
VFPVDF145500		0.45	12	
Cat.No.	Membrane Material	Filter Pore Size(μm)	Qty (pcs/cs)	Case Size(cm)
Filter System (Membrane Diameter : 90mm, Funnel Capacity : 500ml,				
Receiver Bottle Capacity : 500ml)				
VFPPE222500	PES	0.22	12	52.5*38.5*29.5
VFPPE245500		0.45	12	
VFPMCE222500	MCE	0.22	12	
VFPMCE245500		0.45	12	
VFPCA222500	CA	0.22	12	
VFPCA245500		0.45	12	
VFPNY222500	Nylon	0.22	12	
VFPNY245500		0.45	12	
VFPVDF222500	PVDF	0.22	12	
VFPVDF245500		0.45	12	

Filter System (Membrane Diameter : 90mm, Funnel Capacity : 500ml, Receiver Bottle Capacity:1000ml)					
VFPPEs2221000	PES	0.22	12	52.5*38.5*29.5	
VFPPEs2451000		0.45	12		
VFPMCE2221000	MCE	0.22	12		
VFPMCE2451000		0.45	12		
VFPCA2221000	CA	0.22	12		
VFPCA2451000		0.45	12		
VFPNY2221000	Nylon	0.22	12		
VFPNY2451000		0.45	12		
VFPVPDF2221000	PVDF	0.22	12		
VFPVPDF2451000		0.45	12		
Cat.No.	Membrane Material	Filter Pore Size(μm)	Qty (pcs/cs)	Case Size(cm)	
Funnel Capacity : 250ml, Membrane Diameter : 50mm					
VFPPEs122250F	PES	0.22	24	33*33*30	
VFPPEs145250F		0.45	24		
VFPMCE122250F	MCE	0.22	24		
VFPMCE145250F		0.45	24		
VFPCA122250F	CA	0.22	24		
VFPCA145250F		0.45	24		
VFPNY122250F	Nylon	0.22	24		
VFPNY145250F		0.45	24		
VFPVPDF122250F	PVDF	0.22	24		
VFPVPDF145250F		0.45	24		
Funnel Capacity : 500ml, Membrane Diameter : 90mm					
VFPPEs222500F	PES	0.22	24	52.5*38.5*29.5	
VFPPEs245500F		0.45	24		
VFPMCE222500F	MCE	0.22	24		
VFPMCE245500F		0.45	24		
VFPCA222500F	CA	0.22	24		
VFPCA245500F		0.45	24		
VFPNY222500F	Nylon	0.22	24		
VFPNY245500F		0.45	24		
VFPVPDF222500F	PVDF	0.22	24		
VFPVPDF245500F		0.45	24		
Cat.No.	Description		Qty(pcs/cs)	Case Size(cm)	
Receiver Bottle					
VFP250B	250ml polystyrene receiver bottle, sterile		24	33*33*30	
VFP500B	500ml polystyrene receiver bottle, sterile		24	52.5*38.5*29.5	
VFP1000B	1000ml polystyrene receiver bottle, sterile		24	52.5*38.5*29.5	

Note: 1) All of the price listed below are **Ex work price** below are Ex work price.

2) All the price listed below are for Membrane Solutions' Brand label.

3) **Order Qty** means **total items** in one order.

4) Membrane Solutions reserves itself the right to modify prices and content of this document without previous notice.

	Item#	Description	Packaging
Glass Vacuum Filter	VFG030010	1) 300ml Glass Funnel with cover	1 set
		2) 1000ml Glass Solvent Collection	
		3) Stainless Steel Clamp	
		4) Solid Glass Frit	
		Packaging: 30cmX30cmX25cm	
Glass Vacuum Filter	VFG030020	1) 300ml Glass Funnel with cover	1 set
		2) 2000ml Glass Solvent Collection	
		3) Stainless Steel Clamp	
		4) Solid Glass Frit	
		Packaging: 30cmX30cmX25cm	
Glass Vacuum Filter	VFG050020	1) 500ml Glass Funnel with cover	1 set
		2) 2000ml Glass Solvent Collection	
		3) Stainless Steel Clamp	
		4) Solid Glass Frit	
		Packaging: 30cmX30cmX25cm	
Support Screen	VFG147SS	Pyrex glass support screen (1000ML)	1pcs
	VFG047SS	Pyrex glass support screen (2000ML)	1pcs
	VFG010RB	Solvent Collection Bottle, 1000ml	1pcs
	VFG020RB	Solvent Collection Bottle, 2000ml	1pcs

Vacuum Pump list 2013

	Item#	Description	Packaging
Vacuun Pump	VPJ0201	Vacuum Pump, 300mbar, 12L/min GM-0.20((Oil free))	1set
Vacuun Pump	VPJ0332	Vacuum Pump, 250mbar, 20L/min GM-0.33A(Oil free),	1set
Vacuun Pump	VPJ0333	Vacuum Pump, 50mbar, 20L/min	1set
Vacuun Pump	VPJ0501	Vacuum Pump, 250mbar, 30L/min	1set
Vacuun Pump	VPJ0502	Vacuum Pump ,50mbar 30L/min	1set
Vacuun Pump	VPJ1001	Vacuum Pump, 250mbar, 60L/min GM-1.0A(Oil free),	1set

MS Glass Chromatography Vials Product List 2014

MS Items	Type	Products Name/Size	Pcs/pk
LBSV012C	8-425 & 9-425 2ml Clear Vial	2ml Clear vial, 8-425 screw top	100pcs
LBSV032C		2ml Clear vial, 8-425 screw top, graduated with writing area	100pcs
LBSV002C	12x32mm	2ml Clear vial, 9-425 screw top	100pcs
LBSV022C		2ml Clear vial, 9-425 screw top, graduated with writing area	100pcs
LBSV042C		2ml Clear vial, crimp top	100pcs
LBSV062C		2ml Clear vial, crimp top, graduated with writing area	100pcs
LBSV052C		2ml Clear vial, snap top	100pcs
LBSV072C		2ml Clear vial, snap top, graduated with writing area	100pcs
LBSV012A	8-425 & 9-425 2ml Amber Vial	2ml Amber vial, 8-425 screw top	100pcs
LBSV032A		2ml Amber vial, 8-425 screw top, graduated with writing area	100pcs
LBSV002A	12x32mm	2ml Amber vial, 9-425 screw top	100pcs
LBSV022A		2ml Amber vial, 9-425 screw top, graduated with writing area	100pcs
LBSV042A		2ml Amber vial, crimp top	100pcs
LBSV062A		2ml Amber vial, crimp top, graduated with writing area	100pcs
LBSV052A		2ml Amber vial, snap top	100pcs
LBSV072A		2ml Amber vial, snap top, graduated with writing area	100pcs
LBSV102CSS	Cap & Septa for 2ml Sample Vials	Red PTFE/white silicone septa + Blue screw cap with hole, for 2ml 9-425 screw top vial	100pcs
LBSV002CRS		White PTFE/red silicone septa + Blue screw cap with hole, for 2ml 9-425 screw top vial	100pcs
LBSV012CSS	12x32mm	Red PTFE/white silicone septa + Black screw cap with hole, for 2ml 8-425 screw top vial	100pcs
LBSV012CRS		White PTFE/red silicone septa + Black screw cap with hole, for 2ml 8-425 screw top vial	100pcs
LBSV222CSS		Blue PTFE/white silicone septa, Pre-slit + Blue screw cap with hole, 2ml 9-425 screw top vial	100pcs
LBSV221CSS		Pre-slit, Blue PTFE/White silicone septa+ Black screw cap with hole, for 8-425 screw top vial	100pcs
LBSV122CRS		Red PTFE/white silicone septa + Aluminium cap, for 2ml crimp top vial	100pcs
LBSV112CRS		White PTFE/red silicone septa+ Aluminium cap, for 2ml crimp top vial	100pcs
LBSV212CRS		Pre-silt, White PTFE/red silicone septa+ Aluminium cap, for 2ml crimp top vial	100pcs
LBSV222CRS		Pre-silt, Red PTFE/white silicone septa + Aluminium cap, for 2ml crimp top vial	100pcs
LBSV202CSS		White PTFE/Red silicone septa + snap cap with hole, for 2ml snap top vial	100pcs
LBSV302CSS		Red PTFE/White silicone septa + snap cap with hole, for 2ml snap top vial	100pcs
LBSV212CSS		Pre-slit, White PTFE/Red silicone septa + snap cap with hole, for 2ml snap top vial	100pcs
LBSV242CSS		Pre-slit, Red PTFE/white silicone septa + snap cap with hole, for 2ml snap top vial	100pcs
LBSV004C	4ml, 13-425 Sample Vials	4ml Clear vial, 15x45mm, screw top	100pcs
LBSV004A		4ml Amber vial, 15x45mm, screw top	100pcs
LBSV014CRS	Caps & Septa for 13-425 Sample Vials	Red PTFE/White silicone septa + Black cap with hole, for 4ml screw top vial	100pcs
LBSV104CSS		Red PTFE/White silicone septa + Black cap without hole, for 4ml screw top vial	100pcs
LBSV210C	Headspace Vials	10ml Clear vial, 22.5x46mm, crimp top, flat bottom	100pcs
LBSV220C		20ml Clear vial, 22.5x75mm, crimp top, flat bottom	100pcs
LBSV320C		20ml Clear vial, 22.5x75mm, crimp top, round bottom	100pcs
LBSV010C		10ml Clear vial, 18mm screw top, round bottom	100pcs
LBSV120C		20ml Clear vial, 18mm screw top, round bottom	100pcs
LBSV220CSS	Caps & Septa for Headspace Vials	White PTFE/White Silicone + Aluminium cap with hole, for 10ml/20ml crimp top vial	100pcs
LBSV322CRS		Red PTFE/White Silicone + Aluminium cap with hole, for 10ml/20ml crimp top vial	100pcs
LBSV120CSS		Blue PTFE/ White silicone septa + Silver screw cap with hole, 18mm screw top vial	100pcs
LBSV020C	24-400 Vials	20ml Clear vial, 24 400 screw top	100pcs
LBSV020A		20ml Amber vial, 24 400 screw top	100pcs
LBSV040C		40ml Clear vial, 24 400 screw top	100pcs
LBSV040A		40ml Amber vial, 24 400 screw top	100pcs
LBSV040CSS	Caps & Septa for 24-400 Vials	Nature PTFE/Nature silicone septa + Black screw cap without hole, 24 400 screw top vial	100pcs
LBSV140CSS		Nature PTFE/Nature silicone septa + Black screw cap with hole, 24 400 screw top vial	100pcs