


96-Well Cap Mats

Advantage™ 96-Well Cap Mats & Film

- Thin Polypropylene
- Prevents cross-contamination between wells
- Economical

Cat No.	Description	UOM
96057	1mL Autosampler Compatible Round Cap Mat, Polypropylene with 'X' on Surface	50/Pkg
1726	1mL Polypropylene Round Cap Mat	20/Pkg
2726	2mL Autosampler Compatible Square Cap Mat, Pierceable with 'X' on surface	20/Pkg
3726	2mL Square Cap Mat, low temperature (-80°), For long term storage	20/Pkg



96-Well Cap Mats

Resealable Round Cap Mats

- Compatible with all current autosampler needles
- No coring or cross-contamination between wells
- Fits all standard 96-Position Blocks
- Reseals after multiple injections
- Available pre-slit for ease of insertion and reduced vacuum formation

Catalog No.	Description	UOM
Molded Green Silicone - Good for long-term storage at -80°C		
963675	Round Well with Molded Silicone Liner	5
963670	Pre-Slit Round Well with Molded Silicone Liner	5
Ultra-Thin, Soft Clear Silicone, PTFE Coated		
964007	Ultra Thin Round Well w/ Molded Soft Clear Silicone	5
964075	Pre-Slit, Ultra Thin Round Well w/ Molded Soft Clear Silicone	5
Ultra-Thin, Purple Silicone		
965007	Ultra Thin Round Well w/ Molded Soft Purple Silicone	5
965075	Pre-Slit, Ultra Thin Round Well w/ Molded Soft Purple Silicone	5
Gray, PTFE Coated		
96FSS0	Round Well with PTFE/Silicone	5
96FPS0	Pre-Slit, Round Well with PTFE/Silicone	5



Thicker Round Cap Mats

- Thicker Silicone, PTFE Coated Cap Mats
- Recommended for Shimadzu autosamplers
- Better suited for 2° to 8°C autosampler temperatures

Cat. No.	Description	UOM
967702	Round Well with Molded Soft Clear Silicone	5
967701	Pre-Slit, Round Well with Molded Soft Clear Silicone	5



Resealable Square Cap Mats

- Compatible with all current autosampler needles
- No coring or cross-contamination between wells
- Fits all standard 96-position blocks
- Reseals after multiple injections
- Available pre-slit for ease of insertion and reduced vacuum formation

Cat. No.	Description	UOM
Molded Green Silicone		
963725	Square Well with Molded Silicone Liner	5
963724	Pre-Slit, Square Well with Molded Silicone Liner	5
Ultra-Thin, Soft Clear Silicone, PTFE Coated		
964008	Ultra Thin Square Well with Molded Soft Clear Silicone	5
964085	Pre-Slit, Ultra Thin Sq. Well w/ Molded Soft Clear Silicone	5
Ultra-Thin, Purple Soft Silicone, PTFE Coated, Mighty Mat		
965008	Purple Ultra Thin Square Well with Molded Soft Silicone	5
965085	Purple Pre-Slit, Ultra Thin Sq. Well w/ Molded Soft Silicone	5
Gray, PTFE Coated		
29660	Square Well with PTFE/Silicone	5
29665	Pre-Slit, Square Well with PTFE/Silicone	5

