

**Give us a call if you don't find
what you need in the table below**

Replacement Needles for Gas Syringes

Slip-on needles for Series A-2, C, and D syringes feature a needle extension through the hub that fits inside the PTFE tip on the syringe. This design eliminates the dead volume inherent in conventional hubbed needles.

Needles are manufactured to exacting tolerances from high precision stainless steel tube stock. The point styles and lengths listed below are standard, but we are willing to supply others. Let us know your requirements.



1. Blunt tip, open end



2. Bevel, open end



3. Side port, taper

Pressure-Lok	OD	ID	Length	Type	Pkg. of 3
Series A	.029"	x .012"	x 2.25"	2	913050
	.029"	x .012"	x 2.25"	3	913052
Series A-2	.028"	x .005"	x 2"	2	943050
	.029"	x .012"	x 2"	2	943051
	.029"	x .012"	x 2"	3	943052
Series C	.028"	x .005"	x 2"	2	943050
	.029"	x .012"	x 2"	2	943051
	.029"	x .012"	x 2"	3	943052
Series D	.029"	x .012"	x 2"	2	943051
	.029"	x .012"	x 2"	3	943052
Magnum	.035"	x .023"	x 2"	2	953050
	.035"	x .023"	x 2"	3	953052
Luer	OD	ID	Length	Type	Pkg. of 3
Series A-2	.028"	x .006"	x 2"	2	943060
	.028"	x .016"	x 2"	2	943061
	.028"	x .016"	x 2"	3	943062
Series C	.028"	x .006"	x 2"	2	943060
	.028"	x .016"	x 2"	2	943061
	.028"	x .016"	x 2"	3	943062
Series D	.028"	x .016"	x 2"	2	943061
	.028"	x .016"	x 2"	3	943062
Magnum	.028"	x .016"	x 2"	2	943061
	.028"	x .016"	x 2"	3	943062