

# Hard Drawn T304 Stainless Steel Tubing

## - Hypodermic Quality

Part No.	SIZE (mm)	METRIC O.D.	METRIC I.D.	W.T. (mm)	IMPERIAL O.D	IMPERIAL I.D.	W.T. ( " )	NEAREST SWG/OD
7001	0.4	0.40	0.22	0.09	0.016	0.009	0.0035	27
7002	0.45	0.45	0.25	0.10	0.018	0.010	0.004	26
7003	0.5	0.50	0.20	0.15	0.020	0.008	0.006	25
7004	0.6	0.60	0.30	0.15	0.032	0.020	0.006	23
7005	0.7	0.70	0.40	0.15	0.028	0.015	0.007	22
7006	0.8	0.80	0.50	0.15	0.032	.020	0.006	21
7011	1.0	1.0	0.6	0.2	0.040	0.024	0.008	19
7012	1.2	1.2	0.8	0.2	0.048	0.032	0.008	18
7013	1.4	1.4	1.0	0.2	0.056	0.040	0.008	17
7014	1.6	1.6	1.2	0.2	0.063	0.047	0.008	16
7015	1.8	1.8	1.4	0.2	0.071	0.063	0.008	15
7021	2.0	2.0	1.5	0.25	0.078	0.058	0.010	14
7022	2.2	2.2	1.7	0.25	0.087	0.067	0.010	13
7023	2.5	2.5	2.0	0.25	0.098	0.078	0.010	12
7024	2.8	2.8	2.3	0.25	0.110	0.090	0.010	11
7031	3.0	3.0	2.5	0.25	0.118	0.098	0.010	10
7032	3.5	3.5	3.0	0.25	0.138	0.118	0.010	9
7041	4.0	4.0	3.5	0.25	0.157	0.137	0.010	8
7042	4.5 x 0.25	4.5	4.0	0.25	0.178	0.157	0.010	7
7043	4.5 x 0.50	4.5	3.5	0.50	0.178	0.137	0.020	7
7051	5.0	5.0	4.5	0.25	0.198	0.178	0.010	6
7061	6.0	6.0	5.5	0.25	0.237	0.217	0.010	4
7075	8.0	8.0	7.5	0.25	0.316	0.266	0.010	2

<b>Soft Annealed T316 Stainless Steel Capillary Tubing</b>							
Part No.	SIZE (mm)	METRIC O.D.	METRIC I.D.	W.T. (mm)	IMPERIAL O.D	IMPERIAL I.D.	W.T. ( " )
7080	1/8" x 30g	3.18	2.54	0.33	0.125	0.099	0.013
7081	1/8" x 27g	3.18	2.36	0.41	0.125	0.093	0.016
7082	1.8" x 22g	3.18	1.74	0.72	0.125	0.069	0.028
7083	1/8" x 20g	3.18	1.38	0.90	0.125	0.054	0.035

## Stainless Steel Hypo Tubing Price list updated **2019+ per metre**

Part No.	SIZE (mm)	GRADE S.S.	LENGTH		
			0.9M LONG	1.5M LONG	2.0M LONG
7001	0.4	T304			\$ 36.00
7002	0.45	T304		\$25.00	
7003	0.5	T304	\$15.00		
7004	0.6	T304		\$ 4000	
7005	0.7	T304			\$ 40.00
7006	0.8	T304			\$ 42.00
7011	1.0	T304			\$ 65.00
7012	1.2	T304			\$ 68.00
7013	1.4	T304			\$ 50.00
7014	1.6	T304			\$62.00
7015	1.8	T304			\$68.00
7021	2.0	T304			\$68.00
7022	2.2	T304			\$72.00
7023	2.5	T304			\$51.00
7024	2.8	T304			\$75.. 00
7031	3.0	T304			\$78.00
7032	3.5	T304			\$78.00
7041	4.0	T304			\$85.00
7042	4.5 X.25	T304	Thin wall		\$85.00
7043	4.5X .50	T304	Heavy Wall		\$88.00
7051	5.0	T304			\$98.00
7061	6.0	T304			\$98.00
7075	8.0	T304			\$98.00

<b>Soft Annealed T316 Stainless Steel Coiled Capillary Tubing</b>		
SIZE (mm)	PRICING APPLIES TO ALL SIZES.	
7080	1/8" x 30g	FULL 30M COILS ; \$60.00 per Metre, or \$1500.00 per Coil
7081	1.8" x 27g	CUT LENGTHS ; \$ 60.00 per Metre.
7082	1.8" x 22g	
7083	1.8" x 20g	

Annealed(Soft) SS #304 and #316 in 1/4OD(x3/16), 1/8(X2.1MMID) 1/16"OD (1.0, 0.75, 0.5, 0.25, 0.18, 0.13mmID)  
 Any Length > > 200 foot : Unique Silcosteel/Sulfinert : . . . Polymeric > Teflon, LDPE etc various sizes **ASK !**  
 see Flip Catalog <http://flip.chromalytic.net.au/books/ptit/>

**MINIMUM INVOICE VALUE IS \$50.00 EXCLUDING 10% GST AND FREIGHT**

	<b>+61(0)3 9762 2034</b>	<b>Australian Distributors</b> Importers & Manufacturers <a href="http://www.chromtech.net.au">www.chromtech.net.au</a>
<b>ECH</b> nology Pty Ltd Website NEW : <a href="http://www.chromalytic.net.au">www.chromalytic.net.au</a> E-mail : <a href="mailto:info@chromtech.net.au">info@chromtech.net.au</a> Tel: 03 9762 2034 . . . in AUSTRALIA		