



Sulfinert® Treated Ultra-High Pressure Sample Cylinders and Valves

- Stable storage of samples containing sulfur compounds and mercury.
- Cylinders manufactured by Swagelok® and U.S. D.O.T. rated to 5000psig.
- Valves rated to 6000psig.
- 316 grade stainless steel with 1/4" female NPT threads on both ends.

Ideal for collecting samples at gas wellhead or other applications requiring sampling at extremely high pressures.

Description	Fittings	Volume	qty.	Cat.#
Sulfinert Sample Cylinder, 1/4" Female NPT Thread		150cc	ea.	22111
Sulfinert Sample Cylinder, 1/4" Female NPT Thread		300cc	ea.	22112
Sulfinert Sample Cylinder, 1/4" Female NPT Thread		500cc	ea.	22113
Sulfinert Sample Cylinder Valve	1/4" Male NPT (both ends)		ea.	22109
Sulfinert Sample Cylinder Valve	1/4" Male NPT x 1/4" Female NPT		ea.	22110

HROMalytic Chromatography
Australian Distributors ECHnology Products '08

www.chromtech.net.au E-mail : info@chromtech.net.au Tel : +61 3 9762 2034 Fax : +61 3 9761 1169

