

BPCL220 (-10 -81)
Blood Pressure Adapter with 5/16-24 UNF Thread and 6% Taper to 200 Series Barb, 3/32” (2.4 mm) ID Tubing



BPC220 (-10)
Blood Pressure Adapter with 5/16-24 UNF Thread to 200 Series Barb, 3/32” (2.4 mm) ID Tubing



BPCL230 (-10 -81)
Blood Pressure Adapter with 5/16-24 UNF Thread and 6% Internal Taper to 200 Series Barb, 1/8” (3.2 mm) ID Tubing



BPC230 (-10)
Blood Pressure Adapter with 5/16-24 UNF Thread to 200 Series Barb, 1/8” (3.2 mm) ID Tubing



BPCL025 (-10)
Blood Pressure Adapter with 5/16-24 UNF Thread 6% Internal Taper to 500 Series Barb, 5/32” (4.0 mm) ID Tubing



BPC025 (-10)
Blood Pressure Adapter with 5/16-24 UNF Thread to 500 Series Barb, 5/32” (4.0 mm) ID Tubing



BPCL240 (-10 -81)
Blood Pressure Adapter with 5/16-24 UNF Thread and 6% Internal Taper to 200 Series Barb, 5/32” (4.0 mm) ID Tubing



BPC240 (-10)
Blood Pressure Adapter with 5/16-24 UNF Thread to 200 Series Barb, 5/32” (4.0 mm) ID Tubing



BPCL250 (-10 -81)
Blood Pressure Adapter with 5/16-24 UNF Thread and 6% Internal Taper to 200 Series Barb, 3/16” (4.8 mm) ID Tubing



BPC250 (-10)
Blood Pressure Adapter with 5/16-24 UNF Thread to 200 Series Barb, 3/16” (4.8 mm) ID Tubing



CF229 (-1)
Neonatal Blood Pressure Slip Fit (non luer) with 200 Series Barb, 3/32” (2.4 mm) ID Tubing



CF230 (-1)
Neonatal Blood Pressure Slip Fit (non luer) with 200 Series Barb, 5/32” (4.0 mm) ID Tubing



Did you know?

Value Plastics has an Interactive PDF available online at www.valueplastics.com. Access part drawings, reference charts, data sheets, and more from just one location!

BLOOD
PRESSURE

- 001 (various, see part description)
- 06 Orange Nylon
- 07 Yellow Nylon
- 0 Natural Nylon
- 1 White Nylon
- 2 Black Nylon
- 3 Red Nylon
- 4 Green Nylon
- 5 Blue Nylon
- 7 Glass Filled Natural Nylon
- 8 Gray Nylon
- 8B Gray Nylon
- 9 Clear Polycarbonate
- 10 Acetal
- 40 Polysulfone
- 81 Lace White ABS
- 91 White Polycarbonate
- 1006 Natural Acetal
- 1007 Black Acetal
- 6005 Animal-Free Polypropylene
- 6006 White Polypropylene
- 8003 Natural ABS
- 8012 White ABS
- 9002 Rad. Stable Polycarbonate
- 9010 Rad. Stable Polycarbonate
- B Buna-N O-ring
- CM030 Dark Blue Nylon
- E PTEE/EPDM
- EP EPDM O-ring
- J1A Natural Kynar PVDF
- J1A2 Black Kynar PVDF
- PX Platinum-Cured Silicone
- U71 White Class VI Glass Filled Nylon
- U75 Blue Class VI Glass Filled Nylon
- VP1 Animal-free Polyethylene
- X2 Buna-N Rubber
- XO Stainless Steel

Resin Reference

- Pg 82 Chemical Resistance
- Pg 83 Regulatory Status
- Pg 84 Sterilization Stability
- Pg 85 Typical Properties