



**VELOX**<sup>™</sup>  
Solvent Waste System

[www.zodiaclifesciences.com](http://www.zodiaclifesciences.com)



**HROMalytic** +61(0)3 9762 2034  
**ECH**nology Pty Ltd

**Australian Distributors**  
Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)


Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

# Velox Solvent waste systems

## Introduction:

Fluids in HPLC waste containers are overwhelmingly dangerous containing potent volatile organic compounds (VOC) posing potent risk to health of Lab workers & lab safety.

Zodiac came with an innovative solution of introducing **VELOX Solvent waste system** with activated carbon filters not only creating an organised containment for HPLC solvent waste but also providing an efficient way for collection and disposal of volatile vapours.



## Advantages

- Filters system against harmful vapours
- Capillary fixations with PFA Fittings
- Several systems can be connected to 1 container
- Different size containers compatible to any space condition

[www.zodiaclifesciences.com](http://www.zodiaclifesciences.com)







90 Adaptor included  
for maximum  
storage flexibility



Normal cap is additionally  
provided for convenient  
transport of the wastes



Blind fittings are provided  
to block the remaining  
ports when not in use

*Zodiac Velox Solvent Waste System allows 99 percent of the volatile vapours to be trapped into its Powerful exhaust filters keeping the laboratory environment safe and clean , This system is more organised than traditional techniques protects against spills and leaks. This is a necessary addition to every lab for improved solvent waste handling and disposal*

## ORDERING INFORMATION

| Item Description                    | Part Number  |
|-------------------------------------|--------------|
| Velox Solvent Waste System - 5Ltrs  | ZLS.332.2419 |
| Velox Solvent Waste System - 10Ltrs | ZLS.332-3413 |
| Velox Solvent Waste System - 20Ltrs | ZLS.332-4413 |

[www.zodiaclifesciences.com](http://www.zodiaclifesciences.com)

