

Verder QCUT 150M Metal Cutter pdf

www.chromtech.net.au/pdf2/Wet-Abrasive-Cut-Off-Machinee28094QCUT-150-A.pdf

Video QCUT 150A

<https://www.youtube.com/watch?v=W2kC5Vop-LE>

my email to > 19Nov22

To Verder (Germany) QATM.com

Price on QCUT 150M or Custom cut-off machine

Specifically for cutting SS capillary tubing 1/4,1/8,1/16,(1/32inch OD; various IDs'

Fine grit > #500 grit

Ultrathin blade 0.4-0.5mm with spare blades readily available ! ~1-2 metre MAX length

Water Cooling essential

WE purchased elsewhere years ago but need a replacement machine specifically for this type of job

Cut-off but also for cross-slotting of appropriate size tubing and section of tube along short length

~10mm for ID visualisation.

Custom machine but maybe you have something already similar to the QCUT150

50-75mm Diam blade Is OK