

“Pharma” Products / Vial Filling - Liquid : re Production Upgrade : Chromalytic Technology Pty Ltd

Price in AUD ex GST varies as to XR USD/AUD at date of Ordering and payment

Equipment : **Payment** >\$2000 75% up front at date of ordering


ex Freight depends on shipping method

ASK for a QUOTE or details /Spec Sheets (some demo videos are) available

Stage 1 Manual Filling with Syringe Manual Insert Stopper & Cap Crimp with Hand Crimper	Stage 2 Fill with BottleTop Dispenser BTD-3 (10ml) \$195 other size (2.5, 10, ,25, 50ml)	Stage 3 Fill with Peristaltic Pump NE-9000 (programmable or via foot pedal) Manual move head from vial to vial in a rack (50each) \$695.00 PP Rack 50ea \$45	Stage 4 Auto Fill vials Robotic XY Table & Tank 50Kg bench Top \$11-12k Batches of 300/500ea manual load into rack Fills from 2l of 5liter SS tank per hour ??? 2 persons Availability 5 to 6 weeks by Air	Stage 5 Auto feed of vials, auto fill of vials; semi manual insert of caps & stoppers semi auto Crimping BenchTop 60kg \$120-150K max 1800ea/hour one person Availability ???	Stage 6 Complete vial filling / auto feed of vials Auto insertion of stoppers & caps AutoCrimping estimated 400KgV; 180Kg Shipping/ \$20-25k max 2000ea/hour Availability; 3/3 weeks By Sea Freight FIRM more details ~\$1250 for second vial size mods
Stoppers & Caps 13mm / 20mm Chlorobutyl Rubber	Manual	Manual with Vial Racks (PP 50ea \$45)			
Crimping All size heads –S / -FO	Manual Hand Crimper New SS Hand Crimper Sterilisable	XL Bench Crimper Vertical Hand Air Crimper (in stock) Bench Air Crimper Counter-Balanced Air Crimper (large vials / Heavy Duty)			
Filtration (Membrane Solutions) Syringe Filter (in stock) Many Types 0.22/0.45um 25mmD MultiLayer prefilters GC/X Whatman	1 VacUFill (sterile plastic) 500/500ml PES/ PVDF 0.2 and 0.45um 2 Glass Filtration 47mm Diam Disks - Glass Fibre preFilters Hand Vac Pump Motorised Vac Pumps	VacUCap PES 0.,2/0.45um - high(er) Speed Filtration 90mm Diam	Filter Cartridges (ASK)	Syringe Pump Syringe Heater Kit / Syringe Filters (Viscous Solutions)	
Sterilising / Autoclave in-stock	Now available 24liter TableTop Auto-clave Auto Steam/Steriliser “Pro grade” CE Auto steam/vacuum cycled \$3500 60Kg				
Ampoules : 1.2,5, 10ml Clear/amber Ampoule Filling	Manual Liquid Fill (1ml or 2,5,10ml) Auto Flame seal fro a Turntable (Manual) fill 24 at time) \$10k BenchTop 10Kg	AutoFill of Vial (N2 Purge) Auto Flame seal (Stage 1) \$25-30K BenchTop 20Kg			
Accessory : Snap-It tool					
Spray Caps to Order Only (some stock)	Manual / Heavy Duty CrimpPress				
			CT Popular “pharma” Vials	2ml Clear CTSC-2V-13T-15D32H	371 per slab
				4ml Clear CTSC-4V-13T-16D38H	368 per slab
				10ml Clear CTSC-10V-20T-22D45H	324 per slab
				13ml Clear CTSC-13V-20T-24D49H	145 per slab
				7ml Amber CTSC-7V-20T-22D38H	150 per slab

CT Popular “pharma” Vials	2ml Clear CTSC-2V-13T-15D32H	371 per slab
	4ml Clear CTSC-4V-13T-16D38H	368 per slab
	10ml Clear CTSC-10V-20T-22D45H	324 per slab
	13ml Clear CTSC-13V-20T-24D49H	145 per slab
	7ml Amber CTSC-7V-20T-22D38H	150 per slab


“Pharma” Products / Vial Filling - Liquid : re Production Upgrade : Chromalytic Technology Pty Ltd




icanCLAVE
 "Icanclave" autoclaves sterilization



“Pharma” Products / Vial Filling - Liquid : *re Production Upgrade* : Chromalytic Technology Pty Ltd





Syringe Filters



NE-300 Syringe Pump



Pump Foot



Syringe Filtration / Heater Kit

Vacufil™ Disposable Vacuum Filtration




Vacuum Pumps



CT-TD20W



VPJ0332



Vacuum Filtration



Vac-U-Caps

