

Distilled Water

- SS Semi-Automatic / Glass receiver

Tapwater : Melbourne : Conductivity 80megohms; Sydney Adelaide >200megohms
after Water Still <2.0megohms . . . slowly increases as water absorbs CO2 from ambient Air (& harmless !)

Wash Lab Glass ware; PC screens; windows & dry spotless !

"even good enough for your car"

GOOD ENOUGH TO DRINK . . . even ! . . . and what's more beneficial to you !

1 to 2 litres / hour . . . only needs manual topping up with water . . . auto turn-off !

Highlights Of The Countertop Water Distiller



Countertop Water Distiller

- Includes 1 x 4 litre **glass Collector Bottle / Jug**
- Build / componentry quality - **excellent / ultra**

- Stainless Steel Under-Lid** - extra thick for longevity
- Stainless Steel Boiling Chamber** - extra thick for longevity
- Stainless Steel Condensing Coil** - not aluminium
- Non rust poly-carbonate outer casing
- Quiet operation / fan cooled
- Power **on** indicator / light
- Easy to use & clean
- Power consumption only **6-12¢** per litre

- Water purification by steam distillation & carbon filtration
- Semi automatic** - i.e. manual fill - **auto shut off** after each cycle
- Provides enough pure drinking water for a family of **1-5 people**
- Provides up to **24 litres** in **24 hours!**
- Produces approximately **4 litres** in each **4 hour cycle**
- Removes in excess of **99+%** of water impurities
- Almost no water loss during process - saves precious water
- Produces clean water for sparkling ice cubes
- Enhances the natural flavour of juices, soups, coffee, tea and other beverages. Oh, and it **Tastes great!**
- More**

"Very popular ~ very robust ~ ultra reliable!"

Easy to use semi-automatic Water Distillation unit!

HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA