

MS Syringe Filters Product List 2014

(Ctrl F)
> Cat.#s

> Current System AUD Prices



SFPVDF033045S	Sterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:33(mm)	
SFPVDF013022NP	Nonsterile PVDF Syringe Filter, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPVDF013045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFPVDF025022NP	Nonsterile PVDF Syringe Filter, Pore: 0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPVDF025045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(µm), Diameter: 25(mm), PP Prefilter	100
SFPVDF050022NP	Nonsterile PVDF Syringe Filter, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFPVDF050045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:50(mm), PP Prefilter	12
SFPVDF013022SP	Sterile PVDF Syringe Filter, Pore: 0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPVDF013045SP	Sterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFPVDF025022SP	Sterile PVDF Syringe Filter, Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPVDF025045SP	Sterile PVDF Syringe Filter, Pore: 0.45(µm), Diameter:25(mm), PP Prefilter	100
PES Syringe Filters : : ALL Sizes		
SFPES013002NG	Nonsterile PES Syringe Filters, Pore: 0.02(µm), Diamter:13(mm), GF Prefilter	100
SFPES013022NG	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF prefilter	100
SFPES013045NG	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter: 13(mm), GF Prefilter	100
SFPES025002NG	Nonsterile PES Syringe Filters, Pore:0.02(µm), Diameter:25(mm), GF Prefilter	100
SFPES025022NG	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm), GF Prefilter	100
SFPES025045NG	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFPES050022NG	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:50(mm), GF Prefilter	12
SFPES050045NG	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:50(mm), GF prefilter	12
SFPES013022SG	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFPES013045SG	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFPES025022SG	Sterile PES Syringe Filters, Pore:0.02(µm), Diameter:25(mm), GF Prefilter	100
SFPES025045SG	Sterile PES Syringe Filters, Pore:0.02(µm), Diameter:25(mm), GF Prefilter	100
SFPES013010N	Nonsterile PES Syringe Filters, Pore:0.1(µm), Diameter:13(mm)	500
SFPES013022N	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm)	500
SFPES013045N	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:13(mm)	500
SFPES017022N	Nonsterile PES Syringe Filter, Pore:0.22(µm), Diameter:17(mm)	500
SFPES017045N	Nonsterile PES Syringe Filter, Pore:0.45(µm), Diameter:17(mm)	500
SFPES025010N	Nonsterile PES Syringe Filters, Pore:0.1(µm), Diameter:25(mm)	200
SFPES025022N	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm)	200
SFPES025045N	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm)	200
SFPES02500N	Nonsterile PES Syringe Filters, Pore:2.0(µm), Diameter:25(mm)	200
SFPES030022N	Nonsterile PES Syringe Filter, Pore:0.22(µm), Diameter:30(mm)	200
SFPES030045N	Nonsterile PES Syringe Filter, Pore:0.45(µm), Diameter:30(mm)	200
SFPES030200N	Nonsterile PES Syringe Filter, Pore:2.0(µm), Diameter:30(mm)	200
SFPES033022N	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:33(mm)	200
SFPES033045N	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:33(mm)	200
SFPES050022N	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:50(mm)	12
SFPES050045N	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:50(mm)	12
SFPES013022S	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm)	200
SFPES013045S	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:13(mm)	200
SFPES025010S	Sterile PES Syringe Filters, Pore:0.1(µm), Diameter:25(mm)	100
SFPES025022S	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm)	200
SFPES025045S	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm)	100
SFPES030022S	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:30(mm)	100
SFPES033022S	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:33(mm)	200
SFPES033045S	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:33(mm)	200
SFPES013022NP	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100