

MS Syringe Filters Product List 2014

(Ctrl F)
> Cat.#s



> Current System AUD Prices

SFPE013045NP	Nonsterile PES Syringe Filters,Pore: 0.45(µm), Diameter:13(mm), PP Prefilter	
SFPE025022NP	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPE025045NP	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP prefilter	100
SFPE050022NP	Nonsterile PES Syringe Filters, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFPE050045NP	Nonsterile PES Syringe Filters, Pore:0.45(µm), Diameter:50(mm), PP Prefilter	12
SFPE013022SP	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPE013045SP	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP prefilter	100
SFPE025022SP	Sterile PES Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP prefilter	100
SFPE025045SP	Sterile PES Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP prefilter	100

Nylon Syringe Filters : : ALL Sizes

SFNY013010NG	Nonsterile Nylon Syringe Filters, Pore:0.10(µm), Diameter:13(mm), GF Prefilter	100
SFNY013022NG	Nonsterile Nylon Syringe Filters,Pore: 0.22(µm), Diameter: 13(mm), GF Prefilter	100
SFNY013045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFNY025022NG	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm), GF Prefilter	100
SFNY025045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFNY025500NG	Nonsterile Nylon Syringe Filters, Pore:5.0(µm), Diameter:25(mm), GF Prefilter	100
SFNY030045NG	Nonsterile Nylon Syringe Filter, Pore:0.45(µm), Diameter:30(mm),GF Prefilter	100
SFNY050022NG	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:50(mm), GF Prefilter	12
SFNY050045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:50(mm), GF Prefilter	12
SFNY013022SG	Sterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFNY013045SG	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFNY025022SG	Sterile Nylon Syringe Filters, Pore:0.22(µm),Diameter: 25(mm), GF Prefilter	100
SFNY025045SG	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SF1322N1000	Non-sterile Nylon Syringe Filter, Pore:0.22(µm),Diameter:13(mm)	1000
SF1345N1000	Non-sterile Nylon Syringe Filter, Pore:0.45(µm),Diameter:13(mm)	1000
SFNY004020N	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:4(mm)	100
SFNY013010N	Nonsterile Nylon Syringe Filters, Pore:0.1(µm), Diameter:13(mm)	500
SFNY013022N	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:13(mm)	500
SFNY013045N	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm)	500
SFNY013080N	Nonsterile Nylon Syringe Filters, Pore:0.8(µm), Diameter:13(mm)	500
SFNY013100N	Nonsterile Nylon Syringe Filters, Pore:1.0(µm), Diameter:13(mm)	500
SFNY013300N	Nonsterile Nylon Syringe Filters, Pore:3.0(µm), Diameter:13(mm)	500
SFNY013500N	Nonsterile Nylon Syringe Filters, Pore:5.0(µm), Diameter:13(mm)	500
SFNY017022N	Nonsterile Nylon Syringe Filter, Pore:0.22(µm), Diameter:17(mm)	500
SFNY017045N	Nonsterile Nylon Syringe Filter, Pore:0.45(µm), Diameter:17(mm)	500
SFNY025022N	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm)	200
SFNY025045N	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:25(mm)	200
SFNY025080N	Nonsterile Nylon Syringe Filters, Pore:0.8(µm), Diameter:25(mm)	200
SFNY025100N	Nonsterile Nylon Syringe Filters, Pore:1.0(µm), Diameter:25(mm)	200
SFNY025200N	Non-sterile Nylon Syringe Filters, Pore:2.0(µm), Diameter:25(mm)	200
SFNY025300N	Nonsterile Nylon Syringe Filters, Pore:3.0(µm), Diameter:25(mm)	200
SFNY025500N	Nonsterile Nylon Syringe Filters, Pore:5.0(µm), Diameter:25(mm)	200
SFNY030022N	Nonsterile Nylon Syringe Filter, Pore:0.22(µm), Diameter:30(mm)	200
SFNY030045N	Nonsterile Nylon Syringe Filter, Pore:0.45(µm), Diameter:30(mm)	200

06/2014



HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA