

MS SyringeFilters Product List 2014

(Ctrl F)
> Cat.#s



> Current System AUD Prices

SFNY033022N	Nonsterile Nylon Syringe Filters, Pore:0.22(µm),Diameter:33(mm)	
SFNY033045N	Nonsterile Nylon Syringe Filters, Pore:0.45(µm),Diameter:33(mm)	200
SFNY050022N	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:50(mm)	12
SFNY050045N	Nonsterile Nylon Syringe Filters, Pore:0.45µm, Diameter:50(mm)	12
SFNY013010S	Sterile Nylon Syringe Filters, Pore:0.1(µm), Diameter:13(mm)	200
SFNY013022S	Sterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:13(mm)	200
SFNY013045S	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm)	200
SFNY013080S	Sterile Nylon Syringe Filters, Pore:0.8(µm), Diameter:13(mm)	200
SFNY013100S	Sterile Nylon Syringe Filters, Pore:1.0(µm), Diameter:13(mm)	200
SFNY013300S	Sterile Nylon Syringe Filters, Pore:3.0(µm), Diameter:13(mm)	200
SFNY013500S	Sterile Nylon Syringe Filters, Pore:5.0(µm), Diameter:13(mm)	200
SFNY025010S	Sterile Nylon Syringe Filters, Pore:0.1(µm), Diameter:25(mm)	200
SFNY025022S	Sterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm)	200
SFNY025045S	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:25(mm)	200
SFNY025080S	Sterile Nylon Syringe Filters, Pore:0.8(µm), Diameter:25(mm)	200
SFNY025100S	Sterile Nylon Syringe Filters, Pore:1.0(µm), Diameter:25(mm)	200
SFNY025200S	Sterile Nylon Syringe Filters, Pore:2.0(µm), Diameter:25(mm)	200
SFNY025300S	Sterile Nylon Syringe Filters, Pore:3.0(µm), Diameter:25(mm)	200
SFNY025500S	Sterile Nylon Syringe Filters, Pore:5.0(µm), Diameter:25(mm)	200
SFNY033022S	Sterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:33(mm)	200
SFNY033045S	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:33(mm)	200
SFNY013022NP	Nonsterile Nylon Syringe Filters,Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFNY013045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFNY025022NP	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFNY025045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), PP Prefilter	100
SFNY050022NP	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFNY050045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:50(mm), PP Prefilter	12
SFNY013022SP	Sterile Nylon Syringe Filters, Pore:0.22(µm),Diameter: 13(mm), PP Prefilter	100
SFNY013045SP	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFNY025022SP	Sterile Nylon Syringe Filters,Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFNY025045SP	Sterile Nylon Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), PP Prefilter	100
Pre Filter Syringe Filters : ALL Sizes		
SFCA025045NG	Nonsterile CA Syringe Filter, Pore:0.45(µm), Diameter:25(mm), GF prefilter	200
SFNY013010NG	Nonsterile Nylon Syringe Filters, Pore:0.10(µm), Diameter:13(mm), GF Prefilter	100
SFNY013022NG	Nonsterile Nylon Syringe Filters,Pore: 0.22(µm), Diameter: 13(mm), GF Prefilter	100
SFNY013022NP	Nonsterile Nylon Syringe Filters,Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFNY013022SG	Sterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFNY013022SP	Sterile Nylon Syringe Filters, Pore:0.22(µm),Diameter: 13(mm), PP Prefilter	100
SFNY013045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFNY013045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFNY013045SG	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFNY013045SP	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFNY025022NG	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm), GF Prefilter	100
SFNY025022NP	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100

SEARCH [cat#] Write Down Cat#

6

where "SUPPLIER" = restek, vicivalco, vicijour, striges, valueplastics (MAIN Suppliers ONLY > Latest Prices

7 [Ctrl "F"] to SEARCH for Item > Description/Cat# etc

06/2014



HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Go To www.chromtech.net.au/fliip/catalogs/"SUPPLIER" THEN add "/index.php" to URL