

MS Syringe Filters Product List 2014

(Ctrl F)
> Cat.#s



> Current System AUD Prices

SFPTFE025045NP	Nonsterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP Prefilter	
SFPTFE025045SG	Sterile PTFE Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), GF Prefilter	100
SFPTFE025045SP	Sterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:25(mm), PP Prefilter	100
SFPTFE050022NG	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:50(mm), GF Prefilter	12
SFPTFE050022NP	Nonsterile PTFE Syringe Filters, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFPTFE050045NG	Nonsterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:50(mm), GF prefilter	12
SFPTFE050045NP	Nonsterile PTFE Syringe Filters, Pore:0.45(µm), Diameter:50(mm), PP Prefilter	12
SFPVDF013022NG	Nonsterile PVDF Syringe Filter, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFPVDF013022NP	Nonsterile PVDF Syringe Filter,Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPVDF013022SG	Sterile PVDF Syringe Filter,Pore: 0.22(µm), Diameter:13(mm), GF Prefilter	100
SFPVDF013022SP	Sterile PVDF Syringe Filter,Pore: 0.22(µm), Diameter:13(mm), PP Prefilter	100
SFPVDF013045NG	Nonsterile PVDF Syringe Filter,Pore: 0.45(µm),Diameter: 13(mm), GF Prefilter	100
SFPVDF013045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFPVDF013045SG	Sterile PVDF Syringe Filters, Pore:0.45(µm),Diameter: 13(mm), GF Prefilter	100
SFPVDF013045SP	Sterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFPVDF025022NG	Nonsterile PVDF Syringe Filter, Pore:0.22(µm),Diameter: 25(mm), GF Prefilter	100
SFPVDF025022NP	Nonsterile PVDF Syringe Filter,Pore: 0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPVDF025022SG	Sterile PVDF Syringe Filter,Pore: 0.22(µm),Diameter: 25(mm), GF Prefilter	100
SFPVDF025022SP	Sterile PVDF Syringe Filter,Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFPVDF025045NG	Nonsterile PVDF Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFPVDF025045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(µm),Diameter: 25(mm), PP Prefilter	100
SFPVDF025045SG	Sterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFPVDF025045SP	Sterile PVDF Syringe Filter,Pore: 0.45(µm), Diameter:25(mm), PP Prefilter	100
SFPVDF050022NG	Nonsterile PVDF Syringe Filter, Pore:0.22(µm), Diameter:50(mm), GF Prefilter	12
SFPVDF050022NP	Nonsterile PVDF Syringe Filter, Pore:0.22(µm), Diameter:50(mm), PP Prefilter	12
SFPVDF050045NG	Nonsterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:50(mm), GF Prefilter	12
SFPVDF050045NP	Nonsterile PVDF Syringe Filter, Pore:0.45(µm), Diameter:50(mm), PP Prefilter	12
SFCA025045NG	Nonsterile CA Syringe Filter, Pore:0.45(µm), Diameter:25(mm), GF prefilter	200
SFNY013010NG	Nonsterile Nylon Syringe Filters, Pore:0.10(µm), Diameter:13(mm), GF Prefilter	100
SFNY013022NG	Nonsterile Nylon Syringe Filters,Pore: 0.22(µm), Diameter: 13(mm), GF Prefilter	100
SFNY013022NP	Nonsterile Nylon Syringe Filters,Pore:0.22(µm), Diameter:13(mm), PP Prefilter	100
SFNY013022SG	Sterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:13(mm), GF Prefilter	100
SFNY013022SP	Sterile Nylon Syringe Filters, Pore:0.22(µm),Diameter: 13(mm), PP Prefilter	100
SFNY013045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFNY013045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFNY013045SG	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), GF Prefilter	100
SFNY013045SP	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:13(mm), PP Prefilter	100
SFNY025022NG	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm), GF Prefilter	100
SFNY025022NP	Nonsterile Nylon Syringe Filters, Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFNY025022SG	Sterile Nylon Syringe Filters, Pore:0.22(µm),Diameter: 25(mm), GF Prefilter	100
SFNY025022SP	Sterile Nylon Syringe Filters,Pore:0.22(µm), Diameter:25(mm), PP Prefilter	100
SFNY025045NG	Nonsterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFNY025045NP	Nonsterile Nylon Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), PP Prefilter	100
SFNY025045SG	Sterile Nylon Syringe Filters, Pore:0.45(µm), Diameter:25(mm), GF Prefilter	100
SFNY025045SP	Sterile Nylon Syringe Filters, Pore:0.45(µm),Diameter: 25(mm), PP Prefilter	100

where "SUPPLIER" = restek, vicivalco, vicijour, striges, valueplastics (MAIN Suppliers ONLY > Latest Prices

SEARCH [cat#] Write Down Cat#

06/2014



HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA