

MS Plastic Vacuum Filter System ProductList 2

[Ctrl F]
 > Cat.#s

> Current System AUD Prices



Filter System (Membrane Diameter : 90mm, Funnel Capacity : 500ml, Receiver Bottle Capacity:1000ml)				
VFPPEs2221000	PES	0.22	12	52.5*38.5*29.5
VFPPEs2451000		0.45	12	
VFPMCE2221000	MCE	0.22	12	
VFPMCE2451000		0.45	12	
VFPCA2221000	CA	0.22	12	
VFPCA2451000		0.45	12	
VFPNY2221000	Nylon	0.22	12	
VFPNY2451000		0.45	12	
VFPVPDF2221000	PVDF	0.22	12	
VFPVPDF2451000		0.45	12	
Cat.No.	Membrane Material	Filter Pore Size(μm)	Qty (pcs/cs)	Case Size(cm)
Funnel Capacity : 250ml, Membrane Diameter : 50mm				
VFPPEs122250F	PES	0.22	24	33*33*30
VFPPEs145250F		0.45	24	
VFPMCE122250F	MCE	0.22	24	
VFPMCE145250F		0.45	24	
VFPCA122250F	CA	0.22	24	
VFPCA145250F		0.45	24	
VFPNY122250F	Nylon	0.22	24	
VFPNY145250F		0.45	24	
VFPVPDF122250F	PVDF	0.22	24	
VFPVPDF145250F		0.45	24	
Funnel Capacity : 500ml, Membrane Diameter : 90mm				
VFPPEs222500F	PES	0.22	24	52.5*38.5*29.5
VFPPEs245500F		0.45	24	
VFPMCE222500F	MCE	0.22	24	
VFPMCE245500F		0.45	24	
VFPCA222500F	CA	0.22	24	
VFPCA245500F		0.45	24	
VFPNY222500F	Nylon	0.22	24	
VFPNY245500F		0.45	24	
VFPVPDF222500F	PVDF	0.22	24	
VFPVPDF245500F		0.45	24	
Cat.No.	Description		Qty(pcs/cs)	Case Size(cm)
Receiver Bottle				
VFP250B	250ml polystyrene receiver bottle, sterile		24	33*33*30
VFP500B	500ml polystyrene receiver bottle, sterile		24	52.5*38.5*29.5
VFP1000B	1000ml polystyrene receiver bottle, sterile		24	52.5*38.5*29.5