



Note: 1) All of the price listed below are **Ex work price** below are Ex work price.

2) All the price listed below are for Membrane Solutions' Brand label.

3) **Order Qty** means **total items** in one order.

4) Membrane Solutions reserves itself the right to modify prices and content of this document without previous notice.

	Item#	Description	Packaging
<b>Glass Vacuum Filter</b>	VFG030010	1) 300ml Glass Funnel with cover	1 set
		2) 1000ml Glass Solvent Collection	
		3) Stainless Steel Clamp	
		4) Solid Glass Frit	
		Packaging: 30cmX30cmX25cm	
<b>Glass Vacuum Filter</b>	VFG030020	1) 300ml Glass Funnel with cover	1 set
		2) 2000ml Glass Solvent Collection	
		3) Stainless Steel Clamp	
		4) Solid Glass Frit	
		Packaging: 30cmX30cmX25cm	
<b>Glass Vacuum Filter</b>	VFG050020	1) 500ml Glass Funnel with cover	1 set
		2) 2000ml Glass Solvent Collection	
		3) Stainless Steel Clamp	
		4) Solid Glass Frit	
		Packaging: 30cmX30cmX25cm	
<b>Support Screen</b>	VFG147SS	Pyrex glass support screen (1000ML)	1pcs
	VFG047SS	Pyrex glass support screen (2000ML)	1pcs
	VFG010RB	Solvent Collection Bottle, 1000ml	1pcs
	VFG020RB	Solvent Collection Bottle, 2000ml	1pcs