

# MS Glass Chromatography Vials Product

[Ctrl F]  
> Cat.#s

> Current System AUD Prices



MS Items	Type	Products Name/Size	Pcs/pk
LBSV012C	<b>8-425 &amp; 9-425 2ml Clear Vial</b>	2ml Clear vial, 8-425 screw top	100pcs
LBSV032C		2ml Clear vial, 8-425 screw top, graduated with writing area	100pcs
LBSV002C	<b>12x32mm</b>	2ml Clear vial, 9-425 screw top	100pcs
LBSV022C		2ml Clear vial, 9-425 screw top, graduated with writing area	100pcs
LBSV042C		2ml Clear vial, crimp top	100pcs
LBSV062C		2ml Clear vial, crimp top, graduated with writing area	100pcs
LBSV052C		2ml Clear vial, snap top	100pcs
LBSV072C		2ml Clear vial, snap top, graduated with writing area	100pcs
LBSV012A	<b>8-425 &amp; 9-425 2ml Amber Vial</b>	2ml Amber vial, 8-425 screw top	100pcs
LBSV032A		2ml Amber vial, 8-425 screw top, graduated with writing area	100pcs
LBSV002A	<b>12x32mm</b>	2ml Amber vial, 9-425 screw top	100pcs
LBSV022A		2ml Amber vial, 9-425 screw top, graduated with writing area	100pcs
LBSV042A		2ml Amber vial, crimp top	100pcs
LBSV062A		2ml Amber vial, crimp top, graduated with writing area	100pcs
LBSV052A		2ml Amber vial, snap top	100pcs
LBSV072A		2ml Amber vial, snap top, graduated with writing area	100pcs
LBSV102CSS	<b>Cap &amp; Septa for 2ml Sample Vials</b>	Red PTFE/white silicone septa + Blue screw cap with hole, for 2ml 9-425 screw top vial	100pcs
LBSV002CRS		White PTFE/red silicone septa + Blue screw cap with hole, for 2ml 9-425 screw top vial	100pcs
LBSV012CSS	<b>12x32mm</b>	Red PTFE/white silicone septa + Black screw cap with hole, for 2ml 8-425 screw top vial	100pcs
LBSV012CRS		White PTFE/red silicone septa + Black screw cap with hole, for 2ml 8-425 screw top vial	100pcs
LBSV222CSS		Blue PTFE/white silicone septa, Pre-slit + Blue screw cap with hole, 2ml 9-425 screw top vial	100pcs
LBSV221CSS		Pre-slit, Blue PTFE/White silicone septa+ Black screw cap with hole, for 8-425 screw top vial	100pcs
LBSV122CRS		Red PTFE/white silicone septa + Aluminium cap, for 2ml crimp top vial	100pcs
LBSV112CRS		White PTFE/red silicone septa+ Aluminium cap, for 2ml crimp top vial	100pcs
LBSV212CRS		Pre-silt, White PTFE/red silicone septa+ Aluminium cap, for 2ml crimp top vial	100pcs
LBSV222CRS		Pre-silt, Red PTFE/white silicone septa + Aluminium cap, for 2ml crimp top vial	100pcs
LBSV202CSS		White PTFE/Red silicone septa + snap cap with hole, for 2ml snap top vial	100pcs
LBSV302CSS		Red PTFE/White silicone septa + snap cap with hole, for 2ml snap top vial	100pcs
LBSV212CSS		Pre-slit, White PTFE/Red silicone septa + snap cap with hole, for 2ml snap top vial	100pcs
LBSV242CSS		Pre-slit, Red PTFE/white silicone septa + snap cap with hole, for 2ml snap top vial	100pcs
LBSV004C	<b>4ml, 13-425 Sample Vials</b>	4ml Clear vial, 15x45mm, screw top	100pcs
LBSV004A		4ml Amber vial, 15x45mm, screw top	100pcs
LBSV014CRS	<b>Caps &amp; Septa for 13-425 Sample Vials</b>	Red PTFE/White silicone septa + Black cap with hole, for 4ml screw top vial	100pcs
LBSV104CSS		Red PTFE/White silicone septa + Black cap without hole, for 4ml screw top vial	100pcs
LBSV210C	<b>Headspace Vials</b>	10ml Clear vial, 22.5x46mm, crimp top, flat bottom	100pcs
LBSV220C		20ml Clear vial, 22.5x75mm, crimp top, flat bottom	100pcs
LBSV320C		20ml Clear vial, 22.5x75mm, crimp top, round bottom	100pcs
LBSV010C		10ml Clear vial, 18mm screw top, round bottom	100pcs
LBSV120C		20ml Clear vial, 18mm screw top, round bottom	100pcs
LBSV220CSS	<b>Caps &amp; Septa for Headspace Vials</b>	White PTFE/White Silicone + Aluminium cap with hole, for 10ml/20ml crimp top vial	100pcs
LBSV322CRS		Red PTFE/White Silicone + Aluminium cap with hole, for 10ml/20ml crimp top vial	100pcs
LBSV120CSS		Blue PTFE/ White silicone septa + Silver screw cap with hole, 18mm screw top vial	100pcs
LBSV020C	<b>24-400 Vials</b>	20ml Clear vial, 24 400 screw top	100pcs
LBSV020A		20ml Amber vial, 24 400 screw top	100pcs
LBSV040C		40ml Clear vial, 24 400 screw top	100pcs
LBSV040A		40ml Amber vial, 24 400 screw top	100pcs
LBSV040CSS	<b>Caps &amp; Septa for 24-400 Vials</b>	Nature PTFE/Nature silicone septa + Black screw cap without hole, 24 400 screw top vial	100pcs
LBSV140CSS		Nature PTFE/Nature silicone septa + Black screw cap with hole, 24 400 screw top vial	100pcs