



VICI-Jour PRODUCT LIST 2012

> Current System AUD Prices

JR-68252T1	Column PEEK guard cartridge holder complete with 1 x guard column JR-6824:	pc		
JR-68252T12	Column PEEK guard cartridge holder complete with 1 x guard column JR-6824:	pc		
JR-F-197-062-2TINC	Frit Ti 0.197 diam. 0.062in thick 2 um no chamfer	pc		
JR-F-197-062-5T1	Frit Ti 0.197 diam. 0.062in thick 5 um	pc		
JR-26501	Coupler PEEK universal 0.13 mm ID red	pc		
JR-26502	Coupler PEEK universal 0.17 mm ID yellow	pc		
JR-26503	Coupler PEEK universal 0.25 mm ID blue	pc		
JR-26504	Coupler PEEK universal 0.50 mm ID orange	pc		
JR-4071-10	Plug Nylon column endstopper white 10-32 10/pkg	pkg	10	pc
JR-4072-10	Plug Nylon column endstopper black 10-32 10/pkg	pkg	10	pc
JR-4073-10	Plug Nylon column endstopper red 10-32 10/pkg	pkg	10	pc
JR-4074-10	Plug Nylon column endstopper green 10-32 10/pkg	pkg	10	pc
JR-4075-10	Plug Nylon column endstopper blue 10-32 10/pkg	pkg	10	pc
JR-4081-10	Plug Nylon column endstopper white 10-32 10 /pkg	pkg	10	pc
JR-4082-10	Plug Nylon column endstopper black 10-32 10 /pkg	pkg	10	pc
JR-4083-10	Plug Nylon column endstopper red 10-32 10 /pkg	pkg	10	pc
JR-4084-10	Plug Nylon column endstopper green 10-32 10 /pkg	pkg	10	pc
JR-4085-10	Plug Nylon column endstopper blue 10-32 10 /pkg	pkg	10	pc
JR-5005-5	Single Ferrule PEEK 1/16in 5/pkg	pkg	5	pc
JR-SLC-85-10	Sample Loop PEEK Cheminert 10 ul	pc		
JR-SLC-85-100	Sample Loop PEEK Cheminert 100 ul	pc		
JR-SLC-85-1000	Sample Loop PEEK Cheminert 1 ml	pc		
JR-SLC-85-20	Sample Loop PEEK Cheminert 20 ul	pc		
JR-SLC-85-200	Sample Loop PEEK Cheminert 200 ul	pc		
JR-SLC-85-250	Sample Loop PEEK Cheminert 250 ul	pc		
JR-SLC-85-300	Sample Loop PEEK Cheminert 300 ul	pc		
JR-SLC-85-5	Sample Loop PEEK Cheminert 5 ul	pc		
JR-SLC-85-50	Sample Loop PEEK Cheminert 50 ul	pc		
JR-SLC-85-500	Sample Loop PEEK Cheminert 500 ul	pc		
JR-SLR-2190	Sample Loop PEEK Rheodyne 1 uL	pc		
JR-SLR-2191	Sample Loop PEEK Rheodyne 2 uL	pc		
JR-SLR-2192	Sample Loop PEEK Rheodyne 3 uL	pc		
JR-SLR-2193	Sample Loop PEEK Rheodyne 4 uL	pc		
JR-SLR-2194	Sample Loop PEEK Rheodyne 5 uL	pc		
JR-SLR-2195	Sample Loop PEEK Rheodyne 20 uL	pc		
JR-SLR-2196	Sample Loop PEEK Rheodyne 10 uL	pc		
JR-SLR-2200	Sample Loop PEEK Rheodyne 30 uL	pc		
JR-SLR-2201	Sample Loop PEEK Rheodyne 40 uL	pc		
JR-SLR-2202	Sample Loop PEEK Rheodyne 50 uL	pc		
JR-SLR-2204	Sample Loop PEEK Rheodyne 100 uL	pc		
JR-SLR-2206	Sample Loop PEEK Rheodyne 200 uL	pc		
JR-SLR-2207	Sample Loop PEEK Rheodyne 250 uL	pc		
JR-SLR-2210	Sample Loop PEEK Rheodyne 500 uL	pc		
JR-SLR-2222	Sample Loop PEEK Rheodyne 1 ml	pc		
JR-SLV-75-10	Sample Loop PEEK Valco 10 ul	pc		
JR-SLV-75-100	Sample Loop PEEK Valco 100 ul	pc		
JR-SLV-75-1000	Sample Loop PEEK Valco 1 ml	pc		
JR-SLV-75-20	Sample Loop PEEK Valco 20 ul	pc		
JR-SLV-75-250	Sample Loop PEEK Valco 250 ul	pc		
JR-SLV-75-5	Sample Loop PEEK Valco 5 ul	pc		
JR-SLV-75-50	Sample Loop PEEK Valco 50 ul	pc		
JR-SLV-75-500	Sample Loop PEEK Valco 500 ul	pc		
JR-SLC-86-10	Sample Loop SS Cheminert 10 ul	pc		
JR-SLC-86-100	Sample Loop SS Cheminert 100 ul	pc		
JR-SLC-86-1000	Sample Loop SS Cheminert 1 ml	pc		
JR-SLC-86-20	Sample Loop SS Cheminert 20 ul	pc		
JR-SLC-86-250	Sample Loop SS Cheminert 250 ul	pc		