



VICI-Jour PRODUCT LIST 2012

> Current System AUD Prices

JR-SLC-86-5	Sample Loop SS Cheminert 5 ul	pc
JR-SLC-86-50	Sample Loop SS Cheminert 50 ul	pc
JR-SLC-86-500	Sample Loop SS Cheminert 500 ul	pc
JR-SLV-76-10	Sample Loop SS Valco 10 ul	pc
JR-SLV-76-100	Sample Loop SS Valco 100 ul	pc
JR-SLV-76-1000	Sample Loop SS Valco 1 ml	pc
JR-SLV-76-20	Sample Loop SS Valco 20 ul	pc
JR-SLV-76-250	Sample Loop SS Valco 250 ul	pc
JR-SLV-76-5	Sample Loop SS Valco 5 ul	pc
JR-SLV-76-50	Sample Loop SS Valco 50 ul	pc
JR-SLV-76-500	Sample Loop SS Valco 500 ul	pc
JR-SLR-111	Sample Loop SS Rheodyne 5 uL	pc
JR-SLR-111-77	Sample Loop SS Rheodyne 7725 5 uL	pc
JR-SLR-112	Sample Loop SS Rheodyne 10 uL	pc
JR-SLR-112-77	Sample Loop SS Rheodyne 7725 10 uL	pc
JR-SLR-113	Sample Loop SS Rheodyne 20 uL	pc
JR-SLR-113-77	Sample Loop SS Rheodyne 7725 20 uL	pc
JR-SLR-114	Sample Loop SS Rheodyne 50 uL	pc
JR-SLR-114-77	Sample Loop SS Rheodyne 7725 50 uL	pc
JR-SLR-115	Sample Loop SS Rheodyne 100 uL	pc
JR-SLR-115-77	Sample Loop SS Rheodyne 7725 100 uL	pc
JR-SLR-116	Sample Loop SS Rheodyne 200 uL	pc
JR-SLR-116-77	Sample Loop SS Rheodyne 7725 200 uL	pc
JR-SLR-117	Sample Loop SS Rheodyne 500 uL	pc
JR-SLR-117-77	Sample Loop SS Rheodyne 7725 500 uL	pc
JR-SLR-118	Sample Loop SS Rheodyne 1 ml	pc
JR-SLR-118-77	Sample Loop SS Rheodyne 7725 1ml	pc
JR-201531	Easy-Flange clamp pair.Spare part for JR-201540B	pc
JR-201532	Easy-Flange flanger and distance gauge	pc
JR-201536	Easy-Flange flanging disc with 1.3 mm polymer pin	pc
JR-201537	Easy-Flange flanging disc with 1.3 mm Titanium pin	pc
JR-201539	Easy-Flange combi kit in plastic case	pc
JR-201540	Easy-Flange kit in plastic case	pc
JR-201540L	Box for JR-201540 with foam insert	pc
JR-201541	Easy-Flange flanging disc with 0.8 mm polymer pin	pc
JR-201554	Easy-Flange flanging disc with 0.8 mm Titanium pin	pc
JR-202235	Easy-Flange flanging disc with 0.5 mm SS pin for PEEK tubing	pc
JR-357090	Tubing Elbow 90°	pc
JR-357180	Tubing Elbow 180°	pc
JR-794	Tool Guillotine Cutter for polymer tubing	pc
JR-795	Tool replacement blade for JR-794	pc
JR-797	Tool Clean-Cut Tubing Cutter	pc
JR-798	Tool replacement blade for JR-797	pc
JR-9001	Tubing Clip	pc
JR-792	Tool Tubing Cutter for metal tubing	pc
JR-793	Tool replacement cutter wheel for JR-792	pc
JR-796	Tool Tubing Plier for metal tubing	pc
JR-802	Tool file	pc
JR-800	Tool ValvTool 1/4in and 5/16in open wrench	pc
JR-BPR1	Backpressure Regulator variable 15 - 300psimax. 50C complete for liquids only	pc
JR-BPR1-KF	Backpressure Regulator variable 15 - 300 psimax. 50C complete for liquids on	pc
JR-BPR2	Backpressure Regulator variable 300 - 1500psimax. 50C complete for liquids o	pc
JR-BPR2-KF	Backpressure Regulator variable 300 - 1500psimax. 50C complete for liquids c	pc
JR-804	Tool wrench 1/4in x 5/16in	pc
JR-805	Tool wrench 3/8in x 7/16in	pc
JR-806	Tool Tweezers Polypropylen Pinzette	pc
JR-FRK1	Ferrule removal kit	pc