

Typical Vial & Accessories PRODUCT LIST

exGST prices as per
June 2014

(Ctrl F)
> Cat.#s




AUD

Serum Pharma Vials CLEAR; 13mm or 20mm CrimpTop

> Current System AUD Prices

CTSC-2V-13T-15D32H	2ml Clear Neutral/Schott Borosilicate Type I Glass Vial 13mm Crimp Top
CTSC-4V-13T-16D38H	4ml Clear Neutral/Schott Borosilicate Type I Glass Vial 13mm Crimp Top
CTSC-7V-20T-22D38H	7ml Clear Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top
CTSC-10V-13T-22D45H	10ml Clear Neutral/Schott Borosilicate Type I Glass Vial 13mm Crimp Top
CTSC-10V-20T-22D45H	10ml Clear Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top *
CTSC-13V-20T-24D49H	13ml Clear Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top
CTSC-13V-13T-24D49H	13ml Clear Neutral/Schott Borosilicate Type I Glass Vial 13mm Crimp Top
CTSC-20V-20T-26D66H	20ml Clear Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top
CTSC-30V-20T-29D72H	30ml Clear Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top
CTSMCW-50V-20T-43D73H	50ml Clear Moulded Type II Glass Vial Sulphur Treated and Washed 20mm Crimp Top
CTSMCW-100V-20T-52D95H	50ml Clear Moulded Type II Glass Vial Sulphur Treated and Washed 20mm Crimp Top



371	0.39
425	0.39
150	0.55
198	0.66
324	0.56
155	0.66
155	0.66
145	0.66
113	0.77
113	0.85
81	0.95

Serum Pharma Vials AMBER; 13mm or 20mm CrimpTop

CTSA-2V-13T-15D32H	2ml Amber Neutral/Schott Borosilicate Type I Glass Vial 13mm Crimp Top
CTSA-3V-13T-15D35H	3ml Amber Neutral/Schott Borosilicate Type I Glass Vial 13mm Crimp Top
CTSA-7V-20T-22D38H	7ml Amber Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top
CTSA-10V-20T-22D48H	10ml Amber Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top
CTSA-13V-20T-24D50HTP	13ml Amber Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top
CTSA-23V-20T-29D60HTP	23ml Amber Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top *
CTSAW-30V-20T-36D63H	30ml Amber Moulded Type II Glass Vial Sulphur Treated and Washed 20mm Crimp Top
CTSAW-50V-20T-43D73H	50ml Amber Moulded Type II Glass Vial Sulphur Treated and Washed 20mm Crimp Top
CTSAW-100V-20T-52D95H	100ml Amber Moulded Type II Glass Vial Sulphur Treated and Washed 20mm Crimp Top



371	0.39
425	0.39
210	0.55
155	0.77
155	0.77
113	0.77
92	0.89
84	0.95
66	0.95


Serum Stoppers; 13mm or 20mm Crimp Top ... also 32mm

CT-20BNRB-C	20mm Chlorobutyl Rubber Stopper
CT-20BNRB-FD2	20mm Grey Butyl Lyophilization Freeze Drying Rubber Stopper(2-Leg)
CT-13BNRB-C	13mm Chlorobutyl Rubber Stopper
CT-13BNRB-FD2	13mm Grey Butyl Lyophilization Freeze Drying Rubber Stopper(2-Leg)
CT-20BNRB	20mm Grey butyl rubber stopper (softer Light Gray type)
CT-20BNRB-R	20mm Butyl Rubber Stoppers - Red



0.20
0.24
0.20
0.24
0.20
0.20

CrimpTop Caps Plain Aluminium -X



0.16

CT-20CTO-Crimp-X	Aluminum Vial Seals 20mm Centre Tear-out : Color-X
CT-13CTO-Crimp-X	Aluminum Vial Seals 13mm Centre Tear-out : Color-X


0.24
0.24

Plastic FlipOff / WPFO -X "FLIP OFF" writing



0.20
0.20

Plastic FlipOff / KFO -X (plain Finish)



0.20
0.20



Hand Crimpers; Air, BenchPress, Hand Pneumatic Crimper with detachable Head ... also available

CT-Crimp-20S & -13S	HandCrimper 20mm (&13mm) Standard Plain AlumCaps
CT-Crimp-20FO & -13FO	HandCrimper 20mm (&13mm) FlipOff Caps
CT-Decap-20 & -13	Hand Decapper 20mm (&13mm) both Plain and FO Cpas



259.00
259.00
259.00

Chromatography Vials 2ml 12x32ml : Autosampler; 8mm Universal : 9mm Agilent

ScrewCap LBSV012CSS LBSV012C LBSV032C LBSV012A LBSV032A LBSV002C LBSV022C LBSV002A LBSV022A LBSV102CSS CrimpTop LBSV042C LBSV062C LBSV042A LBSV062A LBSV112CRS	<p>Red PTFE/white silicone septa + Black screw cap with hole, /2ml 8-425 screw top vial,100-pk</p> <p>2ml Clear vial, 8-425 screw top,100-pk</p> <p>2ml Clear vial, 8-425 screw top, graduated with writing area,100-pk *</p> <p>2ml Amber vial, 8-425 screw top,100-pk</p> <p>2ml Amber vial, 8-425 screw top, graduated with writing area,100-pk</p> <p>2ml Clear vial, 9-425 screw top,100-pk</p> <p>2ml Clear vial, 9-425 screw top, graduated with writing area,100-pk *</p> <p>2ml Amber vial, 9-425 screw top,100-pk</p> <p>2ml Amber vial, 9-425 screw top, graduated with writing area,100-pk</p> <p>Red PTFE/white silicone septa + Blue screw cap with hole, /2ml 9-425 screw top vial,100-pk</p> <p>2ml Clear vial, crimp top,100-pk</p> <p>2ml Clear vial, crimp top, graduated with writing area,100-pk *</p> <p>2ml Amber vial, crimp top,100-pk</p> <p>2ml Amber vial, crimp top, graduated with writing area,100-pk</p> <p>White PTFE/red silicone septa+ Aluminium cap, for 2ml crimp top vial,100-pk</p>	<p>* * *</p>  <p>* * *</p> 
--	---	---

ASK
for
current
prices

Ask f