

# Typical Vial & Accessories PRODUCT LIST

exGST prices as per June 2014

[Ctrl F]  
> Cat.#s



94

## Limited Volume Inserts for Chromatography Autosampler Vials

> Current System AUD Prices

AUD

Style : Shell or Conical with Bottom Spring Size : 100ul to 300ul  
Microvolume glass inserts for chromatography autosampler Vials ;  
Microvolume conical polypropylene insert with bottom spring for autosampler Vials

V0401100	250ul,5.8x31mm,conical glass insert
V0402100	200ul,5.8x28.5mm,conical glass insert w/bottom polypropylene spring
V0403100	200ul,5.0x31mm,conical glass insert
V0404100	150ul,5.0x28.5mm,conical glass insert,w/bottom polypropylene spring
V0405100	100ul,4.6x28.5mm,conical glass insert ,w/bottom polypropylene spring
V0411100	300ul,5.8x31mm,shell style glass insert
V0413100	250ul,5.0x31mm,shell style glass insert
V0507100	200ul,5.0x31mm,conical glass insert w/bottom spring for 1ml shell vials
V0421100	250ul,5.8x28.5mm,conical polypropylene insert w/bottom polypropylene spring



ASK

## 4ml 13-425 Chromatography Autosampler Vials; 15x45mm

4ml, 15x45mm, 13-425 screw top vials for Autosampler, clear glass or amber glass, with marking spot or without marking spot, PP Screw top caps, PTFE/silicone septa, for chromatography autosampler, 100pcs/pk.

V0501100	Vials only,4ml ,15x45mm,13-425 Clear Screw top ; 100-pk	26.00
V0503100	Vials only,4ml ,15x45mm,13-425 Amber Screw top; 100-pk	26.00
V0521100	Septa,Red PTFE/White Silicone, 12mm for 13-425 Caps ,0.06"thick; 100-pk	15.00
V0511100	Solid top caps,Black,with silicone Septa for 13-425 Screw Vials; 100-pk	20.00
V0512100	Solid top caps,white,with PP Septa for 13-425 Screw Vials; 100-pk	15.00
V0513100	Open top caps,Black,for 13-425 Screw Vials; 100-pk	13.00
V0514100	Solid top caps,with pre-assembled Red PTFE/White Silicone Septa; 100-pk	10.50
V0515100	Open top caps,Black, with pre-assembled Red PTFE/White Silicone Septa; 100-pk	24.00
V0516100	Solid top caps,with pre-assembled white PTFE Septa; 100-pk	ASK



## Screw Top Headspace Vials

Volume : 10ml, 20ml

Size : 23\*46mm, 23\*75mm

Screw Top Headspace vials for chromatography autosampler. 20ml, 23\*75mm, 10ml, 23\*46mm, round bottom, clear glass, 18mm magnetic screw top caps with PTFE/silicone septa, 100pcs/pk

V0901100	Vials only,20ml,23x75mm,Screw top , round bottom, clear glass;100-pk
V0902100	Vials only,10ml,23x46mm,Screw top, round bottom, clear glass;100-pk
V0911100	Open top Screw Magnetic cap;100-pk
V0921100	Septa, Blue PTFE/white silicone, 17.5mm ,0.06"thickness;100-pk
V0912100	Open top Screw magnetic caps with pre-assembled Septa;100-pk



ASK

## Crimp Top Headspace Vials

20ml,23\*75mm, 10ml, 23\*46mm, flat bottom, clear glass, 20mm aluminum seals with PTFE/silicone septa or 20mm Bimetal magnetic seals with PTFE/silicone septa, 100pcs/pk

V0801100	Vials only,Crimp top ,20ml,23x75mm,flat bottom,clear glass;100-pk
V0802100	Vials only,Crimp top ,10ml,23x46mm,flat bottom,clear glass;100-pk
V0811100	20mm open top,Silver Aluminum Seals;100-pk
V0813100	Magnetic Bimetal cap ,silver color ,20mm,8mm hole;100-pk
V0821100	Septa, Transparent PTFE/Transparent silicone, 20mm ,0.125"thickness;100-pk
V0822100	Septa, Tan PTFE/white silicone, 20mm ,0.125"thickness;100-pk
V0814100	20mm Silver aluminum Seal w/Transparent PTFE/Transparent silicone septa;100-pk
V0815100	20mm Silver Aluminum Seal w/Tan PTFE/white silicone Septa;100-pk
V0816100	20mm magnetic Bimetal seal w/Transparent PTFE/Transparent silicone Septa;100-pk
V0817100	20mm magnetic Bimetal seal w/Tan PTFE/white silicone Septa;100-pk
V1101001	20mm hand crimper; CT-CRIMP_20S



ASK

## Pre-Cleaned Volatile Organic Analyte (VOA) Sampling Vials

• Container, liner, and closure cleaned and assembled. • Clear or amber.  
• Open top caps. • Teflon® faced 0.125" silicone septa. • Each case lot numbered.

21798	Precleaned VOA Vials 20mL Clear Glass 24mm/400 72-pk.
21799	Precleaned VOA Vials 20mL Amber Glass 24mm/400 72-pk.
21796	Precleaned VOA Vials 40mL Clear Glass 24mm/400 72-pk.
21797	Precleaned VOA Vials 40mL Amber Glass 24mm/400 72-pk.
26121	Collection Vials for ASE 200 60mL Clear Glass 24mm/400 72-pk.
26122	Collection Vials for ASE 200 60mL Amber Glass 24mm/400 72-pk.



ASK

## ClearBorosilicateType1glassvials, 24mmScrewFinish

CT-24RS-SCR-WH	24mm Ring Seal Screw Cap White, ea
CT-24RS-SCR-BLK	24mm Ring Seal Screw Cap White, ea
CT502105B	20ml 24mm Screw Cap Vial Clear, ea
CT40MA	40ml 24mm Screw Cap Vial Clear, ea



ASK

Ask for

**CHROMalytic** +61(0)3 9762 2034  
**ECH**nology Pty Ltd  
Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

**Australian Distributors**  
Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)