

SRIGC PRODUCTS LIST 2014

(Ctrl F)
> Cat.#s



> Current System

	8670-1247	PID HV POWER SUPPLY CIRCUIT BOARD, 8610C, 310 MODEL
	8670-1248	PID DETECTOR INLET BULKHEAD FITTING
	0110-0040	PID MOUNTED ON MODEL 110 CHASSIS
	0110-4010	PID-FID MOUNTED ON MODEL 110 CHASSIS
	0110-2426	PID-FID-DELCD MOUNTED ON MODEL 110 CHASSIS
	8690-0040	PID DETECTOR

FPD	0110-0080	FPD MOUNTED ON MODEL 110 CHASSIS
	0110-0085	DFPD MOUNTED ON MODEL 110 CHASSIS
	0110-1080	FPD-FID MOUNTED ON MODEL 110 CHASSIS
	8670-0080	PHOTOMULTIPLIER (PMT) TUBE
	8670-0081	FPD DETECTOR BODY
	8670-0082	SULFUR WAVELENGTH FILTER
	8670-0083	PHOSPHORUS WAVELENGTH FILTER
	8670-0084	FPD FILTER ASSEMBLY WITH BUILT-IN SULFUR FILTER
	8670-0085	FPD FILTER ASSEMBLY WITH BUILT-IN PHOSPHORUS FILTER
	8670-0086	FPD PMT SOCKET ASSEMBLY
	8670-2085	DFPD-FID DETECTOR BODY
	8690-0080	FPD DETECTOR
	8690-0085	SIMULTANEOUS S/P FPD DETECTOR (DFPD)
	8690-1080	FPD-FID DUAL DETECTOR (DUAL ELECTRONICS)
	8690-2085	DFPD-FID DETECTOR (3 AMPLIFIERS)

TCD	8670-0007	TCD CELL w/ FILAMENTS, HEATING ELEMENT, SRI MODEL
	8670-0011	TCD GAIN SWITCH
	8670-9007	UPGRADE TCD TO CURRENT HIGH TEMPERATURE MODEL (HARDWARE ONLY)
	8670-9100	TCD CELL STAINLESS STEEL, GOWMAC-TYPE
	8670-9110	O-RING FOR TCD FILAMENT, EACH (MIN. QTY 10)
	8670-9120	TCD FILAMENT, TUNGSTEN-RHENIUM (SPECIFY COPPER OR VITON SEAL)
	8670-9121	TCD FILAMENT, GOLD PLATED, SET OF FOUR (AuW 13-300)
	8670-9122	O-RING, FOR TCD FILAMENT, PACK OF 10
	8670-9123	NUT, TCD, FILAMENT RETAINING (WAS FILAMENT W/ CU SEAL)
	8670-9124	COPPER SEAL, FOR TCD FILAMENT, PACK OF 10
	8670-9125	TCD FILAMENT PROTECTION CIRCUIT (RETROFIT)
	0110-0007	TCD MOUNTED ON MODEL 110 CHASSIS
	8690-0007	TCD DETECTOR
	8690-0008	TCD DETECTOR, 1 FILAMENT

HID	8670-0030	HID COLLECTOR ELECTRODE ASSEMBLY
	8670-0031	HID DETECTOR BODY
	8690-0030	HID DETECTOR, SRI DESIGN
	0110-0030	HID MOUNTED ON MODEL 110 CHASSIS

TID	8670-0017	TID DETECTOR BEAD
	8670-1017	TID DETECTOR BODY
	8690-0017	TID DETECTOR
	0110-0017	TID MOUNTED ON MODEL 110 CHASSIS

ASD	8690-0006	AROMATIC SELECTIVE DETECTOR
	0110-0006	ASD DETECTOR MOUNTED ON MODEL 110 CHASSIS

INJECTORS		
Autosamplers	8640-0001	LIQUID AUTOSAMPLER, 100 VIAL TRAY, COBRA (VERTICAL, FOR 8610V ONLY)
	8640-0002	MOUNTING BRACKET KIT, FOR VERTICAL AUTOSAMPLER (8610V ONLY)
	8640-0003	SYRINGE, FOR COBRA AUTOSAMPLER, 5 CM. NEEDLE
	8640-0050	AUTOCALIBRATION SYSTEM, GAS
	8640-1006	SINGLE BAY AUTOINJECTOR
	8670-0068	SYRINGE NEEDLE, REPLACEMENT, FOR 42 VIAL LIQUID AUTOSAMPLER
	8670-2368	INJECTION PORT COUPLER TO LIQUID AUTOSAMPLER CAROUSEL
	8690-0068	LIQUID AUTOSAMPLER CAROUSEL, 42 VIAL TRAY
	8690-0053	10 STATION PURGE-AND-TRAP SAMPLE SELECTOR
	8690-0075	ON-LINE LIQUID SAMPLER, FOR P&T
	8690-0087	DGA-TOGA PERMEATION LOOP ACCESSORY
	8690-1068	ON-LINE SYRINGE SAMPLER (NO CAROUSEL)
	8690-2253	10 STATION PURGE-AND-TRAP SAMPLE SELECTOR 220 V
Direct Injection	8670-0024	HEATED INJECTOR BODY