

# SRIGC PRODUCTS LIST 2014

(Ctrl F)  
> Cat.#s



> Current System AUD Prices

8680-0052	METHOD 5035 COMPLIANT VOA VIAL SPARGE ACCESSORY
8680-0352	SPARE DESSICANT CANISTER
8690-0006	AROMATIC SELECTIVE DETECTOR (ASD)
8690-0007	TCD DETECTOR
8690-0008	TCD DETECTOR, 1 FILAMENT
8690-0009	RGD DETECTOR (WITH EMPTY REACTOR)
8690-0010	FID DETECTOR
8690-0011	BLANK REACTION TUBE W/FRITS
8690-0015	NPD DETECTOR
8690-0017	TID DETECTOR
8690-0020	ECD DETECTOR
8690-0023	ADDITIONAL ON-COLUMN INJECTOR WITH EPC CIRCUIT
8690-0025	HEATED INJECTION PORT UPGRADE
8690-0030	HID DETECTOR, SRI DESIGN
8690-0034	SPLIT / SPLITLESS INJECTOR, AUTOMATED VENT
8690-0035	GAS PHASE SPLITTER
8690-0037	ON-COLUMN PTV INJECTOR UPGRADE
8690-0040	PID DETECTOR
8690-0045	HEATED STATIC HEADSPACE INJECTOR
8690-0046	HEATED STATIC HEADSPACE INJECTOR UPGRADE, TO EXISTING VALVE, VALVE OVEN
8690-0051	AUTOMATED PURGE AND TRAP CONCENTRATOR, METHOD 5030
8690-0052	AUTOMATED PURGE AND TRAP CONCENTRATOR, METHOD 5030/35
8690-0053	10 SAMPLE PURGE AND TRAP AUTOSAMPLER, w/4 RELAY BD
8690-0054	10 SAMPLE METHOD 5035 PURGE-AND-TRAP AUTOSAMPLER
8690-0060	10 PORT GAS SAMPLING VALVE, MULTIGAS 3 PLUMBING, WITH SOLENOID
8690-0061	10 PORT GAS SAMPLING VALVE, PNEUMATIC ACTUATOR
8690-0062	4 PORT LIQUID SAMPLING VALVE, ELECTRIC ACTUATOR
8690-0063	10 PORT GAS SAMPLING VALVE, MANUAL ACTUATOR
8690-0064	4 PORT LIQUID SAMPLING VALVE, MANUAL ACTUATOR
8690-0065	AUTOMATED GAS SAMPLE VALVE, 10 PORT
8690-0066	10 INPUT HEATED STREAM SELECTOR W/ TRANSFER LINE
8690-0067	BLOCK HEATER w/ CARTRIDGE FOR GAS VALVE
8690-0068	LIQUID AUTOSAMPLER CAROUSEL, 42 VIAL TRAY
8690-0069	SOLENOID-CONTROLLED GAS VALVE
8690-0070	INTERNAL WHISPER-QUIET AIR COMPRESSOR, 115V
8690-0071	PERISTALTIC PUMP, REPLACEMENT
8690-0072	SPLIT VENT VALVE
8690-0073	VACUUM PUMP INTERFACE, DATA SYSTEM CONTROLLED
8690-0074	THERMOSTATTED VALVE OVEN, HIGH TEMP. (300°C MAX.)
8690-0075	ON-LINE LIQUID SAMPLER FOR PURGE-AND-TRAP
8690-0077	22 PORT, 10 STREAM SELECTOR VALVE, AUTOMATED, STD. ALONE
8690-0078	22 PORT, 10 STREAM SELECTOR VALVE, ELECTRIC ACTUATOR, PLUMBED
8690-0079	22 PORT, 10 STREAM SELECTOR VALVE, ELECTRIC ACTUATOR, UNPLUMBED
8690-0080	FPD DETECTOR
8690-0081	METHANIZER ACCESSORY BUILT ONTO EXISTING VALVE OVEN
8690-0082	METHANIZER ACCESSORY, IN JET
8690-0083	TO-14 CRYOTRAP LOOP, AUTOMATED CONTROL
8690-0084	HEATED / FAST COOLING ADSORBENT TRAP FOR EXISTING GAS VALVE
8690-0085	SIMULTANEOUS S/P FPD DETECTOR (DFPD)
8690-0086	PELTIER COOLING w/ ELECTRONICS, FOR EXISTING TRAP
8690-0087	DGA-TOGA PERMEATION LOOP ACCESSORY
8690-0088	THERMOSTATTED VALVE OVEN
8690-0089	UNHEATED INSULATED VALVE OVEN
8690-0090	PELTIER COOLED TRAP, FOR GAS SAMPLING VALVE
8690-0091	THERMOSTATTED VALVE OVEN, NARROW, FOR 8610D GC
8690-0092	ENRICHMENT COIL ACCESSORY
8690-0093	ONLINE MEMBRANE SAMPLER W/ SPARGER FOR COOLING TOWER WATER
8690-0094	OXIDIZER ACCESSORY
8690-0095	REUSEABLE SHIPPING CONTAINER
8690-0096	DGA-TOGA EXTRACTOR ASSEMBLY
8690-0097	TRAP PACKING ACCESSORY
8690-0151	NAFION MEMBRANE REPLACEMENT FOR HYDROGEN GEN.
8690-0152	NAFION TUBING, DRYER, FOR HYDROGEN GENERATOR
8690-0153	WATER FILTER, FOR HYDROGEN GENERATOR
8690-0160	TURBO-COOL FAST OVEN COOLING
8690-0230	230VAC 50Hz CONVERSION KIT
8690-1015	NPD CONVERSION KIT FOR FID