

21117	Split Liner	2mm x 2.75 x 120 for Thermo GCs IP
21117-200.5	Split Liner	2mm x 2.75 x 120 for Thermo GCs IP
21117-202.5	Split Liner	2mm x 2.75 x 120 for Thermo GCs IP
21117-209.5	Split Liner	2mm x 2.75 x 120 for Thermo GCs IP Deact., w/CarboFrit, 5pk
21117-210.5	Split Liner	2mm x 2.75 x 120 for Thermo GCs Base Deact., 5pk
21117-211.5	Split Liner	2mm x 2.75 x 120 for Thermo GCs Base Deact., w/Base Deact. Wool, 5pk
21117-213.5	Split Liner	2mm x 2.75 x 120 for Thermo GCs Siltek Deact., w/Deact. Wool, 5pk
21117-214.5	Split Liner	2mm x 2.75 x 120 for Thermo GCs Siltek Deact., 5pk
21117-216.5	Split Liner	2mm x 2.75 x 120 for Thermo GCs Siltek Deact., w/CarboFrit, 5pk
21118	Syringe Adapter Kit	For use with VOC single stage regulators
21121	Shimadzu Ferrules	5mm, Pack of 10
21122	Ferrule, Teflon	1/16", Pack of 10
21123	Ferrule, Teflon	1/16" x 0.4mm Pack of 10
21124	Ferrule, Teflon	1/16" x 0.5mm Pack of 10
21125	Ferrule, Teflon	1/16" x 0.8mm Pack of 10
21126	Ferrule, Teflon	1/8", Pack of 10
21127	Ferrule, Teflon	3/16", Pack of 10
21128	Ferrule, Teflon	1/4", Pack of 10
21131	FID Collector Assembly	Siltek, for Agilent 5890 GCs
21132	FID Collector Assembly	Siltek, for Agilent 6890/6850/7890 GCs
21135	FID Castle	Siltek, for Agilent 5890/6890/6850/7890 GCs
21136	FID Collector Nut	and Washer for Agilent 5890/6890/6850/7890 GCs
21137	Ignitor Castle	for Agilent 5890/6890/6850/7890 GCs
21138	FID Collector	for Agilent 5890 GCs
21139	FID Collector Assembly	Assembly Kit for Agilent 6890/6850/7890 GCs
21140	Vials	2ml Clear 9mm Thread w/Grad Marking Spot 100PK
21141	Vials	2ml Clear 9mm Thread w/Grad Marking Spot 1000PK
21142	Vials	2ml Amber 9mm Thread w/Grad Marking Spot 100PK
21142-220	Vials	2ml Amber Thread w/Grad Marking Spot, IP Deact. 100pk
21143	Vials	2ml Amber 9mm Thread w/Grad Marking Spot 1000PK
21147	Vial Cap, Septa	8mm x 0.065" Pack of 1000 Red PTFE/White Silicone
21148	Vial Cap, Septa	10mm x 0.060" 1000Pk Red PTFE/White Silicone
21149	Vial Closure	8-425 Blk Closure 0.065" Red PTFE/White Silic 100PK
21150	Fritted sparger,5mL,EA	
21151	Fritted sparger,25mL,EA	
21152	Vial, Crimp Top	2ml Clr 12x32mm 11mm Finish Pack of 100 (no marking spot)
21153	Vial, Crimp Top	2ml Clr 12x32 11mm Finish Pack of 1000 (no marking spot)
21154	Vial, Short Cap	2ml Clr 12x32 9mm Screw Thread Finish Pack of 100
21154-220	Vial, Short Cap	2mL Clr 12x32 9mm Screw Thread Finish Pack of 100, IP Deact
21155	Vial, Short Cap	2ml Clr 12x32mm 9mm Screw Thread Finish Pack of 1000
21156	CIS4/PTV Liner	2mm x 3.0 x 71 for Agilent GCs IP Deact., w/Deact. Wool,10pk
21156-211.10	CIS4/PTV Liner	2mm x 3.0 x 71 for Agilent GCs Base Deact., w/Base Deact. Wool, 10pk
21156-231.10	CIS4/PTV Liner	2mm x 3.0 x 71 for Agilent GCs IP Deact., w/ SV Wool, 10pk
21157	Straight Glass PTV Liner	2mm x 3.0 x 71 for Agilent GCs IP Deact., 10pk
21157-209.10	Straight Glass PTV Liner	2mm x 3.0 x 71 for Agilent GCs IP Deact., w/CarboFrit, 10pk
21157-210.10	Straight Glass PTV Liner	2mm x 3.0 x 71 for Agilent GCs Base Deact., 10pk
21157-211.10	Straight Glass PTV Liner	2mm x 3.0 x 71 for Agilent GCs Base Deact., Base Wool, 10pk
21157-213.10	Straight Glass PTV Liner	2mm x 3.0 x 71 for Agilent GCs Siltek Deact., w/Deact. Wool, 10pk
21157-214.10	Straight Glass PTV Liner	2mm x 3.0 x 71 for Agilent GCs Siltek Deact., 10pk
21157-216.10	Straight Glass PTV Liner	2mm x 3.0 x 71 for Agilent GCs Siltek Deact., w/CarboFrit, 10pk
21159	Shimadzu Injector Wrench	For Shimadzu 17A, 2010 & 2014 GCs
21160	Vials, Headspace	27ml Clear 30x60mm Pack of 100
21161	Vials, Headspace	27ml Clear 30x60mm Pack of 1000
21162	Vials, Headspace	20ml Clear 23x75mm Rounded Bottom Pack of 100
21163	Vials, Headspace	20ml Clear 23x75mm Rounded Bottom Pack of 1000
21163-221	Vials, Headspace	20mL Clear 23x75mm Rounded Bottom Pack of 1000 IP Deactivated
21164	Vials, Headspace	10ml Clear 23x46mm Rounded Bottom Pack of 100
21165	Vials, Headspace	10ml Clear 23x46mm Rounded Bottom Pack of 1000
21165-221	Vials, Headspace	10ml Clear 23x46mm Deact. Rounded Bottom Pack of 1000 Deactivated
21166	Vials, Headspace	6ml Clear Pack of 100
21167	Vials, Headspace	6ml Clear Pack of 1000
21172	Injector Mounting Post	for Agilent 7683 Series Autosamplers
21174	Aluminum Crimp Seals	11mm, PTFE/Rubber Pack of 100
21175	Aluminum Crimp Seals	11mm, PTFE/Rubber Pack of 1000
21176	Vial Cap	8-425 Hole, Polypropylene Pack of 100
21177	Vial Cap	8-425 Hole, Polypropylene Pack of 1000

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034



**HROMalytic** +61(0)3 9762 2034  
**ECH**nology Pty Ltd

**Australian Distributors**  
 Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)

Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

(Ctrl F)  
> Cat.#s

> Current System AUD Prices

