

21178	Vial Cap, Septa	8mm Red PTFE/Silicone Pack of 100
21179	Vial Cap, Septa	8mm Red PTFE/Silicone Pack of 100
21192	Vial Kits, Screw Cap	2.0ml Clear, Untreated Pack of 100
21193	Vial Kits, Screw Cap	2.0ml Clear, Untreated Pack of 1000
21194	Vial Kits, Screw Cap	2.0ml Amber, Untreated Pack of 100
21195	Vial Kits, Screw Cap	2.0ml Amber, Untreated Pack of 1000
21196	Vial Kits, Crimp Cap	2.0ml Clear, Untreated Pack of 100
21197	Vial Kits, Crimp Cap	2.0ml Clear, Untreated Pack of 1000
21198	Vial Kits, Crimp Cap	2.0ml Amber, Untreated Pack of 100
21199	Vial Kits, Crimp Cap	2.0ml Amber, Untreated Pack of 1000
21202	Syringe, SGE	Varian CP 8400, CP 8410 10ul/ASN-F/26/50mm/Bevel
21205	Syringe	Dynatech Hvy Duty P&T 2ml
21206	Syringe	Dynatech Hvy Duty P&T 5ml
21207	Syringe	H.D. P&T 10ml w/Sample Lok
21208	Syringe	H.D. P&T 5ml w/Sample Lok
21209	Syringe	Dynatech Hvy Duty P&T 10ml
21210	Syringe, SGE	5F-HP-0.63/.47 5ul/23-26s/42mmCone/F HP 7673 Autosampler
21211	Syringe, SGE	5R-HP-.063/.47 5ul/23-26s/42mmCone/RN HP7673
21212	Syringe, SGE	10F-HP-.63/.47 10ul/23-26s/42mmCone/RN
21213	Syringe, SGE	10R-HP-.63/.47 10ul/23-26s/42mmCone/RN HP7673
21214	Syringe, SGE	SK-5F-HP.63/.47 5UL/23-26s/42mmCone/F HP7673 Autosampler Pack of 6
21215	Syringe, SGE	SK-10F-HP.63/.47 10ul/23-26s/42mmCone/F HP7673 Autosampler Pack of 6
21216	Needle	N5-HP.63/.47 5ul/23-26s/42mmCone HP7673 Autosampler Pack of 2 -
21217	Needle	N10-HP-.63/.47 10ul/23-26s/42mmCone HP7673 Autosampler Pack of 2
21218	Plunger	P10F-HP-GT 10ul/GT/F HP7673 Autosampler Pack of 2
21220	Syringe, SGE	10R-HP-GT-.63/.47 10ul/23-26s/42mmCone HP7673 Autosampler
21221	Syringe, Hamilton	.5ul/26s/ASRN/HP7673
21222	Syringe, Hamilton	5ul/26s/ASRN/HP 7673
21223	Syringe, Hamilton	5ul/23-26s/ASRN/HP 7673
21224	Syringe, Hamilton	10ul/23s/ASRN/GT/HP 7673
21228	Syringe, Hamilton	10ul/26s/GT/ASRN/HP7673
21229	Syringe, Hamilton	10ul/23-26s/GT/ASRN/HP 7673
21230	Syringe, Hamilton	5ul/26s/GT/F/HP 7673 Pack of 6
21237	Injector Mounting Post	for Agilent 7673 Series Autosamplers, for use with 6890 GCs
21238	SS Inlet Seal, 0.8mm ID	Dual Vespel Ring, 2pk
21239	SS Inlet Seals	Dual Vespel Ring, 0.8mmID 10-pk
21240	Gold Plated Inlet Seal	Dual Vespel Ring, 0.8mmID 2pk
21241	Gold Plated Inlet Seals	Dual Vespel Ring, 0.8mmID 10pk
21242	Siltek Inlet Seal	Dual Vespel Ring, 0.8mmID 2pk
21243	Siltek-Treated Inlet Seal	Dual Vespel Ring, 0.8mmID 10pk
21244	SS Inlet Seal 1.2mm ID	Dual Vespel Ring 2pk
21245	SS Inlet Seal 1.2mm ID	Dual Vespel Ring 10pk
21246	Gold Plated Inlet Seal	Dual Vespel Ring 1.2mm ID 2pk
21247	Gold Plated Inlet Seal	Dual Vespel Ring 1.2mm ID 10pk
21248	Siltek Inlet Seal	Dual Vespel Ring 1.2mm ID 2pk
21249	Siltek Inlet Seal	Dual Vespel Ring 1.2mm ID 10pk
21250	Syringe, Hamilton	701 10ul LC Syringe Solid Needle for Rheodyne
21251	Syringe, Hamilton	702 25ul (22s/2"/3) LC Syringe Solid Needle for Rheodyne
21252	Syringe, Hamilton	705 50ul 22s LC Syringe Solid Needle for Rheodyne
21253	Syringe, Hamilton	710SNR 100ul LC Syringe (22s/2"/3) Solid Needle for Rheodyne
21254	Syringe, Hamilton	725SNR 250ul LC Syringe (22/2"/3) Solid Needle for Rheodyne
21255	Syringe, Hamilton	801 10ul LC Syringe Removable Needle for Waters
21256	Syringe, Hamilton	802 25ul LC Syringe Removable Needle for Waters
21257	Syringe, Hamilton	805 50ul LC Syringe Removable Needle for Waters
21258	Syringe, Hamilton	810 100ul LC Syringe Removable Needle for Waters
21259	Syringe, Hamilton	825 250ul LC Syringe Removable Needle for Waters
21260	Syringe, Hamilton	1701 10ul LC Syringe Gastight for Rheodyne
21261	Syringe, Hamilton	1702 25ul LC Syringe Gastight for Rheodyne
21262	Syringe, Hamilton	1705 50ul LC Syringe Gastight for Rheodyne
21263	Syringe, Hamilton	1710 100ul LC Syringe Gastight for Rheodyne
21264	Syringe, Hamilton	1725 250ul LC Syringe Gastight for Rheodyne
21265	Syringe, Hamilton	1010W Priming 10ml Gastight for Waters
21270	Syringe Needles, Hamilton	25S/1.97"/3 Repl. Needle 10-100ul Syringe Volume Pack of 6
21271	Syringe Needles, Hamilton	25S/1.97"/3 Repl. Needle 250ul-10mL Syringe Volume Pack of 6
21272	Splitless Liner	1mm x 6.20 x 92.1 for PE GCs IP Deact., ea
21272-214.1	Splitless Liner	1mm x 6.20 x 92.1 for PE GCs Siltek Deact., ea

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

 <b>HROMalytic</b> +61(0)3 9762 2034 <b>ECHnology</b> Pty Ltd		<b>Australian Distributors</b> Importers & Manufacturers <a href="http://www.chromtech.net.au">www.chromtech.net.au</a>
Website NEW : <a href="http://www.chromalytic.net.au">www.chromalytic.net.au</a> E-mail : <a href="mailto:info@chromtech.net.au">info@chromtech.net.au</a> Tel: 03 9762 2034 . . . in AUSTRALIA		

(Ctrl F)  
> Cat.#s

> Current System AUD Prices

