


		Page 71	(Ctrl F) > Cat.#s	
21340	Flexible Hose, SS	18" Long		
21347	Regulator Manifold	for CGA 580 (N2,He,Ar) Chrome Plated Brass Protocol Station for Single or Dual Stage		
21348	Regulator Manifold	for CGA 350 (H2,P5) Chrome Plated Brass Protocol Station for Single or Dual Stage		
21349	Regulator Manifold	for CGA 590 (Air), Chrome Plated Brass Protocol Station for Single or Dual Stage		
21353	Bore Brushes	Replacements for Injection Port Repair Tool, 2pk		
21356	Tool	HPLC Piston Seal Insertion Tool		
21360-1	Silcosteel Regulator	Dual Stage CGA 240		
21360-10	Silcosteel Regulator	Dual Stage CGA 580		
21360-11	Silcosteel Regulator	Dual Stage CGA 590		
21360-12	Silcosteel Regulator	Dual Stage CGA 660		
21360-13	Silcosteel Regulator	Dual Stage CGA 670		
21360-14	Silcosteel Regulator	Dual Stage CGA 678		
21360-15	Silcosteel Regulator	Dual Stage CGA 705		
21360-16	Silcosteel Regulator	Dual Stage CGA 300		
21360-2	Silcosteel Regulator	Dual Stage CGA 330		
21360-4	Silcosteel Regulator	Dual Stage CGA 300		
21360-5	Silcosteel Regulator	Dual Stage CGA 320		
21360-6	Silcosteel Regulator	Dual Stage CGA 326		
21360-7	Silcosteel Regulator	Dual Stage CGA 350		
21360-8	Silcosteel Regulator	Dual Stage CGA 510		
21360-9	Silcosteel Regulator	Dual Stage CGA 540		
21361-1	Silcosteel Regulator	Single Stage CGA 240		
21361-10	Silcosteel Regulator	Single Stage CGA 590		
21361-11	Silcosteel Regulator	Single Stage CGA 660		
21361-12	Silcosteel Regulator	Single Stage CGA 670		
21361-13	Silcosteel Regulator	Single Stage CGA 678		
21361-14	Silcosteel Regulator	Single Stage CGA 705		
21361-15	Silcosteel Regulator	Single Stage CGA 580		
21361-2	Silcosteel Regulator	Single Stage CGA 300		
21361-3	Silcosteel Regulator	Single Stage CGA 320		
21361-4	Silcosteel Regulator	Single Stage CGA 326		
21361-5	Silcosteel Regulator	Single Stage CGA 330		
21361-6	Silcosteel Regulator	Single Stage CGA 350		
21361-7	Silcosteel Regulator	Single Stage CGA 510		
21361-8	Silcosteel Regulator	Single Stage CGA 540		
21361-9	Silcosteel Regulator	Single Stage CGA 580		
21384	Connector Kit, Siltek	MXT Connector for 0.53mmID Fused Silica Columns		
21385	Connector Kit, Siltek	MXT Connector for 0.32mmID Fused Silica Columns		
21386	Connector Kit, Siltek	MXT Connector for 0.25mmID Fused Silica Columns		
21387	Connector Kit, Siltek	MXT Y-Connector Kits for 0.53mmID Fused Silica Columns		
21388	Connector Kit, Siltek	MXT Y-Connector Kits for 0.32mmID Fused Silica Columns		
21389	Connector Kit, Siltek	MXT Y-Connector Kits for 0.25mmID Fused Silica Columns		
21390-214.1	Drilled Uniliner Liner	Hole on Top, 1mm x 6.3 x 78.5 for Agilent GCs, Siltek Deact., ea		
21391-214.5	Drilled Uniliner Liner	Hole on Top, 1mm x 6.3 x 78.5 for Agilent GCs, Siltek Deact., 5pk		
21392	Inlet Liner Seal	for Trace PTV, 2pk		
21393	Tool	Injection Port Repair Tool for Agilent GC`s		
21394	Sample Cylinder,Sulfinert	2250cc, 1800psig, 304L, SS, 1/4" Female NPT Ends		
21394-PI	Sample Cylinder	Sulfinert, TPED, 2250cc, 1450psig, 304L, SS, 1/4" Female NPT Ends		
21399	Installation Tool	1/16" Capillary Installation Tool		
21400	Sample Cylinder Valve	Sulfinert, 1/4" Male NPT x 1/4" Male NPT		
21401	Sample Cylinder Valve	Sulfinert, 1/4" NPT x 1/4" Compression		
21402-1	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 1-1/2 inches		
21402-10	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 5-1/4 inches		
21402-11	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 4-1/2 inches		
21402-12	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 3-5/8 inches		
21402-13	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 3.0 inches		
21402-14	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 3-5/8 inches		
21402-15	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 4-5/8 inches		
21402-16	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 1-3/4 inches		
21402-17	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 3-7/8 inches		
21402-18	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 4-1/2 inches		
21402-19	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 2-3/8 inches		
21402-2	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 1-5/8 inches		
21402-21	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 4-3/8 inches		
21402-22	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 3-7/8 inches		
21402-23	Sample Cylinder Valve	Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 2-3/8 inches		

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

	+61(0)3 9762 2034	<b>Australian Distributors</b> Importers & Manufacturers www.chromtech.net.au
Website NEW : <a href="http://www.chromalytic.net.au">www.chromalytic.net.au</a> E-mail : <a href="mailto:info@chromtech.net.au">info@chromtech.net.au</a> Tel: 03 9762 2034 . . . in AUSTRALIA		