

22080	Retaining Nut	Heater Block Retaining Nut for Agilent
22081	O-rings	for High Capacity Oxygen Traps 10pk
22082	Reducing Union Kit	for Transfer Lines, 1/8"-1/16" SS w/1/4"
22083	Inlet Seals	Dual Vespel Ring, 0.8mm Gold Plated Cross Disk, 2pk
22084	Inlet Seals	Dual Vespel Ring, 0.8mm Gold Plated Cross Disk, 10pk
22085	Inlet Seals	Dual Vespel Ring, 0.8mm Siltek Treated Cross Disk, 2pk
22086	Inlet Seals	Dual Vespel Ring, 0.8mm Siltek Treated Cross Disk, 10pk
22087	Inlet Seals	Discontinue When Stock is Depleted Dual Vespel Ring, 0.8mm Stainless Steel Cr
22088	Inlet Seals	Discontinue When Stock is Depleted Dual Vespel Ring, 0.8mm Stainless Steel Cr
22089	Gold Plated Liner Cap	for Split/Splitless Injector
22090	SilcoCan Canister	1L SilcoCan, No Valve
22091	SilcoCan Canister	3L SilcoCan, No Valve (cat #563646 is needed if connecting to a sampling kit)
22093	SilcoCan Canister	15L SilcoCan, No Valve
22094	TO Canister	1L TO-Can, No Valve
22095	TO Canister	3L TO-Can, No Valve (cat #563647 is needed if connecting to a sampling kit)
22097	TO Canister	15L TO-Can, No Valve
22098	Orifice	Replacement, 80-340 sccm, Non-Siltek
22099	Orifice	Replacement, 80-340 sccm Siltek Treated
22103	Veriflo Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Replacement Orifice, 80-340 sc
22104	Tedlar Bag	Replacement Septa for Tedlar Bag Valves, 4mm Diameter, Teflon Faced Silicone,
22105	TO-Can Canister	1 Liter, 1/4" Swagelok SS4H Bellows Valve 2-port only, Stainless Steel only
22106	TO-Can Canister	3 Liter, 1/4" Swagelok SS4H Bellows Valve 2-port only, Stainless Steel only (cat #
22107	TO-Can Canister	6 Liter, 1/4" Swagelok SS4H Bellows Valve 2-port only, Stainless Steel only
22108	TO-Can Canister	15 Liter, 1/4" Swagelok SS4H Bellows Valve 2-port only, Stainless Steel only
22110	Sample Cylinder Valve	Discontinue when Stock is Depleted Sulfinitert, 6000psig, 1/4in Male NPT x 1/4in F
22111	Sample Cylinder,Sulfinitert	150cc, 5000psig, 316L, SS, 1/4" Female NPT Ends
22111-PI	Sample Cylinder	Sulfinitert, TPED, 150cc, 4350psig, 316L, SS, 1/4" Female NPT Ends
22112	Sample Cylinder,Sulfinitert	300cc, 5000psig, 316L, SS, 1/4" Female NPT Ends
22112-PI	Sample Cylinder	Sulfinitert, TPED, 300cc, 4350psig, 316L, SS, 1/4" Female NPT Ends
22113	Sample Cylinder,Sulfinitert	500cc, 5000psig, 316L, SS, 1/4" Female NPT Ends
22113-PI	Sample Cylinder	Sulfinitert, TPED, 500cc, 4350psig, 316L, SS, 1/4" Female NPT Ends
22114	PUF Cartridge	Large, 65mm OD x 125mm L, 75mm PUF
22115	PUF Cartridge	Large PUF/XAD Cartridge, 65mm OD x 125mm L, 25mm PUF/10g XAD-2/50mm I
22116	PUF Cartridge	Small, 22mm OD x 100mm L, 76mm PUF
22117	PUF Cartridge	Small PUF/XAD Cartridge, 22mm OD x 100mm L, 30mm PUF/1.5g XAD-2/30mm
22118	Sampling Pump	Vacuum Bag Sampler Model 1062, 1-10L Gas Sampling Bag Capacity, Includes F
22119	Sampling Pump	Replacement Battery for Vacuum Bag Sampler Model 1062
22120	Sampling Pump	Universal Battery Charger for Vacuum Bag Sampler Model 1062
22121	Ashcroft Gauge Connector	to Canister Valve, SS, Connects 1/4" Male NPT to 1/4" Male Compression Fitting
22122	Mini-Can Accessory	Sampling Belt, Adjustable 50in Long, Secures Mini-Can or Other Device w/2 Velcr
22123	Mini-Can Accessory	Personal Sampling Inlet, 3ft x 1/16in OD Teflon Tubing, Clip, Teflon Reducing Fer
22124	Mini-Can Accessory	Mini-Can Stand, Fits 2 3/4in Diameter Cans, Collapsible, 3-Leg, SS
22125	Humidifier	For Wasson Can Cleaner
22126	3 L Option	Includes everything for 12, 3 L Air Cans
22127	Mini-Can Option	Includes everything for either 48, 400cc or 48, 1000cc mini-cans
22128	Column Hanger	for Agilent GCs For 7890, 6890, 5890, and 5880A GCs
22129	Generator, Nitrogen	NitroFlow Lab Model for LC/MS 32 L/min Generating Capacity 12 month warranty
22130	Generator, Nitrogen	N2-04 Model for ELSD 8 L/min Generating Capacity 12 month warranty from date
22131	Generator, Nitrogen	N2-22 Model, 44L/min Generating Capacity 12 month warranty from date of purch
22132	Generator, Nitrogen	N2-22A Model, 29 L/min Generating Capacity 12 month warranty from date of pur
22133	Generator, Nitrogen	N2-35 Model, 44 L/min Generating Capacity 12 month warranty from date of purch
22134	Generator, Nitrogen	N2-35A Model, 44 L/min Generating Capacity 12 month warranty from date of pur
22135	Generator, Nitrogen	Maintenance Kit w/Carbon Filter
22171	Split Precision Liner	2.0mm x 5.0 x 95 for Shimadzu GCs IP Deact., w/Deact. Wool, ea
22172	Split Precision Liner	2.0mm x 5.0 x 95 for Shimadzu GCs IP Deact., w/Deact. Wool, 5-pk
22173	Single Gooseneck	Precision Liner, 3.5mm x 5.0 x 95 for Shimadzu GCs, IP Deact., w/Deact. Wool, e
22174	Single Gooseneck	Precision Liner, 3.5mm x 5.0 x 95 for Shimadzu GCs, IP Deact., w/Deact. Wool, 5
22175	EZ Twist Top	Replacement Shell Weldement for Agilent 7890 GCs, DO NOT SHIP TO JAPAN (
22176	EZ Twist Top	Replacement Shell Weldement for Agilent 7890 GCs, Siltek Treated, DO NOT SH
22177	EZ Twist Top	Complete Injection Port for Agilent 7890 GCs, DO NOT SHIP TO JAPAN OR GE
22178	EZ Twist Top	Complete Injection Port for Agilent 7890 GCs, Siltek Treated, DO NOT SHIP TO J
22179	FID Kit	Detector Maintenance/Start-up Kit for Agilent 6850/6890/7890 GCs
22180	FID Kit	Detector Maintenance/Start-up Kit for Agilent 5890 GCs
22181	Inlet Kit	Injector Maintenance/Start-up Kit for Agilent 5890/6850/6890/7890 GCs
22182	MLE Capillary Tool Kit	for Shimadzu GCs
22183	MLE Capillary Tool Kit	for Thermo Electron GCs

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034



HROMalytic +61(0)3 9762 2034

ECHnology Pty Ltd

Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors

Importers & Manufacturers

www.chromtech.net.au

