

22781	Merlin Microseal
22782	Merlin Microseal
22783	Luer Syringe Needles
22784	Luer Syringe Needles
22788	Certified Vials
22789	Certified Vials
22790	Certified Vials
22791	Certified Vials
22792	Certified Vials
22793	Certified Vials
22800	Split Vent Line
22808	Merlin Microseal
22809	Merlin Microseal
22810	Merlin Microseal
22811	Merlin Microseal
22812	Merlin Microseal
22813	Merlin Microseal
22814	Merlin Microseal
22815	Merlin Microseal Septum
22820	Trap, Chemical
22831	Magnetic Seal w/Septa
22832	Magnetic Seal w/Septa
22833	Magnetic Seal w/Septa
22834	Magnetic Seal w/Septa
22835	Magnetic Seal w/Septa
22836	Magnetic Seal w/Septa
22837	Magnetic Seal w/Septa
22839	Restek Electronic
22839-R	Leak Detector Conformance
22840	Sample Loop, Sulfinert
22841	Sample Loop, Sulfinert
22842	Sample Loop, Sulfinert
22843	Sample Loop, Sulfinert
22844	Sample Loop, Sulfinert
22845	Sample Loop, Sulfinert
22846	Sample Loop, Sulfinert
22847	Sample Loop, Sulfinert
22848	Sample Loop, Sulfinert
22849	Sample Loop, Sulfinert
22850	Sample Loop, Sulfinert
22855	Turret Tray Assembly
22856	Vial Storage Rack
22857	Vial Storage Rack
22858	Capillary Column Caps
22863	Silcosteel-CR, Fitting
22864	Silcosteel-CR, Fitting
22865	Silcosteel-CR, Fitting
22866	Silcosteel-CR, Fitting
22867	Silcosteel-CR, Fitting
22868	Silcosteel-CR, Fitting
22869	Silcosteel-CR, Fitting
22870	Silcosteel-CR, Fitting
22871	Silcosteel-CR, Fitting
22872	Silcosteel-CR, Fitting
22873	Silcosteel-CR, Fitting
22875	Silcosteel-CR, Fitting
22876	Silcosteel-CR, Fitting
22878	Silcosteel-CR, Fitting
22879	Silcosteel-CR, Fitting
22880	Silcosteel-CR, Fitting
22881	Silcosteel-CR, Fitting
22882	Silcosteel-CR, Fitting
22883	Silcosteel-CR, Fitting
22884	Silcosteel-CR, Fitting
22885	Silcosteel-CR, Fitting
22886	Silcosteel-CR, Fitting

Adapter Kit for Perkin Elmer Inlets (pc
Replacement Septum for SPME Appli
for 100uL A2 Gas-Tight Syringes, 22ga, 3-pk, (0.028"OD x 0.012"ID x 2")
for 500uL A2 Gas-Tight Syringes, 22ga, 3-pk, (0.028"OD x 0.012"ID x 2")
Assembled 2mL Amber Vial w/Marking Spot and Bonded PTFE/Silicone Slit Septa
Assembled 2mL Clear Vial w/Marking Spot and Bonded PTFE/Silicone Slit Septa
Unassembled 2mL Amber Vial w/Marking Spot and Bonded PTFE/Silicone Septa
Unassembled 2mL Amber Vial w/Marking Spot and Bonded PTFE/Silicone Slit Se
Unassembled 2mL Clear Vial w/Marking Spot and Bonded PTFE/Silicone Septa, 1
Unassembled 2mL Clear Vial High Recovery, w/Marking Spot and Bonded PTFE/
32" Replacement For HP GC 2pk
Replacement PTFE Washers Pack of 2
High Pressure Nut Kit (1 nut) Fits 300/400 series septa
High Pressure Septa Std Kit (nut, 2 high pressure septa)
High Pressure Septa Starter Kit (nut, 1 high-pressure septum)
High Pressure Replacement Septa (1 septum) 400 series
Standard Kit, Nut, 2 Septa
Septa Starter Kit, Nut, 1 Septa
Replacement Septum 300 series (1 septum)
for HP 6890 GCs
20mm Magnetic Crimp Cap with Teflon/Silicone Septum (tan/white) w/8mm Hole,
20mm Magnetic Crimp Cap with Teflon/Silicone Septum (tan/white) w/8mm Hole,
20mm Magnetic Crimp Cap with Teflon/Silicone Septum (tan/white) w/5mm Hole,
20mm Magnetic Crimp Cap with Teflon/Silicone Septum (tan/white) w/5mm Hole,
20mm Magnetic Crimp Cap with PTFE/Butyl Septum w/8mm Hole, 100pk
20mm Magnetic Crimp Cap with PTFE/Butyl Septum w/8mm Hole, 1000pk
Discontinue when stock is depleted 20mm Magnetic Crimp Cap with PTFE/Butyl S
Leak Detector. Includes Storage/Carry Case and Universal Battery Charger Set (L
Check: Probe, Internal Tubing and Batteries Replaced, Unit Tested and Certified,
5ul
10ul
20ul
25ul
50ul
100ul
250ul
500ul
1ml
2ml
5ml
Autosampler Turret Tray for the Agilent 5890/6890 GC
for 12 x 32mm Vials, 50 Vial Capacity, Polypropylene, ea
for 12 x 32mm Vials, 50 Vial Capacity, Polypropylene, 5-pk
Silicone Caps, 10-pk
1/16" Union
1/8" Union
1/4" Union
1/16" Tee
1/8" Tee
1/4" Tee
1/8" to 1/16" Reducer
1/4" to 1/16" Reducer
1/4" to 1/8" Reducer
1/8" Cross
1/4" Cross
1/8" Elbow
1/4" Elbow
1/8" Plug
1/4" Plug
1/8" to 1/16" Tube End Reducer
1/4" to 1/16" Tube End Reducer
1/8" to 1/4" Tube End Reducer
1/4" to 1/8" Tube End Reducer
1/8" Port Connector
1/4" Port Connector
1/8" to 1/4" Port Connector

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034



**HROMalytic** +61(0)3 9762 2034  
**ECHnology** Pty Ltd  
Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

**Australian Distributors**  
Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)

(Ctrl F)  
> Cat.#s



> Current System AUD Prices