

23986	Syringe Filter, 30mm	0.22um PVDF Coffee Luer Lock 100-pk
23987	Syringe Filter, 30mm	0.45um PVDF Red Luer Lock 100-pk
23987	Syringe Filter, 30mm	0.45um PVDF Red Luer Lock 100-pk
23987	Syringe Filter, 30mm	0.45um PVDF Red Luer Lock 100-pk
23988	Syringe Filter, 30mm	0.22um PES White Luer Lock 100-pk
23988	Syringe Filter, 30mm	0.22um PES White Luer Lock 100-pk
23988	Syringe Filter, 30mm	0.22um PES White Luer Lock 100-pk
23989	Syringe Filter, 30mm	0.45um PES Blue Luer Lock 100-pk
23989	Syringe Filter, 30mm	0.45um PES Blue Luer Lock 100-pk
23989	Syringe Filter, 30mm	0.45um PES Blue Luer Lock 100-pk
23993	Q-sep Tube Rack	For 15mL Centrifuge Tube Holds 60, Polypropylene, White
23994	Q-sep Tube Rack	For 50mL Centrifuge Tube Holds 24, Polypropylene, White
23995	Q-sep Tube Rack	For 2mL Centrifuge Tube Holds 100, Polypropylene, White
23996	ASE Cell Carrier	Wire Basket 9"w x 5"h x 7"d
23997	Q-sep Teflon Tube	50mL Centrifuge Tube Teflon FEP, 2-pk
23998	ASE Cell Organizer	blue 13-bin unit, 12"L x 12"H x 7.5"D
23999	ASE Cell Cleaning Brushes	For ASE 100, 200 and 300 Cells 3-pk fa
24002	SPE Vacuum/Pressure	Station, 115 volt.
24003	SPE vacuum stat.220V,EA	
24004	SPE C18-47mmdisks,20pk	
24008	Replacement Vacuum	Valve and Gauge Assembly.
24011	Vacuum Manifold, 12 port	2pk, replacement lid gasket
24014	Waste Container, Resprep	Replacement Waste Container for 12-port Manifold
24016	Vacuum tubing,1/4"ID	10ft roll
24017	Vac.manifold,6D	Diskcover tubing adapter 1/8"OD 6pk
24020	Diskcover-47	EA, SPE 47mm disk holders
24021	Diskcover-47	6pk, SPE 47mm disk holders
24025	Vac.manifold,6D	6pk,Tef. Sample tubing,(2' length, 1/8"OD)
24027	Diskcover-47	EA, Teflon O-ring 47mm
24028	Diskcover-47	6pk, Teflon O-ring 47mm
24029	Diskcover-47	EA, Sample reservoir,125mL
24030	Diskcover-47	6pk, Sample reservoir,125mL
24031	SPE Florisil, 3mL/500mg	50pk
24032	SPE Florisil, 3mL/500mg	50pk, Teflon Frit
24034	SPE Florisil, 6mL/1000mg	30pk
24035	SPE Silica, 3mL/500mg	50pk
24036	SPE Silica, 3mL/500mg	Teflon Frit, 50-pk
24038	SPE Silica, 6mL/1000mg	30pk
24048	SPE C8-47mmdisks,24pk	
24049	SPE OrganoTin, 16pk	60mL, mixed bed
24050	SPE C18,3mL/500mg, 50pk	
24051	SPE C18,6mL/1g, 30pk	
24052	SPE C18,6mL/500mg, 30pk	
24053	SDVB Resin	Equivalent to XAD-2, Untreated in 1kg Containers
24054	TDU Tubes, SS	Empty, 10pk
24055	TDU Tubes, Glass	Empty, 10pk
24056	TDU Tubes, SS	Tenax TA, Unconditioned, 10pk
24057	TDU Tubes, SS	Graphitized Carbon, Unconditioned, 10pk
24058	TDU Tubes, SS	Tenax GR/Carbopack B Unconditioned, 10pk
24059	TDU Tubes, SS	Carbopack B/Carbosieve SIII Unconditioned, 10pk
24060	TDU Tubes, SS	Tenax TA/Graphitized Carbon/Carboxen 1000, Unconditioned, 10pk
24061	TDU Tubes, SS	Carbopack C/Carbopack B Carbosieve SIII, Unconditioned 10pk
24062	TDU Tubes, Glass	Tenax TA, Unconditioned, 10pk
24063	TDU Tubes, Glass	Graphitized Carbon, Unconditioned, 10pk
24064	TDU Tubes, Glass	Tenax GR/Carbopack B, Unconditioned, 10pk
24065	TDU Tubes, Glass	Carbopack B/Carbosieve SIII, Unconditioned, 10pk
24066	TDU Tubes, Glass	Tenax TA/Graphitized Carbon/Carboxen 1000, Unconditioned, 10pk
24067	TDU Tubes, Glass	Carbopack C/Carbopack B Carbosieve SIII, Unconditioned 10pk
24068	TDU Tube Accessories	1/4in Brass Cap and PTFE Ferrule, 20pk
24069	TDU Tube Accessories	1/4in PTFE Ferrule, 20pk
24070	TDU Tube Accessories	CapLok Tool
24071	TDU Tube Accessories	Pen Clip, 10pk
24072	TDU Tube Accessories	TubeMate Tool
24073	TDU Tube Accessories	1/4in SS Union and PTFE Ferrule, 10pk
24074	TDU Tube Accessories	Diffusion Cap, 10pk
24075	Calibration TDU Tubes	SS, Tenax TA 1cm Bed 10pk

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

 HROMalytic +61(0)3 9762 2034 ECHnology Pty Ltd		Australian Distributors Importers & Manufacturers www.chromtech.net.au
Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA		

[Ctrl F]
> Cat.#s

> Current System AUD Prices

