

24076 Calibration TDU Tubes
24077 Calibration TDU Tubes
24078 Calibration TDU Tubes
24079 Calibration TDU Tubes
24080 TDU Tubes, SS
24081 TDU Tubes, SS
24082 TDU Tubes, SS
24083 TDU Tubes, SS
24084 TDU Tubes, SS
24085 TDU Tubes, SS
24086 TDU Tubes, Glass
24087 TDU Tubes, Glass
24088 TDU Tubes, Glass
24089 TDU Tubes, Glass
24090 TDU Tubes, Glass
24091 TDU Tubes, Glass
24100 Gauge, Vacuum/Pressure
24104 Gauge, Vacuum/Pressure
24108 Gauge, Vacuum/Pressure
24120 Gauge
24121 Fittings
24123 Air Canister Jacket
24129 Small Cylinder Stand
24130 Sample Cylinder,Sulfinert
24130-PI Sample Cylinder,Sulfinert
24131 Sample Cylinder,Sulfinert
24131-PI Sample Cylinder
24132 Sample Cylinder,Sulfinert
24132-PI Sample Cylinder
24133 Sample Cylinder,Sulfinert
24133-PI Sample Cylinder
24134 Sample Cylinder,Sulfinert
24134-PI Sample Cylinder
24140 SilcoCan Canister
24140-650 SilcoCan Canister
24140-65065 SilcoCan Canister
24140-65065 SilcoCan Canister
24140-651 1 Liter SilcoCan
24140-652 SilcoCan Canister
24141 SilcoCan Canister
24141-650 SilcoCan Canister
24141-65065 SilcoCan Canister
24141-651 SilcoCan Canister
24142 SilcoCan Canister
24142-650 SilcoCan Canister
24142-65065 SilcoCan Canister
24142-65065 SilcoCan Canister
24142-65065 SilcoCan Canister
24142-651 SilcoCan Canister
24142-652 SilcoCan Canister
24143 SilcoCan Canister
24143-650 SilcoCan Canister
24143-65065 SilcoCan Canister
24143-65065 SilcoCan Canister
24143-65065 SilcoCan Canister
24143-651 SilcoCan Canister
24143-652 SilcoCan Canister
24143-653 SilcoCan Canister
24144 Valve
24145 Valve
24146 Valve
24147 Valve
24148 Canister Valve
24151 Air Canister Tripod
24160 Sampling Kit
24161 Sampling Kit

Page 100
Glass, Tenax TA 1cm Bed 10pk
Calibration Solution Loading Rig
9.5mm Replacement Septa 10pk, for Calibration Solution Loading Rig
10pk Certified Reference Std, SS 100mg BTX on Tenax TA
Tenax TA, Conditioned and Capped, 10pk
Graphitized Carbon, Conditioned and Capped, 10pk
10pk, Tenax GR/Carbopack B Conditioned and Capped
10pk, Carbopack B/Carbosieve SIII Conditioned and Capped
10pk, Tenax TA/Graphitized Carbon/Carboxen 1000, Conditioned and Capped
10pk, Carbopack C/Carbopack B Carbosieve SIII, Conditioned and Capped
Tenax TA, Conditioned and Capped, 10pk
10pk, Graphitized Carbon, Conditioned and Capped
10pk, Tenax GR/Carbopack B Conditioned and Capped
10pk, Carbopack B/Carbosieve SIII Conditioned and Capped
10pk, Tenax TA/Graphitized Carbon/Carboxen 1000 Conditioned and Capped
10pk, Carbopack C/Carbopack B Carbosieve SIII, Conditioned and Capped
30" to 15 psig compound 1/8" NPT CBM, 2" diam, SS
30" to 30 psig Compound 1/8" NPT CBM, 2" diam, SS
30" to 60 psig Compound 1/8"NPT CBM, 2" diam, SS
30" to 0 psig Vacuum, 1/8" NPT CBM, 1.5" Diam. 316 SS, Lexan Lens
Quick Connect Stem Protector, SS, For Miniature Air Sampling Canister
Heating Jacket, 110 Volt
Support for Disposable 8cm Diameter Gas Cylinders
75cc, 1/4" Female NPT Threads
TPED, 75cc, 1450psig, 304L SS, 1/4" Female NPT Ends
150cc, 1800psig, 304L, SS, 1/4" Female NPT Ends
Sulfinert, TPED, 150cc, 1450psig, 304L, SS, 1/4" Female NPT Ends
300cc, 1800psig, 304L, SS, 1/4" Female NPT Ends
Sulfinert, TPED, 300cc, 1450psig, 304L, SS, 1/4" Female NPT Ends
500cc, 1800psig, 304L, SS, 1/4" Female NPT Ends
Sulfinert, TPED, 500cc, 1450psig, 304L, SS, 1/4" Female NPT Ends
1000cc, 1800psig, 304L, SS, 1/4" Female NPT Ends
Sulfinert, TPED, 1000cc, 1450psig, 304L, SS, 1/4" Female NPT Ends
1 Liter 1/4" Valve w/Gauge
1 Liter w/Gauge and Siltek Treated 1/4" Valve
1 Liter w/Gauge and 1/4" Silcosteel Valve 2" McDaniels Vacuum Pressure Gauge
1 Liter w/Gauge and Siltek Treated 1/4in Valve ship under vacuum
W/ 30 -15 vacuum gauge
1 Liter 1/4" Valve w/2" McDaniels Vacuum Pressure Gauge 30" HG to 30 PSI
3 Liter 1/4" Valve w/Gauge (cat #563646 is needed if connecting to a sampling kit)
3 Liter w/Gauge and Siltek Treated 1/4" Valve (cat #563646 is needed if connectir
3 Liter w/Gauge and 1/4" Silcosteel Valve Ship under vacuum (cat #563646 is nee
3 Liter 1/4" Valve 2" McDaniels Vacuum Pressure Gauge 30" HG to 15 PSI (cat #
6 Liter 1/4" Valve w/Gauge
6 Liter w/Gauge and Siltek Treated 1/4" Valve
6 Liter, 1/4" Silcosteel Valve, 2" McDaniels Pressure Gauge 30HG to 15 PSI
6Liter 1/4" SilcoSteel Valve w/2" McDaniels Vac Pres Gauge 30"HG to 30PSI
6 Liter 1/4" Silcosteel Valve w/Gauge Ship under Vacuum
6 Liter 1/4" Valve w/Gauge 2" McDaniels Vacuum Pressure Gauge 30"HG to 15 P
6 Liter 1/4" Valve 2" McDaniels Vacuum Pressure Gauge 30" HG to 30 PSI
15 Liter 1/4" Valve w/Gauge
15 Liter w/Gauge and Siltek Treated 1/4" Valve
15 Liter 1/4" Silcosteel Valve 2" McDaniels Vacuum Pressure Gauge 30" HG to 15
15 Liter 1/4" Silcosteel Valve 2" McDaniels Vacuum Pressure Gauge 30" HG to 30
15 Liter 1/4" Silcosteel Valve Ship under vacuum
15 Liter 1/4" Valve 2" McDaniels Vacuum Pressure Gauge 30" HG to 15" PSI
15 Liter 1/4" Valve 2" McDaniels Vacuum Pressure Gauge 30" HG to 30 PSI
15 Liter 1/4" Valve w/Gauge Ship under vacuum
2-Port Siltek Replacement 1/4" Valve for Air Monitoring Canisters
2-Port Replacement 1/4" Valve for Air Monitoring Canisters
3-Port Siltek Replacement 1/4" Valve for Air Monitoring Canisters
3-Port Replacement 1/4" Valve for Air Monitoring Canisters
SS4H Bellows Replacement Valve 2-port only, Stainless Steel only
Holds two 1L, 3L, or 6L size canisters
CONFIRM DELIVERY WITH OPERATIONS > 25 Siltek, Complete Kit .0012 Orific
CONFIRM DELIVERY WITH OPERATIONS > 25 Siltek, Complete Kit .0016 Orific

(Ctrl F)
> Cat.#s

> Current System AUD Prices



Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

**HROMalytic** +61(0)3 9762 2034
ECHnology Pty Ltd
Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au